



Product Certificate

Jørgen Kruuse A/S hereby declare the following details regarding the mentioned product

Item number

291015

Product name

Fasol.

**Product type**

Laboratory equipment for sampling.

Flotation medium used to easy and hygienic test of dog or cat for intestinal parasites.

To be used in conjunction with

- 291010 Fecalyzer or
- 291011 KRUUSE Fecal Ova 3 Step.

Supplementary items

291010



291011

Material

Aqua, Magnesiumsulfate & Chloro solution (preservative).

Colour

Clear watery solution.

Storage condition

Keep containers tightly and closed.

Instruction for use (with use of 291010)

1. Lift cap but do not remove green insert. Fill green vial with fasol to the tip of the arrow embossed on side of vial.
2. Rotate green insert vial back and forth to separate ova from fecal sample. Mix thoroughly.
3. Seat green insert vial firmly in place with tongue depressor (or with thumbs).
4. Fill holder completely to form meniscus with additional fasol.
5. Float 22mm cover slip on meniscus for 15-20 minutes.
6. Transfer cover slip to slide for microscopic examination at 100x magnification.
7. Close cap and dispose of fecalyzer.

Shelf-life

5 years.

Size

1 litre.

Legislation, EEC directives, Standards

Material safety datasheet is available.

Origin

DK

Last revised

24-nov-2011

Version

2

Signed

Inge Merete Jensen

Function

QA Manager