

DIGITAL SOLUTIONS

PRODUCT REFERENCES





Virginia Woolf, 17
25005 · Lleida (Spain)

INTERNATIONAL +34 873 450 709
das@dynamicabutment.com

SPAIN +34 973 289 580
spain@dynamicabutment.com

Dynamic Abutment Solutions is the trademark of the manufacturer Talladium España S.L. The information included in this catalogue is exclusively addressed to professionals in dental sector.

All commercial trademarks mentioned herein are fully registered by their respective companies, and the images that appear are just to provide an orientation. Talladium reserves the right to modify this document without previous notice. Talladium is not responsible for the inadequate execution of these products if the warning indications corresponding to every reference are not contemplated.

Ask about the products available with CE legislation. Some of the products are not authorized for sale and distribution or do not have a sales license in some countries according to other legislations. Information: das@dynamicabutment.com.



INDEX

8	LIST OF COMPATIBILITIES AVAILABLE
9	COMPATIBILITIES AVAILABLE
24	COMPATIBILITIES AND REFERENCES
169	DAS LIBRARIES
171	MULTI-UNIT DAS LIBRARIES
173	DYNAMIC SCREWS TECHNICAL SPECIFICATIONS
200	STRAIGHT SCREWS TECHNICAL SPECIFICATIONS
216	SCREWDRIVERS AND STRAIGHT SCREWS
217	DYNAMIC MILLING TOOL SPECIFICATIONS
221	SCREWDRIVERS AND COMPLEMENTS
222	DAS MU SYSTEM COMPONENTS

MULTI-UNIT
DAS SYSTEM

DYNAMIC
DAS SYSTEM

LIST OF COMPATIBILITIES AVAILABLE

AB	BIOTEC	EASY IMPLANT	JDENTALCARE	NOVA IMPLANTS	STRAUMANN
ACE	BIOTEM	ECKERMANN	KEYSTONE	OSSTEM IMPLANT	SURCAM DENTAL
ADIN	BREIDENT MEDICAL	ELITE MEDICA	KLOCKNER	OSTEOPLUS	SYBRON IMPLANT SOLUTIONS
ALFA-GATE	BTI	EUROTEKNIKA	KUWOTECH	OXY	SYSTHEX
ALPHABIO	BTK	F&B IMPLANT (FIT & BRILLIANT)	LASAK	PALTOP	TBR
ALPHA-DENT	CAMLOG	GALIMPLANT	LEADER	PHIBO	TITANIUM-FIX
ANCLADEN	CONEXÃO SISTEMA DE PRÓTESE	GC TECH	MEDENTIKA	POINT IMPLANT	TREE-OSS
ANKYLOS	CORTEX	GMI (ILERIMPLANT)	MEDENTIS	PROCLINIC	TRI DENTAL IMPLANTS
ANTHOGYR	COWELLMEDI	GT MEDICAL	MEGAGEN	PROTEG IMPLANTS	TRINON
ARDS	C-TECH	HAHN IMPLANT (GLIDEWELL)	MICRODENT	RADHEX	UFIT
ASTRA	DENTAL TECH	HIOSSEN	MIS	REFLECT	VULKAN IMPLANTS
AVINENT	DENTAURUM	HI-TEC	MONOIMPLANT	RITTER	WARANTEC (ONEPLANT)
B&W	DENTEGRIS	IBS	MOZO-GRAU (TICARE)	ROOTT	WIN
BEGO	DENTEM	IDO IMPLANTS	MPI	SEWON MEDIX	XIVE
BIOCONCEPT	DENTIS	IHDE DENTAL (IMBIODENT)	NEOBIOTECH	SIC INVENT	YES IMPLANT
BIOGENESIS	DENTIUM	IMPLANT DIRECT	NEODENT	SIGNO VINCES	ZIACOM
BIOHORIZONS	DIO IMPLANTS	IMPLANT GENESIS	NEOSS	SIN IMPLANTS	ZIMMER
BIOLOK	DITRON	IMPLANTSWISS	NOBEL BIO CARE	SOUTHERN IMPLANTS	
BIOMET 3i	DMI DENTAL SUPPLY	INTRA-LOCK	NORIS MEDICAL	STERI-OSS	
BIONER	DSP BIOMEDICAL		NORMON	STERNGOLD	

COMPATIBILITIES AVAILABLE

AB

Model: I2
Implant Ø: 3.5/3.75/4.2/4.5/
Platform: Standard
Code: 0040

Model: I2
Implant Ø: 5/6
Platform: Standard
Code: 0040

Model: I22
Implant Ø: 3.75/4.22
Platform: Standard
Code: 0040

Model: I5
Implant Ø: 3.5/3.75/4.2/4.5/
Platform: Standard
Code: 0040

Model: I5
Implant Ø: 5/6/7/8
Platform: Standard
Code: 0040

Model: I55
Implant Ø: 3.75/4.2/4.5
Platform: Standard
Code: 0040

Model: I55
Implant Ø: 5/6/7/8
Platform: Standard
Code: 0040

Model: I10
Implant Ø: 4.2/5
Platform: Standard
Code: 0040

Model: I15
Implant Ø: 6/7/8
Platform: Standard
Code: 0040

Model: Multi Unit D1-P64
Implant Ø: Multi Unit D1-P64
Platform: Universal
Code: 0025

ACE

Model: External Hex
Implant Ø: 3.3
Platform: NP 3.5
Code: 0023

Model: External Hex
Implant Ø: 3.75/4
Platform: RP 4.1
Code: 0024

Model: External Hex
Implant Ø: 4.75
Platform: WP 5
Code: 0058

Model: Infinity TRI-CAM
Implant Ø: 3.5
Platform: 3.5
Code: 0026

Model: Infinity TRI-CAM
Implant Ø: 4.3
Platform: 4.3
Code: 0027

Model: Infinity TRI-CAM
Implant Ø: 5
Platform: 5
Code: 0028

Model: Infinity Internal Hex
Implant Ø: 3.7/4.1
Platform: 3.5
Code: 0040

Model: Infinity Internal Hex
Implant Ø: 4.7/5.1
Platform: 4.5
Code: 0041

Model: Infinity Octagon
Implant Ø: 3.3
Platform: RP 4.8
Code: 0037

Model: Infinity Octagon
Implant Ø: 4.1
Platform: RP 4.8
Code: 0037

Model: Infinity Octagon
Implant Ø: 4.8
Platform: RP 4.8
Code: 0037

Model: Infinity Octagon
Implant Ø: 4.8
Platform: WP 6.5
Code: 0096

Model: Multi Unit
Implant Ø: Universal
Platform: Universal
Code: 0025

ADIN

Model: Swell
Implant Ø: 3.3
Platform: 3.45
Code: 0040

Model: Swell
Implant Ø: 3.3
Platform: 3.45
Code: 0042

Model: Swell
Implant Ø: 3.75/4.2
Platform: 3.6
Code: 0040

Model: Swell
Implant Ø: 3.75/4.2
Platform: 3.6
Code: 0042

Model: Swell
Implant Ø: 5
Platform: 4
Code: 0040

Model: Swell
Implant Ø: 5
Platform: 4
Code: 0042

Model: Swell
Implant Ø: 6
Platform: 4.6
Code: 0040

Model: Swell
Implant Ø: 6
Platform: 4.6
Code: 0042

Model: Touareg-S / Touareg-OS
Implant Ø: 3.5
Platform: 3.45
Code: 0040

Model: Touareg-S / Touareg-OS
Implant Ø: 3.5
Platform: 3.45
Code: 0042

Model: Touareg-S / Touareg-OS
Implant Ø: 3.75/4.2
Platform: 3.6
Code: 0040

Model: Touareg-S / Touareg-OS
Implant Ø: 3.75/4.2
Platform: 3.6
Code: 0042

Model: Touareg-S / Touareg-OS
Implant Ø: 5
Platform: 4
Code: 0040

Model: Touareg-S / Touareg-OS
Implant Ø: 5
Platform: 4
Code: 0042

Model: Touareg-S / Touareg-OS
Implant Ø: 6
Platform: 5
Code: 0040

Model: Touareg-S / Touareg-OS
Implant Ø: 6
Platform: 5
Code: 0042

Model: Touareg CloseFit
Implant Ø: 2.75
Platform: UNP
Code: 0188

Model: Touareg CloseFit
Implant Ø: 3
Platform: NP
Code: 0145

Model: Touareg CloseFit
Implant Ø: 3.5
Platform: RP
Code: 0021

Model: Touareg CloseFit
Implant Ø: 4.3/5
Platform: WP
Code: 0022

Model: Triple
Implant Ø: 3.5/3.75/4.2/5/6
Platform: Standard
Code: 0040

Model: Triple
Implant Ø: 3.5/3.75/4.2/5/6
Platform: Standard
Code: 0042

Model: Multi Unit TMA
Implant Ø: Universal
Platform: Universal
Code: 0025

ALFA-GATE

Model: Bioactive/Porous/Trio/MAX
Implant Ø: 3.3/3.75/4.2/4.7/5.2/6
Platform: SP
Code: 0040

Model: Conical
Implant Ø: 3.5
Platform: NP
Code: 0021

Model: Conical
Implant Ø: 4.3/5
Platform: RP
Code: 0022

ALPHABIO

Model: Internal Hex Connection (IH) SPI
Implant Ø: 3.3/3.75/4.2/5/6
Platform: Universal
Code: 0040

Model: Internal Hex Connection (IH) SPI
Implant Ø: 3.7/3.75/4.2/4.65/5.3
Platform: Universal
Code: 0040

Model: Internal Hex Connection (IH) ICE
Implant Ø: 3.3/3.75/4.2/5
Platform: Universal
Code: 0040

Model: Internal Hex Connection (IH) ICE
Implant Ø: 3.3/3.75/4.2/5/6
Platform: Universal
Code: 0040

Model: Internal Hex Connection (IH) DF1
Implant Ø: 3.75/4.2/5
Platform: 3.5
Code: 0040

Model: Conical Hex Connection (CHC) NICE
Implant Ø: 3.2
Platform: Narrow
Code: 0136

Model: Conical Hex Connection (CHC) NICE
Implant Ø: 3.2/3.5
Platform: Narrow
Code: 0136

Model: Conical Standard Connection (CS)
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0169

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0195

COMPATIBILITIES AVAILABLE

ALPHA-DENT

Model: Classic Conus
Implant Ø: 3.3/3.75
Platform: 3
Code: 0265

Model: Active Conus
Implant Ø: 3.25/3.75
Platform: 3
Code: 0265

ANCLADEN

Model: Anclalock
Implant Ø: 3.75/4.25/5
Platform: 3.5
Code: 0040

ANKYLOS

Model: Ankylos
Implant Ø: 3.5
Platform: 3.5
Code: 0075

Model: Ankylos
Implant Ø: 4.5
Platform: 4.5
Code: 0075

Model: Ankylos
Implant Ø: 5.5
Platform: 5.5
Code: 0075

Model: Ankylos
Implant Ø: 7
Platform: 7
Code: 0075

Model: Balance Base Narrow
Implant Ø: Universal
Platform: Universal
Code: 0183

ANTHOGRYR

Model: Axiom BL REG / PX
Implant Ø: 3.4
Platform: 3.4
Code: 0161

Model: Axiom BL REG / PX
Implant Ø: 4
Platform: 4
Code: 0149

Model: Axiom BL REG / PX
Implant Ø: 4.6
Platform: 4.6
Code: 0149

Model: Axiom BL REG / PX
Implant Ø: 5.2
Platform: 5.2
Code: 0162

Model: Anthofit HE
Implant Ø: 3.5/3.75/4
Platform: R (4.1)
Code: 0024

Model: Anthofit HE
Implant Ø: 5
Platform: L (5)
Code: 0058

Model: Ossfit
Implant Ø: 3.5/4.2
Platform: 4.8
Code: 0074

Model: Ossfit
Implant Ø: 3.5/4.2
Platform: 4.8
Code: 0037

Model: Ossfit
Implant Ø: 5
Platform: 6.5
Code: 0096

Model: Multi Unit
Implant Ø: 4.8
Platform: Universal
Code: 0163

ARDS

Model: Smart
Implant Ø: 3.75/4.2/4.5
Platform: 3.75
Code: 0040

Model: Classic
Implant Ø: 3.3/3.75/4.2/5/6
Platform: 3.75
Code: 0040

Model: Premium
Implant Ø: 3.3/3.75/4.2/5/6
Platform: 3.75
Code: 0040

Model: CIT
Implant Ø: 3.3/3.75/4.2/5/6
Platform: 3.75
Code: 0040

ASTRA

Model: Yellow
Implant Ø: 3
Platform: Yellow (X-estricha)
Code: 0109

Model: Aqua
Implant Ø: 3.5/4
Platform: Aqua(Estreicha)
Code: 0004

Model: Lilac
Implant Ø: 4.5/5
Platform: Lilac (Ancha)
Code: 0005

Model: Uniabutment Cono 200
Implant Ø: Regular/Wide
Platform: Regular/Wide
Code: 0066

Model: Evolution (Interno)
Implant Ø: 3
Platform: 3
Code: 0090

Model: Evolution (Interno)
Implant Ø: 3.6
Platform: 3.6
Code: 0006

Model: Evolution (Interno)
Implant Ø: 4.2
Platform: 4.2
Code: 0007

Model: Evolution (Interno)
Implant Ø: 4.8
Platform: 4.8
Code: 0091

Model: Evolution (Interno)
Implant Ø: 5.4
Platform: 5.4
Code: 0092

Model: Evolution
(Cono externo 33°)
Implant Ø:
Platform: Universal
Code: 0008

Model: Multibase Abutment
(SmartFix concept)
Implant Ø:
Platform: Universal
Code: 0258

Model: DS OmniTaper EV
Implant Ø: 3
Platform: XS
Code: 0090

Model: DS OmniTaper EV
Implant Ø: 3.4
Platform: S
Code: 0006

Model: DS OmniTaper EV
Implant Ø: 3.8
Platform: M
Code: 0007

Model: DS OmniTaper EV
Implant Ø: 4.5
Platform: L
Code: 0091

Model: DS OmniTaper EV
Implant Ø: 5.5
Platform: XL
Code: 0092

AVINENT

Model: HE/EC
Implant Ø: 3.3//3.5/4
Platform: 3.5
Code: 0023

Model: HE/EC
Implant Ø:
3.3/3.8/4/4.2/4.8//4.5/5
Platform: 4.1
Code: 0024

Model: HE/EC
Implant Ø: 4.8
Platform: 5.1
Code: 0061

Model: HI/IC
Implant Ø: 3.1//3.5/4
Platform: 3.5
Code: 0040_B

Model: HI/IC
Implant Ø:
3.3/3.8/4/4.2/4.8//4.5/5
Platform: 4.1
Code: 0040_B

Model: Transepitelial
Implant Ø:
Platform: Regular
Code: 0025

B&W

Model: Hexágono Externo
Implant Ø: 3.75/4
Platform: 4.1
Code: 0024

Model: Hexágono Externo
Implant Ø: 5
Platform: 5
Code: 0058

Model: Cónico Hexágono
Interno CIH
Implant Ø: 3.3/4
Platform: 4
Code: 0040

Model: Cónico Hexágono
Interno CIH
Implant Ø: 3.3/4
Platform: 4
Code: 0042

BEGO

Model: RS/RSX
Implant Ø: 3.0
Platform: 3.0
Code: 0049

Model: S/RI/RS/RSX
Implant Ø: 3.25/3.75
Platform: 3.67
Code: 0050

Model: S/RI/RS/RSX
Implant Ø: 4.1
Platform: 4.1
Code: 0051

Model: S/RI/RS/RSX
Implant Ø: 4.5
Platform: 4.5
Code: 0052

Model: S/RI/RS/RSX
Implant Ø: 5.5
Platform: 5.5
Code: 0081

Model: Mini
Implant Ø: 2.7/2.9/3.1
Platform: Mini
Code: 0187

Model: MultiPlus
Implant Ø:
Platform: Universal
Code: 0150

BIOCONCEPT

Model: BC Tissue Level Standard
Implant Ø: 3.3/4.1/4.8
Platform: Regular
Code: 0037

COMPATIBILITIES AVAILABLE

Model: BV Tapered Bone Level
Implant Ø: 3.5
Platform: Narrow
Code: 0029

Model: BV Tapered Bone Level
Implant Ø: 4/4.5/5
Platform: Regular
Code: 0030

Model: BC Bone Level
Implant Ø: 4.1/4.8
Platform: Regular
Code: 0035

Model: BC Tissue Level Standard Plus
Implant Ø: 4.8
Platform: Regular
Code: 0037

Model: BC Tissue Level Tapered Effect
Implant Ø: 4.8
Platform: Regular
Code: 0037

Model: BC Bone Level
Implant Ø: 3.3
Platform: Narrow
Code: 0033

BIOGENESIS

Model: 3icon
Implant Ø: 3.3
Platform: Mini (Pink)
Code: 0023

Model: 3icon
Implant Ø: 3.75/4/4.3/4.5
Platform: Regular (Blue)
Code: 0024

Model: 3icon
Implant Ø: 5/5.5
Platform: Wide (Yellow)
Code: 0058

Model: Aticon
Implant Ø: 3.5/4/4.5/5
Platform: Blue
Code: 0005

Model: Iticon
Implant Ø: 3.5/4.1/4.8
Platform: 4.8
Code: 0037

BIOHORIZONS

Model: External
Implant Ø: 3.5
Platform: 3.7 (Yellow)
Code: 0023

Model: Tapered Internal
Implant Ø: 3/3.4
Platform: 3 (Grey)
Code: 0102

Model: Tapered Internal
Implant Ø: 3.8
Platform: 3.5 (Yellow)
Code: 0040

Model: Tapered Internal
Implant Ø: 4.6
Platform: 4.5 (Green)
Code: 0041

Model: Tapered Internal
Implant Ø: 5.8
Platform: 5.7 (Blue)
Code: 0080

Model: Tapered Plus
Implant Ø: 3.8
Platform: 3 (Grey)
Code: 0102

Model: Tapered Plus
Implant Ø: 4.6
Platform: 3.5 (Yellow)
Code: 0040

Model: Tapered Plus
Implant Ø: 5.8
Platform: 4.5 (Green)
Code: 0041

Model: Tapered Pro Conical
Implant Ø: 3.3/3.8
Platform: Narrow (Grey)
Code: 0119

Model: Tapered Pro Conical
Implant Ø: 4.2/4.6/5.2
Platform: Regular (Yellow)
Code: 0120

Model: Mount-free
Tapered Internal
Implant Ø: 3/3.4
Platform: 3 (Grey)
Code: 0102

Model: Mount-free
Tapered Internal
Implant Ø: 3.8
Platform: 3.5 (Yellow)
Code: 0040

Model: Mount-free
Tapered Internal
Implant Ø: 4.6
Platform: 4.5 (Green)
Code: 0041

Model: Mount-free
Tapered Internal
Implant Ø: 5.8
Platform: 5.7 (Blue)
Code: 0080

Model: Tapered Pro
Implant Ø: 3.8
Platform: 3 (Grey)
Code: 0102

Model: Tapered Pro
Implant Ø: 4.2/4.6
Platform: 3.5 (Yellow)
Code: 0040

Model: Tapered Pro
Implant Ø: 5.2
Platform: 4.5 (Green)
Code: 0041

Model: Tapered Short
Implant Ø: 4.6
Platform: 3.5 (Yellow)
Code: 0040

Model: Tapered Short
Implant Ø: 5.8
Platform: 4.5 (Green)
Code: 0041

Model: Tapered PTG
Implant Ø: 4.2
Platform: 3.5 (Yellow)
Code: 0040

Model: Tapered IM
(Immediate Molar)
Implant Ø: 7/8
Platform: 5.7 (Blue)
Code: 0080

Model: Tapered Tissue Level
Implant Ø: 3/3.8
Platform: 3.5 (Yellow)
Code: 0040

Model: Tapered Tissue Level
Implant Ø: 4.6
Platform: 4.5 (Green)
Code: 0041

Model: Tapered Tissue Level
Implant Ø: 5.8
Platform: 5.7 (Blue)
Code: 0080

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0025

BIOLOK

Model: Hexágono Externo
Implant Ø: 3.45
Platform: 3.45
Code: 0003

BIOMET 3i

Model: Osseotite External Hex
Implant Ø: 3.25
Platform: 3.4
Code: 0003

Model: Osseotite External Hex
Implant Ø: 3.75/4
Platform: 4.1
Code: 0024

Model: Osseotite External Hex
Implant Ø: 5
Platform: 5
Code: 0058

Model: Certain
Implant Ø: 3.25/4
Platform: 3.4
Code: 0001

Model: Certain
Implant Ø: 4/5
Platform: 4.1
Code: 0002

Model: Certain
Implant Ø: 5
Platform: 5
Code: 0057

Model: Low Profile
Implant Ø:
Platform: Universal
Code: 0025

BIONER

Model: Ikelt / Bikelt
Implant Ø: 3.3/3.75/4
Platform: 4.1
Code: 0024

Model: Ikelt
Implant Ø: 5
Platform: 5
Code: 0058

Model: Hikelt
Implant Ø: 3.8
Platform: 3.95
Code: 0040

Model: Hikelt
Implant Ø: 4.7
Platform: 4.9
Code: 0041

Model: TopDM
Implant Ø: 3.5
Platform: 3.5
Code: 0021

Model: TopDM
Implant Ø: 4
Platform: 4
Code: 0021

Model: TopDM
Implant Ø: 3.25/4
Platform: 5
Code: 0021

Model: Short DM
Implant Ø: 4/5/6
Platform: Universal
Code: 0021

Model: Hibikelt
Implant Ø: 4/5
Platform: Universal
Code: 0021

Model: Transepitelial A-5M
Implant Ø: Transepitelial A-5M
Platform: Regular
Code: 0025

BIOTEC

Model: SPR/CIM
Implant Ø: 3.3
Platform: 3.3
Code: 0040

Model: SPR/CIM
Implant Ø: 3.75
Platform: 3.75
Code: 0040

Model: SPR/SPTT/CIM
Implant Ø: 4.2
Platform: 4.2
Code: 0040

Model: SPR/SPTT/CIM
Implant Ø: 5
Platform: 5
Code: 0040

COMPATIBILITIES AVAILABLE

BIOTEM

Model: AR Fixture
Implant Ø: 3.7/4/4.5
Platform: Regular
Code: 0030

BREDENT MEDICAL

Model: Narrow Sky
Implant Ø: 3.5
Platform: NP 3.5
Code: 0110

Model: Blue Sky
Implant Ø: 3.5/4/4.5/5.5
Platform: 4
Code: 0111

Model: Blue Sky Classic
Implant Ø: 3.5/4/4.5
Platform: 4
Code: 0111

Model: Copa Sky
Implant Ø: 3.5/4/4.5/5/6
Platform: 3.3
Code: 0251

BTI

Model: Externa Tiny
Implant Ø: 2.5/3/3.3/3.5/3.75
Platform: Tiny 3.5
Code: 0009

Model: Externa Universal
Implant Ø: 3.75/4
Platform: Universal 4.1
Code: 0024

Model: Externa Universal Plus
Implant Ø: 4.5/5
Platform: Universal Plus 4.1
Code: 0024

Model: Externa
Implant Ø: 4.5/5/5.5
Platform: Ancha 5.5
Code: 0060

Model: Interna
Implant Ø: 3.3/3.5/3.75
Platform: 3.5
Code: 0257

Model: Interna Universal
Implant Ø: 3.3/3.5/3.75/4/4.25
Platform: Universal 4.1
Code: 0010

Model: Interna Universal Plus
Implant Ø: 4.5/5/5.5
Platform: Universal Plus 4.1
Code: 0010

Model: Interna Ancha
Implant Ø: 5.5/6/6.25
Platform: Ancha 5.5
Code: 0059

Model: Multi-IM
Implant Ø: converter 4.1
Platform: Universal 4.1
Code: 0151

Model: Multi-IM
Implant Ø: converter 5.5
Platform: Ancha 5.5
Code: 0177

BTK

Model: Klassic / Konic
Implant Ø: 3.25
Platform: 3.4 EN
Code: 0003

Model: Klassic / Konic/Line Plus
Implant Ø: 3.25PL/3.75/4
Platform: 4.1 ER
Code: 0024

Model: IS +/Line Plus IS+
Implant Ø: 3.3/3.7/4.1/4.8/6
Platform: DR
Code: 0029

Model: Klassic / Konic
Implant Ø: 3.3/3.7/4.1/4.8
Platform: KR
Code: 0029

Model: Klassic / Konic
Implant Ø: 3.25/4
Platform: 3.5 IR
Code: 0040

Model: Line Plus IC+
Implant Ø: 3.25/3.75/4.25/5
Platform: LR
Code: 0040

Model: Klassic / Konic
Implant Ø: 3.25/4
Platform: 3.5 IR
Code: 0042

CAMLOG

Model: Camlog Screw-Line/
Progressive-Line
Implant Ø: 3.3
Platform: 3.3
Code: 0087

Model: Camlog Screw-Line/
Progressive-Line
Implant Ø: 3.8
Platform: 3.8
Code: 0011

Model: Camlog Screw-Line/
Progressive-Line
Implant Ø: 4.3
Platform: 4.3
Code: 0012

Model: Camlog Screw-Line/
Progressive-Line
Implant Ø: 5
Platform: 5
Code: 0088

Model: Conelog Screw-Line/
Progressive-Line
Implant Ø: 3.3
Platform: 3.3
Code: 0119

Model: Conelog Screw-Line/
Progressive-Line
Implant Ø: 3.8
Platform: 3.8
Code: 0120

Model: Conelog Screw-Line/
Progressive-Line
Implant Ø: 4.3
Platform: 4.3
Code: 0121

CONEXÃO SISTEMA DE PRÓTESE

Model: NP 24° Flesh Gold
Implant Ø: 3.5/3.75/4.3/5
Platform: Universal
Code: 0021

Model: NP 24° Flash
Implant Ø: 3.5/4.3/5
Platform: Universal
Code: 0021

Model: NP 24° Pterytrans EX
Implant Ø: 3.5
Platform: Universal
Code: 0021

Model: NP 24° Short® NP BLT
Implant Ø: 3.75/4.3
Platform: Universal
Code: 0021

Model: NP 24° Torq
Implant Ø: 3.5/3.75/4
Platform: Universal
Code: 0021

Model: NP 24° Expand
Implant Ø: 3.75/4/5
Platform: Universal
Code: 0021

CORTEX

Model: Internal Hex Dynamix
Implant Ø: 3.3/3.8/4.2/5/6
Platform: 3.75
Code: 0040

Model: Internal Hex Classix
Implant Ø: 3.3/3.8/4.2/5/6
Platform: 3.75
Code: 0040

Model: Internal Hex Saturn
Implant Ø: 3.8/4.2
Platform: 3.5
Code: 0040

Model: Conical Dynamix
Implant Ø: 3
Platform: NP
Code: 0109

Model: Conical Dynamix
Implant Ø: 3.3/3.8/4.2
Platform: RP
Code: 0004

Model: Conical Dynamix
Implant Ø: 5/6
Platform: WP
Code: 0005

Model: Conical Classix
Implant Ø: 3.3/3.8/4.2
Platform: RP
Code: 0004

Model: Conical Classix
Implant Ø: 5/6
Platform: WP
Code: 0005

Model: Conical Magix
Implant Ø: 3.3/3.8/4.2
Platform: RP
Code: 0004

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0025

COWELLMEDI

Model: Inno - External Type
Implant Ø: 5.0/6.0
Platform: 5.1
Code: 0061

Model: Multi S&A Abutment Ø 4.5 mm
Implant Ø: Multi S Abutment
Platform: Universal
Code: 0193

C-TECH

Model: EL Esthetic Line
Implant Ø: 3.8/4.3/5.1
Platform: 4
Code: 0246

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0245

DENTAL TECH

Model: Implogix
Implant Ø: 4.5
Platform: 4.5 (Blue)
Code: 0041

DENTAURUM

Model: Tiologic
Implant Ø: 3.3
Platform: Small
Code: 0130

Model: Tiologic
Implant Ø: 3.7/4.2
Platform: Medium
Code: 0131

Model: Tiologic
Implant Ø: 4.8/5.5
Platform: Large
Code: 0132

DENTEGRIS

Model: SLS-Straight
Implant Ø: 4.5
Platform: 4.5
Code: 0041_B

Model: Sinus-Lift
Implant Ø: 4.5
Platform: 4.5
Code: 0041_B

COMPATIBILITIES AVAILABLE

Model: S&T Implants
Implant Ø: 4.5
Platform: 4.5
Code: 0041_B

DENTEM

Model: Regular
Implant Ø:
Platform: Regular
Code: 0030

DENTIS

Model: OneQ-SL
Implant Ø: 3
Platform: Narrow
Code: 0014

Model: OneQ-SL
Implant Ø: 3.9/4.2/4.7/5.2
Platform: Regular
Code: 0030

Model: OneQ-SL
Implant Ø: 6/7/8
Platform: Wide
Code: 0030

Model: s-Clean Tapered /
Tapered II
Implant Ø: 3.7
Platform: Mini
Code: 0030

Model: s-Clean Tapered /
Tapered II
Implant Ø: 4.1/4.3
Platform: Regular
Code: 0030

Model: s-Clean Tapered /
Tapered II
Implant Ø: 4.8
Platform: Wide
Code: 0030

Model: s-Clean Straight
Implant Ø: 4.1/4.8
Platform: 4.1/4.8
Code: 0030

Model: s-Clean Save
Implant Ø: 5.5/6
Platform: 5.5/6
Code: 0030

Model: SQ-SL
Implant Ø: 3.5
Platform: Narrow
Code: 0014

Model: SQ-SL
Implant Ø: 4/4.5/5
Platform: Regular
Code: 0030

Model: SQ-SL
Implant Ø: 6/7/8
Platform: Wide
Code: 0030

Model: e-Clean
Implant Ø: 3.5
Platform: Mini
Code: 0023

Model: e-Clean
Implant Ø: 4.1
Platform: Regular
Code: 0024

Model: e-Clean
Implant Ø: 5.1
Platform: Wide
Code: 0061

Model: i-Clean Tapered
Implant Ø: 3.7
Platform: Mini
Code: 0037

Model: i-Clean Tapered
Implant Ø: 4.1/4.3
Platform: Regular
Code: 0037

Model: i-Clean Tapered
Implant Ø: 4.8
Platform: Wide
Code: 0037

Model: i-Clean Straight
Implant Ø: 4.1/4.8
Platform: 4.8
Code: 0037

Model: Octa Abutment
Implant Ø: Universal
Platform: Universal
Code: 0074

DENTIUM

Model: NR Line
Implant Ø: 3.1
Platform: 3.2
Code: 0190

Model: NR Line
Implant Ø: 3.1
Platform: 3.6
Code: 0190

Model: NR Line
Implant Ø: 3.6
Platform: 3.6
Code: 0191

Model: NR Line
Implant Ø: 4.3
Platform: 4.3
Code: 0191

Model: NR Line
Implant Ø: 5
Platform: 5
Code: 0191

Model: NR Line
Implant Ø: 6
Platform: 6
Code: 0191

Model: Multi Unit NR Line
Implant Ø:
Platform: 5
Code: 0192

Model: SimpleLine II
Implant Ø: 3.8/4.3
Platform: 4.8
Code: 0074

Model: SimpleLine II
Implant Ø: 3.8/4.3
Platform: 4.8
Code: 0037

Model: SimpleLine II
Implant Ø: 4.3/4.8
Platform: 6.5
Code: 0096

Model: SuperLine/SuperLine II/
Implantium
Implant Ø: 3.4
Platform: 3.6
Code: 0030

Model: SuperLine/SuperLine II/
Implantium
Implant Ø: 3.8
Platform: 4
Code: 0030

Model: SuperLine/SuperLine II/
Implantium
Implant Ø: 4.3
Platform: 4.5
Code: 0030

Model: SuperLine/SuperLine II/
Implantium
Implant Ø: 4.8
Platform: 5
Code: 0030

Model: SuperLine/SuperLine II/
Implantium
Implant Ø: 4.8
Platform: 6
Code: 0030

Model: Multi Unit Superline and
Implantium
Implant Ø:
Platform: 4.5
Code: 0193

DIO IMPLANTS

Model: SM System
Implant Ø: 4.5/5/5.3
Platform: Regular/Wide
Code: 0033

Model: UF II Narrow
Implant Ø: 3/3.3
Platform: Narrow
Code: 0014

Model: UF II
Implant Ø: 3.8/4/4.5/5/5.5
Platform: Regular
Code: 0030

Model: External
Implant Ø: 3.3/3.8
Platform: Narrow 3.5
Code: 0023

Model: External
Implant Ø: 3.75/4/4.5
Platform: Regular 4.1
Code: 0024

Model: External
Implant Ø: 5/5.3/5.5/6
Platform: Wide 5.1
Code: 0061

Model: Internal OCTA
Implant Ø:
Platform: 4.8
Code: 0074

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0247

DITRON

Model: Ultimate Matrix
Implant Ø: 3.75/4.2
Platform: 3.75
Code: 0040

Model: Ultimate Matrix
Implant Ø: 5
Platform: 4
Code: 0040

Model: Ultimate Matrix
Implant Ø: 6
Platform: 4.6
Code: 0040

Model: MPI Matrix
Implant Ø: 3.5
Platform: 3.5
Code: 0040

Model: MPI Matrix
Implant Ø: 3.75
Platform: 3.75
Code: 0040

Model: MPI Matrix
Implant Ø: 4.2
Platform: 4.2
Code: 0040

Model: MPI Matrix
Implant Ø: 5
Platform: 5
Code: 0040

Model: MPI Matrix
Implant Ø: 6
Platform: 6
Code: 0040

DMI DENTAL SUPPLY

Model: DCI/DSI
Implant Ø: 3.3/3.5/3.75/4.2/5/6
Platform: 3.75
Code: 0040

DSP BIOMEDICAL

Model: Hexágono Externo
Implant Ø: 3.75/4/5//3.5/3.8/4.3
Platform: 4.1
Code: 0024

EASY IMPLANT

Model: Master C
Implant Ø: 3.5
Platform: 3.5 (Ocean)
Code: 0004

Model: Master C
Implant Ø: 4
Platform: 4 (Ocean)
Code: 0004

COMPATIBILITIES AVAILABLE

Model: Master C
Implant Ø: 4.5
Platform: 4.5 (Lilas)
Code: 0030

Model: Master C
Implant Ø: 5
Platform: 5 (Lilas)
Code: 0030

Model: Master S
Implant Ø: 3.3
Platform: 3.3 (Ocean)
Code: 0004

Model: Master S
Implant Ø: 3.75
Platform: 3.75 (Lilas)
Code: 0030

Model: Master S
Implant Ø: 4.25
Platform: 4.25 (Lilas)
Code: 0030

Model: Master S
Implant Ø: 4.75
Platform: 4.75 (Lilas)
Code: 0030

Model: Master L
Implant Ø: 3.3
Platform: 3.3 (Lilas)
Code: 0030

Model: Master L
Implant Ø: 3.75
Platform: 3.75 (Lilas)
Code: 0030

Model: Master L
Implant Ø: 4.25
Platform: 4.25 (Lilas)
Code: 0030

Model: Master L
Implant Ø: 4.75
Platform: 4.75 (Lilas)
Code: 0030

Model: MINI
Implant Ø: 3
Platform: 3
Code: 0176

Model: Hexcel-S
Implant Ø: 3.3
Platform: 3.3
Code: 0003

Model: Hexcel-S
Implant Ø: 3.75
Platform: 4.1
Code: 0024

Model: Hexcel-S
Implant Ø: 4.25
Platform: 4.1
Code: 0024

Model: Hexcel-S
Implant Ø: 4.75
Platform: 5
Code: 0058

Model: Multi Unit Conical
Abutment
Implant Ø:
Platform: Universal
Code: 0025

ECKERMANN

Model: All-Spiral
Implant Ø: 4
Platform: Regular
Code: 0069

Model: Duplo
Implant Ø: 4
Platform: Regular
Code: 0070

Model: Hexagon
Implant Ø: 3/3.5/4/4.5/5
Platform: 4.1
Code: 0024

Model: Winner
Implant Ø: 3/3.5/4
Platform: 3.5
Code: 0040_B

Model: Winner
Implant Ø: 3/3.5/4
Platform: 3.5
Code: 0042

Model: Winner
Implant Ø: 4/4.5/5
Platform: 4.5
Code: 0041_B

Model: Winner
Implant Ø: 4/4.5/5
Platform: 4.5
Code: 0043

ELITE MEDICA

Model: Conexión Externa
Implant Ø: 3.75
Platform: Narrow
Code: 0023

Model: Conexión Externa
Implant Ø: 4
Platform: Regular
Code: 0024

Model: Conexión Externa
Implant Ø: 5
Platform: Wide
Code: 0061

EUROTEKNIKA

Model: Naturactis
Implant Ø: 3.5
Platform: 3.4
Code: 0004

Model: Naturactis
Implant Ø: 4
Platform: 3.8
Code: 0004

Model: Naturactis
Implant Ø: 4.5
Platform: 4.3
Code: 0004

Model: Naturactis
Implant Ø: 5
Platform: 4.8
Code: 0004

Model: Uneva
Implant Ø: 3.6
Platform: 4.1
Code: 0024

Model: Uneva
Implant Ø: 4.1
Platform: 4.1
Code: 0024

Model: Uneva
(Platform Switching)
Implant Ø: 4.8
Platform: 4.1
Code: 0024

Model: Uneva
(Platform Switching)
Implant Ø: 6
Platform: 4.1
Code: 0024

Model: Natea
Implant Ø: 3.6/4.1/4.8
Platform: Narrow
Code: 0004

Model: Natea
Implant Ø: 3.6/4.1/4.8
Platform: Regular
Code: 0004

Model: Natea
Implant Ø: 6
Platform: Wide
Code: 0004

Model: Aesthetica
Implant Ø: 4.1
Platform: 4.8
Code: 0074

Model: Aesthetica
Implant Ø: 4.1
Platform: 4.8
Code: 0037

Model: Aesthetica
Implant Ø: 4.8
Platform: 6.5
Code: 0096

Model: Naturall
Implant Ø: 3.5
Platform: Narrow
Code: 0004

Model: Naturall
Implant Ø: 4/4.5
Platform: Regular
Code: 0004

Model: Naturall
Implant Ø: 5
Platform: Wide
Code: 0004

Model: Ibone E/Ibone S
Implant Ø: 3.8/4.3/4.8
Platform: 3.5
Code: 0004

Model: Ibone E/Ibone S
Implant Ø: 4.8/5.5/6.2
Platform: 4.3
Code: 0004

Model: Ibone G
Implant Ø: 4.8/5.5
Platform: RP
Code: 0037

Model: Ibone G
Implant Ø: 5.5/6.2
Platform: WP
Code: 0096

Model: Multi Unit Tetra
Implant Ø: Universal
Platform: Universal
Code: 0025

F&B IMPLANT (FIT & BRILLANT)

Model: FA Submerged Fixture
Implant Ø: 3.9
Platform: Narrow
Code: 0030

Model: FA Submerged Fixture
Implant Ø: 4.1/4.4
Platform: Regular
Code: 0030

Model: FA Submerged Fixture
Implant Ø: 4.8
Platform: Wide
Code: 0030

Model: FA Submerged Fixture
Implant Ø: 5.3/5.8/6.3/6.8
Platform: Ultra-Wide
Code: 0030

GALIMPLANT

Model: Conexión Externa
Implant Ø: 3.5/4
Platform: 4
Code: 0024

Model: Conexión Interna
Implant Ø: 3.5
Platform: 3.5
Code: 0004

Model: Conexión Interna
Implant Ø: 4
Platform: 4
Code: 0004

Model: Conexión Interna
Implant Ø: 5
Platform: 5
Code: 0004

Model: Pilar multi-posicion recto/
Pilar multi-posicion angulado
Implant Ø: Universal
Platform: Universal
Code: 0025

GC TECH

Model: AADVA Standard /
Tapered Implants
Implant Ø: 3.3
Platform: Narrow
Code: 0196

Model: AADVA Standard /
Tapered Implants
Implant Ø: 4
Platform: Regular
Code: 0197

Model: AADVA Standard /
Tapered Implants
Implant Ø: 5
Platform: Wide
Code: 0198

COMPATIBILITIES AVAILABLE

GMI (ILERIMPLANT)

Model: Phoenix
Implant Ø: 3.3/3.75/4
Platform: Standard 4.1
Code: 0024

Model: Phoenix
Implant Ø: 5
Platform: Wide 5.1
Code: 0061

Model: Frontier
Implant Ø: 3.3/3.75/4.25
Platform: RP 3.3
Code: 0042_B

Model: Frontier
Implant Ø: 4.75/5.75
Platform: WP 4.3
Code: 0041_B

Model: Universal
Implant Ø: Universal
Platform: PS-RP 4.8
Code: 0025

Model: Avantgard
Implant Ø: 3.75/4.25
Platform: RP
Code: 0043

GT MEDICAL

Model: Best Fit Octógono Interno
Implant Ø: 3.7/4.3/4.8
Platform: Regular
Code: 0074

Model: Best Fit Octógono Interno
Implant Ø: 3.7/4.3/4.8
Platform: Regular
Code: 0037

Model: Best Fit Hexágono Interno
Implant Ø: 3.7/4.1/4.3/4.8
Platform: Wide
Code: 0005

Model: Best Fit Hexágono Externo
Implant Ø: 3.5
Platform: Narrow
Code: 0023

Model: Best Fit Hexágono Externo
Implant Ø: 4.1
Platform: Regular
Code: 0024

Model: Best Fit Hexágono Externo
Implant Ø: 5
Platform: Wide
Code: 0061

HAHN IMPLANT (GLIDEWELL)

Model: Hahn Tapered Implant
Implant Ø: 3.5/4.3
Platform: 3.5/4.3
Code: 0021

Model: Hahn Tapered Implant
Implant Ø: 5
Platform: 5
Code: 0022

Model: Hahn Tapered Implant
Implant Ø: 7
Platform: 7
Code: 0124

Model: Multi-Unit Abutment system
Implant Ø: Universal
Platform: Universal
Code: 0025

Model: Hahn Tapered Implant
Implant Ø: 3
Platform: 3
Code: 0159

HIOSSSEN

Model: ETII SA / ETIII SA
Implant Ø: 3.5
Platform: Mini
Code: 0029

Model: ETII SA / ETIII SA / ETIV SA
Implant Ø: 4/4.5/5
Platform: Regular
Code: 0030

Model: ETIII BA
Implant Ø: 3.5
Platform: Mini
Code: 0029

Model: ETIII BA
Implant Ø: 4/4.5/5
Platform: Regular
Code: 0030

HI-TEC

Model: Tapered Self Thread
Implant Ø: 3.3/3.75
Platform: 3.5
Code: 0040

Model: Tapered Self Thread
Implant Ø: 4.2/5
Platform: 4.5
Code: 0041

Model: Logic Plus
Implant Ø: 3.5
Platform: 3.7
Code: 0040

Model: Logic Plus
Implant Ø: 4.3
Platform: 3.9
Code: 0040

IBS

Model: Magic FC
Implant Ø: 4/4.5/5/5.5/6/6.5
Platform: 3.8
Code: 0030

Model: N.R. Fix
Implant Ø: 3/3.5
Platform: 3.8
Code: 0030

IDO IMPLANTS

Model: I Do Implant
Implant Ø: 3.8/4/4.5/5/5.5/6/7
Platform: Universal
Code: 0030

IHDE DENTAL (IMBIODENT)

Model: Bone Level Plus
Implant Ø: 3.3
Platform: 3.3
Code: 0033

Model: Bone Level Plus
Implant Ø: 4.1
Platform: 4.1
Code: 0035

Model: Bone Level Plus
Implant Ø: 4.8
Platform: 4.8
Code: 0035

IMPLANT DIRECT

Model: RePlus / Replant / Reactive
Implant Ø: 3.5/3.7/4.2
Platform: 3.5
Code: 0026

Model: RePlus / Replant / Reactive
Implant Ø: 4.3/4.7
Platform: 4.3
Code: 0027

Model: RePlus / Replant / Reactive
Implant Ø: 5/5.7
Platform: 5
Code: 0028

Model: Legacy
Implant Ø: 3.7/4.2
Platform: 3.5
Code: 0040

Model: Legacy
Implant Ø: 4.7/5.2
Platform: 4.5
Code: 0041

Model: Swishplant / Swishplus
Implant Ø: 4.1/4.8
Platform: 4.8
Code: 0074

Model: Swishplant / Swishplus
Implant Ø: 4.1/4.8
Platform: 4.8
Code: 0037

Model: Swishplant / Swishplus
Implant Ø: 4.8/5.7
Platform: 6.5
Code: 0096

Model: SwishActive
Implant Ø: 3.3
Platform: 3
Code: 0021

Model: SwishActive
Implant Ø: 4.1/4.8
Platform: 3.4
Code: 0022

Model: Interactive
Implant Ø: 3.2/3.7
Platform: 3
Code: 0021

Model: Interactive
Implant Ø: 4.3/5
Platform: 3.4
Code: 0022

Model: Simply Iconic
Implant Ø: 3.2/3.7/4.2
Platform: Purple 3.0
Code: 0021

Model: Simply Iconic
Implant Ø: 4.7/5.2/5.7
Platform: Gold 3.4
Code: 0022

Model: Multi Unit Abutment
Implant Ø:
Platform: 5
Code: 00158

IMPLANT GENESIS

Model: Aktiv System
Implant Ø: 3.5/3.75/4.2/5
Platform: Standard
Code: 0040

IMPLANTSWISS

Model: Bone Level
Implant Ø: 3.3
Platform: 3.3
Code: 0004

Model: Bone Level
Implant Ø: 3.7
Platform: 3.7
Code: 0030

Model: Bone Level
Implant Ø: 4.3
Platform: 4.3
Code: 0030

Model: Bone Level
Implant Ø: 4.8
Platform: 4.8
Code: 0030

Model: Bone Level
Implant Ø: 5.5
Platform: 5.5
Code: 0030

Model: Multi Unit Abutment
Implant Ø: Universal
Platform: 4.8
Code: 0025

INTRA-LOCK

Model: Unihex
Implant Ø: 4
Platform: Regular
Code: 0024

Model: Unihex
Implant Ø: 4.75
Platform: Wide
Code: 0024

COMPATIBILITIES AVAILABLE

Model: IntraHex
Implant Ø: 3.75/4
Platform: 3.5
Code: 0040

Model: IntraHex
Implant Ø: 4.75
Platform: 4.5
Code: 0041

JDENTALCARE

Model: JDEvolution/
JDEvolution Plus
Implant Ø: 3.7
Platform: 3.7
Code: 0040

Model: JDEvolution/
JDEvolution Plus
Implant Ø: 4.3/5
Platform: 4
Code: 0040

Model: JDEvolution/
JDEvolution Plus
Implant Ø: 6
Platform: 5
Code: 0040

Model: JD Pterygo
Implant Ø: 4
Platform: 4
Code: 0040

Model: JD ICON/JD ICON F
Implant Ø: 3.9
Platform: 3.9
Code: 0022

Model: JD ICON/JD ICON F
Implant Ø: 4.3
Platform: 4
Code: 0022

Model: JD ICON/JD ICON F
Implant Ø: 5
Platform: 4.7
Code: 0022

Model: JD ICON Plus
Implant Ø: 3.7
Platform: 3.7
Code: 0015

Model: JD ICON Plus
Implant Ø: 4.3
Platform: 4
Code: 0015

Model: JD ICON Plus
Implant Ø: 5
Platform: 4.8
Code: 0015

Model: JD ICON Plus T
Implant Ø: 3.5
Platform: 3.5
Code: 0015

Model: JD ICON Plus T
Implant Ø: 4
Platform: 3.5
Code: 0015

Model: JD ICON Plus T
Implant Ø: 4.5
Platform: 3.5
Code: 0015

Model: Conical Abutment
Implant Ø:
Platform: Universal
Code: 0025

KEYSTONE

Model: Restore
Implant Ø: 3.75/4
Platform: RD 4.1
Code: 0024

Model: Restore
Implant Ø: 5/6
Platform: WD 5
Code: 0061

Model: Internal TiLobe
PrimaConnex
Implant Ø: 3.3/3.5
Platform: 3.5
Code: 0044

Model: Internal TiLobe
PrimaConnex
Implant Ø: 4/4.1
Platform: 4.1
Code: 0045

Model: Internal TiLobe
PrimaConnex
Implant Ø: 5
Platform: 5
Code: 0046

Model: Internal TiLobe Prima Plus
Implant Ø: 3.5
Platform: 3.5
Code: 0044

Model: Internal TiLobe Prima Plus
Implant Ø: 4.1
Platform: 4.1
Code: 0045

Model: Internal TiLobe Prima Plus
Implant Ø: 5/6
Platform: 5
Code: 0046

KLOCKNER

Model: Essential Cone
Implant Ø: 3.5/4/4.5
Platform: 4.5
Code: 0054

Model: Essential Cone
Implant Ø: 4.8
Platform: 6
Code: 0071

Model: Essential Cone Pilar 25°
Implant Ø: 3.5/4/4.5
Platform: 4.5
Code: 0054

Model: Essential Cone
Octacone 12°
Implant Ø: 3.5/4/4.5
Platform: 4.5
Code: 0054

Model: KL
Implant Ø: 3.3
Platform: Narrow
Code: 0023

Model: KL
Implant Ø: 3.7/4.2
Platform: Regular
Code: 0024

Model: KL
Implant Ø: 4.7
Platform: Wide
Code: 0061

Model: Vega
Implant Ø: 3.5
Platform: NV
Code: 0082

Model: Vega
Implant Ø: 4/4.5
Platform: RV
Code: 0083

Model: Vega+
Implant Ø: 3.6
Platform: NV
Code: 0082

Model: Vega+
Implant Ø: 4.1/4.6
Platform: RV
Code: 0083

Model: Multi Unit Permanent
Implant Ø: 4.2
Platform: Universal
Code: 0173

KUWOTECH

Model: KISPLANT
Implant Ø: Cono Morse +
Hexágono Interno
Platform: 3.5
Code: 0030

Model: KISPLANT
Implant Ø: Cono Morse +
Hexágono Interno
Platform: 4/4.5/5/5.5/6
Code: 0030

LASAK

Model: Bioniq
Implant Ø: 2.9
Platform: QR (Amarillo)
Code: 0166

Model: Bioniq
Implant Ø: 3.5
Platform: QR (Azul)
Code: 0167

Model: Bioniq
Implant Ø: 4
Platform: QR (Azul)
Code: 0167

Model: Bioniq
Implant Ø: 5
Platform: QR (Azul)
Code: 0167

LEADER

Model: Tixos Internal Hex
Implant Ø: 3.3
Platform: 3.5
Code: 0040

Model: Tixos Internal Hex
Implant Ø: 3.75
Platform: 4
Code: 0040

Model: Tixos External Hex
Implant Ø: 3.3/3.75
Platform: 4.1
Code: 0024

Model: Tixos External Hex
Implant Ø: 5
Platform: 5
Code: 0058

MEDENTIKA

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0025

MEDENTIS

Model: Premium/Active Master
Implant Ø: 3.3
Platform: 3.3 (Pink)
Code: 0249

Model: Premium/Active Master
Implant Ø: 3.75
Platform: 3.75 (Red)
Code: 0125

Model: Premium/Active Master
Implant Ø: 4.1
Platform: 4.1 (Green)
Code: 0125

Model: Premium/Active Master
Implant Ø: 4.8
Platform: 4.8 (Blue)
Code: 0125

MEGAGEN

Model: AnyRidge
Implant Ø: 3.5
Platform: Small
Code: 0015

Model: AnyRidge
Implant Ø: 4/4.5
Platform: Regular
Code: 0015

Model: AnyRidge
Implant Ø: 5/5.5
Platform: Wide
Code: 0015

Model: AnyOne Internal
Implant Ø: 3.5/4/4.5/5/6/7
Platform: General
Code: 0030

Model: AnyOne External
Implant Ø: 3.5
Platform: Small 3.5
Code: 0023

COMPATIBILITIES AVAILABLE

Model: AnyOne External
Implant Ø: 4
Platform: Regular 4.1
Code: 0024

Model: AnyOne External
Implant Ø: 4.5
Platform: Regular 4.5
Code: 0024

Model: AnyOne External
Implant Ø: 5
Platform: Regular 5
Code: 0058

Model: AnyOne External
Implant Ø: 6
Platform: SuperWide 5.5
Code: 0058

Model: Cone Abutment
Implant Ø: Universal
Platform: 3.8
Code: 0128

Model: Cone Abutment
Implant Ø: Universal
Platform: 4.8
Code: 0074

Model: Mini Narrow Ridge
Implant Ø: 3/3.4
Platform: Mini
Code: 0014

Model: ExFeel
Implant Ø: 3.5
Platform: Small
Code: 0037

Model: ExFeel
Implant Ø: 4.1
Platform: Regular
Code: 0037

Model: ExFeel
Implant Ø: 4.8/5
Platform: Wide
Code: 0037

Model: Multi Unit N Type
Implant Ø: Multi Unit N Type
Platform: Universal
Code: 0025

Model: Multi Unit S Type
Implant Ø: Multi Unit S Type
Platform: Universal
Code: 0264

Model: Blue Diamond/BD Cuff
Implant Ø: 3.3/3.7
Platform: NC
Code: 0222

Model: Blue Diamond/BD Cuff
Implant Ø: 4.1/4.4/4.8/5.3/5.8/6.3
Platform: RC
Code: 0223

MICRODENT

Model: Universal
Implant Ø: 2.8/3.25
Platform: 3.5
Code: 0003

Model: Universal
Implant Ø: 3.3/3.5/3.75/4
Platform: 4.1
Code: 0024

Model: Universal
Implant Ø: 4.2/5
Platform: 5.1
Code: 0058

Model: System
Implant Ø: 2.8/3.25
Platform: 3.5
Code: 0003

Model: Ektos
Implant Ø: 3.7/4.2
Platform: 3.5
Code: 0040_B

MIS

Model: Lance
Implant Ø: 3.75/4.2
Platform: Standard
Code: 0024

Model: Lance
Implant Ø: 5
Platform: Wide
Code: 0058

Model: Multi-Unit
Implant Ø:
Platform: General
Code: 0020

Model: Seven
Implant Ø: 3.3
Platform: Narrow
Code: 0019

Model: Seven
Implant Ø: 3.75/4.2
Platform: Standard
Code: 0040

Model: Seven
Implant Ø: 5/6
Platform: Wide
Code: 0041

Model: M4
Implant Ø: 3.3
Platform: Narrow
Code: 0019

Model: M4
Implant Ø: 3.75/4.2
Platform: Standard
Code: 0040

Model: M4
Implant Ø: 5/6
Platform: Wide
Code: 0041

Model: C1
Implant Ø: 3.3
Platform: Narrow
Code: 0016

Model: C1
Implant Ø: 3.75/4.2
Platform: Standard
Code: 0017

Model: C1
Implant Ø: 5
Platform: Wide
Code: 0018

Model: V3
Implant Ø: 3.9/4.3/5
Platform: Standard
Code: 0017

MONOIMPLANT

Model: MonoImplant Multi Unit
Implant Ø: 3/3.7/4.1
Platform: 4.8
Code: 0025

MOZO-GRAU (TICARE)

Model: MG Osseous
Implant Ø: 3.3
Platform: 3.4 Mini
Code: 0003

Model: MG Osseous
Implant Ø: 3.4/3.75/4.25
Platform: 4.1 Standard
Code: 0024

Model: MG Osseous
Implant Ø: 5
Platform: 5 Maxi
Code: 0061

Model: MG Inhex
Implant Ø: 3.3
Platform: 2.3 Mini
Code: 0109

Model: MG Inhex
Implant Ø: 3.75/4.25
Platform: 2.8 Standard
Code: 0004

Model: MG Inhex
Implant Ø: 5
Platform: 3.8 Maxi
Code: 0005

MPI

Model: Conexión Externa HE
Privilege
Implant Ø: 3.3
Platform: 3.5
Code: 0009

Model: Conexión Externa HE
Privilege
Implant Ø: 3.3/4
Platform: 4.1
Code: 0024

Model: Conexión Externa HE
Privilege
Implant Ø: 5
Platform: 5
Code: 0058

Model: Privilege CM
Implant Ø: 3.5/4
Platform: Regular
Code: 0004

Model: Privilege CM
Implant Ø: 5
Platform: Wide
Code: 0005

Model: Excellence CM
Implant Ø: 3.5/4
Platform: Regular
Code: 0004

Model: Excellence CM
Implant Ø: 5
Platform: Wide
Code: 0005

NEOBIOTECH

Model: EB External System
Implant Ø: 3.5
Platform: Narrow
Code: 0023

Model: IS Implant System
Implant Ø: 3.2
Platform: S-Narrow
Code: 0029

Model: IS Implant System
Implant Ø: 3.5
Platform: Narrow
Code: 0030

Model: IS Implant System
Implant Ø: 4
Platform: Regular
Code: 0030

Model: IS Implant System
Implant Ø: 4.5
Platform: Regular
Code: 0030

Model: IS Implant System
Implant Ø: 5
Platform: Wide
Code: 0030

Model: IS Implant System
Implant Ø:
Platform: 4.8
Code: 0025

NEODENT

Model: Helix GM/Drive GM/
Titamax GM
Implant Ø: 3.5/3.75/4/4.3/5/6
Platform: Regular
Code: 0186

Model: Smart HE
Implant Ø: 3.75/4
Platform: 4.1
Code: 0024

Model: Helix HE
Implant Ø: 3.75/4/4.3
Platform: 4.1
Code: 0024

Model: Mini Pilar CM /
Mini Pilar Angulado CM
Implant Ø: Mini Pilar CM /
Mini Pilar Angulado CM
Platform: Universal
Code: 0025

COMPATIBILITIES AVAILABLE

NEOSS

Model: ProActive Straight/
Tapered/Edge
Implant Ø: 3.5 Green
Platform: Standard
Code: 0047

Model: ProActive Straight/
Tapered/Edge
Implant Ø: 4 Yellow
Platform: Standard
Code: 0047

Model: ProActive Straight/
Tapered/Edge
Implant Ø: 4.5 Blue
Platform: Standard
Code: 0048

Model: ProActive Straight/
Tapered/Edge
Implant Ø: 5 Peach
Platform: Standard
Code: 0048

Model: ProActive Straight/
Tapered/Edge
Implant Ø: 5.5 Lilac
Platform: Standard
Code: 0048

Model: ProActive Wide
Implant Ø: 6
Platform: Standard
Code: 0048

Model: ProActive Sinus
Implant Ø: 6.5
Platform: Standard
Code: 0048

Model: Short Implant
Implant Ø: 3.5/4/4.5
Platform: Standard
Code: 0047

Model: Short Implant
Implant Ø: 5/5.5/6/6.5
Platform: Standard
Code: 0048

NOBEL BIO CARE

Model: Branemark
Implant Ø: 3.3
Platform: Narrow
Code: 0023

Model: Branemark
Implant Ø: 3.75/4
Platform: Regular
Code: 0024

Model: Branemark
Implant Ø: 5/6
Platform: Wide
Code: 0061

Model: Replace
Implant Ø: 3.5
Platform: Narrow
Code: 0026

Model: Replace
Implant Ø: 4.3
Platform: Regular
Code: 0027

Model: Replace
Implant Ø: 5
Platform: Wide
Code: 0028

Model: Replace
Implant Ø: 6
Platform: Platform 6
Code: 0129

Model: Active/Replace
Conical Connection
Implant Ø: 3
Platform: Mini 3.0
Code: 0159

Model: Active/Replace
Conical Connection
Implant Ø: 3.5
Platform: Narrow
Code: 0021

Model: Active/Replace
Conical Connection
Implant Ø: 4.3/5
Platform: Regular
Code: 0022

Model: Active/Replace
Conical Connection
Implant Ø: 5.5
Platform: Wide
Code: 0124

Model: NobelSpeedy
Implant Ø: 3.3
Platform: Narrow
Code: 0023

Model: NobelSpeedy
Implant Ø: 4/5
Platform: Regular
Code: 0024

Model: NobelSpeedy
Implant Ø: 5/6
Platform: Wide
Code: 0061

Model: NobelParallel
Implant Ø: 3.75
Platform: Narrow
Code: 0021

Model: NobelParallel
Implant Ø: 4.3/5
Platform: Regular
Code: 0022

Model: NobelParallel
Implant Ø: 5.5
Platform: Wide
Code: 0124

NORIS MEDICAL

Model: Tuff
Implant Ø: 3.3/3.75/4.2/5/6
Platform: 3.75
Code: 0040

Model: Tuff TT
Implant Ø: 3.3/3.75/4.2/5/6
Platform: 3.75
Code: 0040

Model: Onix
Implant Ø: 3.3/3.75/4.2/5/6
Platform: 3.75
Code: 0040

Model: Cortical
Implant Ø: 4.0/5/6
Platform: 3.75
Code: 0040

Model: PteryCore
Implant Ø: 4.2
Platform: 3.75
Code: 0040

Model: PteryFit
Implant Ø: 4.2
Platform: 3.75
Code: 0040

NORMON

Model: Normoimplant HE
Implant Ø: 3.25/3.75/4.25/4.75
Platform: 4.1
Code: 0024

Model: Normoimplant HI
Implant Ø: 3.75/4.25/4.75
Platform: 3.5
Code: 0040_B

NOVA IMPLANTS

Model: PSI/PCI
Implant Ø: 3.3/3.75/4.2/5/6
Platform: 3.75
Code: 0040_B

OSSTEM IMPLANT

Model: TS
Implant Ø: 3
Platform: Mini 3.0
Code: 0029

Model: TS
Implant Ø: 3.5
Platform: Mini 3.5
Code: 0029

Model: TS
Implant Ø: 4/4.5/5/6/7
Platform: Regular
Code: 0030

Model: US
Implant Ø: 3.3/3.5
Platform: Mini 3.5
Code: 0023

Model: US
Implant Ø: 3.75/4/4.5
Platform: Regular 4.1
Code: 0024

Model: US
Implant Ø: 5/5.5
Platform: Wide 5.1
Code: 0061

Model: US
Implant Ø: 5/5.5
Platform: Wide PS 5
Code: 0058

Model: Multi Unit Abutment/
Esthetic Low abutment
Implant Ø: Universal
Platform: Mini/Regular
Code: 0025

OSTEOPLUS

Model: She
Implant Ø: 3.45
Platform: 3.45
Code: 0009

Model: She
Implant Ø: 3.75 / 4
Platform: 4
Code: 0024

Model: Shi
Implant Ø: 3.3 / 3.75 / 4.2
Platform: 3.5
Code: 0040

OXY

Model: K1 Line
Implant Ø: 3.5/4/4.5/5
Platform: Regular
Code: 0015

Model: K1 Line
Implant Ø: 5.5/6/6.5
Platform: Wide
Code: 0015

Model: PSK Line
Implant Ø: 3.75/4.5/5
Platform: Regular
Code: 0015

Model: MD Line KONE
Implant Ø: 3.75/4.25/5
Platform: Regular
Code: 0015

Model: MD Line Ext
Implant Ø: 3.75/4.25
Platform: Standard
Code: 0024

Model: Fixo
Implant Ø: Universal
Platform: 4.8
Code: 0242

PALTOP

Model: Advanced classic
Implant Ø: 3.25
Platform: Narrow (blue)
Code: 0229

Model: Advanced classic
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0040_B

Model: Advanced classic
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0042

Model: Advanced classic
Implant Ø: 6
Platform: Wide (purple)
Code: 0041

Model: Advanced *
Implant Ø: 3.25
Platform: Narrow (blue)
Code: 0229

COMPATIBILITIES AVAILABLE

Model: Advanced +
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0040_B

Model: Advanced +
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0042

Model: Advanced +
Implant Ø: 6
Platform: Wide (purple)
Code: 0041

Model: Dynamic
Implant Ø: 3.25
Platform: Narrow (blue)
Code: 0029

Model: Dynamic
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0040_B

Model: Dynamic
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0042

Model: Dynamic
Implant Ø: 6
Platform: Wide (purple)
Code: 0041

Model: PAI
Implant Ø: 3.25
Platform: Narrow (blue)
Code: 0029

Model: PAI
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0040_B

Model: PAI
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0042

Model: PAI
Implant Ø: 6
Platform: Wide (Purple)
Code: 0041

Model: DIVA/ACTIVE
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0040_B

Model: DIVA/ACTIVE
Implant Ø: 3.75/4.2/5
Platform: Standard
Code: 0042

Model: Conical Active
Implant Ø: 3.25/3.75/4.2/5
Platform: Standard
Code: 0029

Model: Universal Multi-Unit
Implant Ø: Multi-Unit
Platform: Universal
Code: 0181

PHIBO

Model: TSH/BNT Serie 3
Implant Ø: 3.6
Platform: 4
Code: 0024

Model: TSH/BNT Serie 4
Implant Ø: 4.2
Platform: 4
Code: 0024

Model: TSH/BNT Serie 5
Implant Ø: 4.8
Platform: 5
Code: 0060

POINT IMPLANT

Model: SM II System UV Implant
Implant Ø: 4/4.5/5
Platform: Regular
Code: 0030

PROCLINIC

Model: Cilindrico Externo
Implant Ø: 3.3
Platform: 3.3 Mini
Code: 0009

Model: Cilindrico Externo/
Cónico Externo
Implant Ø: 3.75/4.25//3.5/4
Platform: 4.1 Estandar
Code: 0024

Model: Cilindrico Externo/
Cónico Externo
Implant Ø: 5
Platform: 5 Maxi
Code: 0058

Model: Cilindrico Interno/
Cónico Interno
Implant Ø:
3.3/3.75/4.25/5/3.5/4/5
Platform: 3.5
Code: 0040

Model: SP Octa
Implant Ø: 3.3/4.1/4.8
Platform: 4.8
Code: 0074

Model: SP Octa
Implant Ø: 3.3/4.1/4.8
Platform: 4.8
Code: 0037

Model: SP Octa
Implant Ø: 4.8
Platform: 6.5
Code: 0096

Model: Aqua CM
Implant Ø: 3.5/4/5
Platform: 2.82
Code: 0004

PROTEG IMPLANTS

Model: PR01
Implant Ø: 3.3/3.75/4.2/5/6
Platform: Narrow
Code: 0265

RADHEX

Model: PHE
Implant Ø: 3.5
Platform: 3.5
Code: 0023

Model: PHE
Implant Ø: 4/4.5/5
Platform: 4.1
Code: 0024

Model: PHI
Implant Ø: 3.75
Platform: 3.5
Code: 0040_B

Model: PHI
Implant Ø: 4.5/5
Platform: 4.5
Code: 0041_B

REFLECT

Model: Rapid
Implant Ø: 3.0
Platform: 3.0
Code: 0159

Model: Rapid
Implant Ø: 3.5
Platform: NP
Code: 0021

RITTER

Model: SB/LA/QSI
Implant Ø: 3.75/4.2/5/6
Platform: Standard
Code: 0040_B

Model: Rapid
Implant Ø: 4.3/5/5.5
Platform: RP
Code: 0022

Model: Recover
Implant Ø: 3.5
Platform: NP
Code: 0021

Model: Recover
Implant Ø: 4.3/5
Platform: RP
Code: 0022

Model: Aspire
Implant Ø: 3.5/4
Platform: Aqual(Estrecha)
Code: 0004

Model: Aspire
Implant Ø: 5
Platform: Lilac (Ancha)
Code: 0005

Model: Tapered Screw
Implant Ø: 3.5
Platform: 3.5
Code: 0040_B

Model: Tapered Screw
Implant Ø: 4.1/4.7
Platform: 4.5
Code: 0041

ROOT

Model: R
Implant Ø: 3.5/3.8/4.2/4.8/5.5
Platform: Universal
Code: 0236

SEWON MEDIX

Model: IH2 SLA SYSTEM
Implant Ø: 3.5
Platform: Mini
Code: 0029

Model: IH2 SLA SYSTEM
Implant Ø: 3.5/4/4.5/5
Platform: Regular
Code: 0030

Model: IH2 RBM SYSTEM
Implant Ø: 3.5
Platform: Mini
Code: 0029

Model: IH2 RBM SYSTEM
Implant Ø: 3.5/4/4.5/5
Platform: Regular
Code: 0030

Model: IH SYSTEM
Implant Ø: 3.5/4/4.5/5
Platform: Universal
Code: 0025

SIC INVENT

Model: HEXAGONAL SYSTEM
SiCace
Implant Ø: 3.4/4
Platform: 3.3
Code: 0170

Model: HEXAGONAL SYSTEM
SiCace
Implant Ø: 4.5/5
Platform: 4.2
Code: 0171

Model: HEXAGONAL SYSTEM
SiCMax
Implant Ø: 3.4/4.2
Platform: 3.3
Code: 0170

Model: HEXAGONAL SYSTEM
SiCMax
Implant Ø: 4.7/5.2
Platform: 4.2
Code: 0171

Model: HEXAGONAL SYSTEM
SiCtapered
Implant Ø: 3.4/4.2
Platform: 3.3
Code: 0170

Model: HEXAGONAL SYSTEM
SiCtapered
Implant Ø: 4.7/5.2
Platform: 4.2
Code: 0171

SIGNO VINCES

Model: Duo
Implant Ø: 4.6
Platform: 4.1
Code: 0024

COMPATIBILITIES AVAILABLE

Model: Integra
Implant Ø: 3.75/4
Platform: 4.1
Code: 0024

Model: Compact
Implant Ø: 4.5
Platform: CM3.8
Code: 0004

Model: Duocon
Implant Ø: 3.8
Platform: CM3.8
Code: 0004

Model: Duocon
Implant Ø: 4.6/5.5
Platform: CM4.6
Code: 0005

Model: Infra
Implant Ø: 3.3/3.8/4.6
Platform: CM
Code: 0004

SIN IMPLANTS

Model: HE
Implant Ø: 3.75
Platform: 4.1
Code: 0024

Model: HI SW
Implant Ø: 3.8
Platform: 3.8
Code: 0039

Model: Tryon
Implant Ø: 3.25/3.75/4
Platform: 4.1
Code: 0024

Model: Tryon CO
Implant Ø: 4
Platform: 4.1
Code: 0024

Model: Revolution
Implant Ø: 3.25/3.75/4
Platform: 4.1
Code: 0024

Model: Stylus
Implant Ø: 4
Platform: 4.1
Code: 0024

Model: Epikut/Epikut Plus
Implant Ø: 4.5
Platform: 4.5
Code: 0024

Model: Epikut/Epikut Plus
Implant Ø: 5
Platform: 5
Code: 0058

Model: Epikut/Epikut Plus
Implant Ø: 3.5/3.8/4.5/5
Platform: Universal
Code: 0020

Model: Mini Pilar
Implant Ø: Universal
Platform: Universal
Code: 0025

SOUTHERN IMPLANTS

Model: Tri-Nex
Implant Ø: 3.5
Platform: 3.5
Code: 0026

Model: Tri-Nex
Implant Ø: 4.3
Platform: 4.3
Code: 0027

Model: Tri-Nex
Implant Ø: 5
Platform: 5
Code: 0028

Model: Tri-Nex
Implant Ø: 6
Platform: 6
Code: 0129

Model: IT Connection
Implant Ø: 3.3/4/4.1/4.9/5
Platform: 4.8
Code: 0037

Model: IT Connection
Implant Ø: 4.9/5/6
Platform: 6.5
Code: 0096

Model: External Hex
Implant Ø: 3.25
Platform: 3.4
Code: 0003

Model: External Hex
Implant Ø: 3.75/4
Platform: 4.1
Code: 0024

Model: External Hex
Implant Ø: 4.7/5
Platform: 5
Code: 0058

Model: External Hex
Implant Ø: 5.7/6
Platform: 6
Code: 0058

Model: Deep Conical
Implant Ø: 3
Platform: 2.45
Code: 0109

Model: Deep Conical
Implant Ø: 3.5/4
Platform: 2.95/3.1
Code: 0004

Model: Deep Conical
Implant Ø: 5
Platform: 4.1
Code: 0005

Model: Internal Hex
Implant Ø: 3.75/4.2/5
Platform: Universal
Code: 0040

Model: Provata
Implant Ø: 4/5/6
Platform: Standard
Code: 0040

Model: Compact Conical
Implant Ø: 4.8
Platform: 4.8
Code: 0025

STERI-OSS

Model: Hex-Loc
Implant Ø: 3.25
Platform: 3.3
Code: 0023

STERNGOLD

Model: STERN EX
Implant Ø: 3.75/4/5
Platform: 4.1
Code: 0024

STRAUMANN

Model: Tissue Level
Implant Ø: 3.3
Platform: 3.5
Code: 0160

Model: Tissue Level
Implant Ø: 3.3/4.1/4.8
Platform: Regular 4.8
Code: 0037

Model: Tissue Level
Implant Ø: 4.8
Platform: Wide 6.5
Code: 0096

Model: Synocta
Implant Ø:
Platform: Regular 4.8
Code: 0074

Model: Synocta
Implant Ø:
Platform: Wide 6.5
Code: 0137

Model: Bone Level Tapered SC
Implant Ø: 2.9
Platform: SC- 2.9
Code: 0135

Model: Bone Level
Implant Ø: 3.3
Platform: NC- 3.3
Code: 0033

Model: Bone Level
Implant Ø: 4.1
Platform: RC-4.1
Code: 0035

Model: Bone Level
Implant Ø: 4.8
Platform: RC-4.8
Code: 0035

Model: Screw-Retained
Implant Ø: NC/RC Ø4.6
Platform: Universal
Code: 0101

Model: BLX
Implant Ø: 3.5/3.75/4/4.5
Platform: RB (Regular Base)
Code: 0207

Model: BLX
Implant Ø: 5/5.5/6.5
Platform: WB (Wide Base)
Code: 0208

Model: BLC
Implant Ø: 3.3/3.75
Platform: RB
Code: 0207

Model: BLC
Implant Ø: 4.5/5.5/6.5
Platform: WB
Code: 0208

Model: TLX / TLX S
Implant Ø: 3.75/4.5
Platform: NT
Code: 0260

Model: TLX / TLX S
Implant Ø: 3.75/4.5
Platform: RT
Code: 0261

Model: TLX / TLX S
Implant Ø: 5.5/6.5
Platform: WT
Code: 0262

Model: TLC
Implant Ø: 3.3
Platform: NT
Code: 0260

Model: TLC
Implant Ø: 3.3/3.75/4.5
Platform: RT
Code: 0261

Model: TLC
Implant Ø: 4.5/5.5/6.5
Platform: WT
Code: 0262

SURCAM DENTAL

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0181

SYBRON IMPLANT SOLUTIONS

Model: Endopore (Innova)
Implant Ø: 4.1
Platform: 4.1
Code: 0024

SYSTHEX

Model: Classic-ci / Estetic-ci
Implant Ø: 3.5/3.75/4
Platform: 4.1
Code: 0024

TBR

Model: Hex-Conic
Implant Ø: 3.5
Platform: Narrow
Code: 0023

Model: Hex-Conic
Implant Ø: 5
Platform: Wide
Code: 0058

COMPATIBILITIES AVAILABLE

Model: Connect / Infinity
Implant Ø: 3.5
Platform: 3.5
Code: 0266

Model: Connect / Infinity
Implant Ø: 4
Platform: 4
Code: 0267

Model: Connect / Infinity
Implant Ø: 5
Platform: 5
Code: 0268

Model: Baby 8
Implant Ø: 4
Platform: 4
Code: 0267

Model: Baby 8
Implant Ø: 5
Platform: 5
Code: 0268

TITANIUM-FIX

Model: b-fix
Implant Ø: 3.5/4
Platform: Regular
Code: 0004

Model: b-fix
Implant Ø: 4.5/5
Platform: Larga
Code: 0005

TREE-OSS

Model: Simple
Implant Ø: 3.3/3.75/5
Platform: 3.75 Amarillo
Code: 0040

Model: Rapid/Anatomic
Implant Ø: 3.3
Platform: 3.5 Rosa
Code: 0023

Model: Rapid/Anatomic
Implant Ø: 3.75/4
Platform: 4.1 Amarillo
Code: 0024

Model: Rapid/Anatomic
Implant Ø: 5
Platform: 5.1 Azul
Code: 0061

Model: Anatomic/HS
Implant Ø: 3.5
Platform: 3.5 Rosa
Code: 0026

Model: Anatomic/HS
Implant Ø: 4.3
Platform: 4.3 Amarillo
Code: 0027

Model: Anatomic/HS
Implant Ø: 5
Platform: 5 Azul
Code: 0028

Model: Multi Unit
Implant Ø:
Platform: Universal
Code: 0025

TRI DENTAL IMPLANTS

Model: TRI-Vent
Implant Ø: 3.75/4.1/4.7
Platform: 3.5
Code: 0040

Model: TRI-Vent
Implant Ø: 3.75/4.1/4.7
Platform: 3.5
Code: 0042

TRINON

Model: Q2
Implant Ø: 3.5/3.75/4.5
Platform: 4
Code: 0024

Model: OK
Implant Ø: 4
Platform: 4.8
Code: 0074

Model: OK
Implant Ø: 4
Platform: 4.8
Code: 0037

UFIT

Model: Gtz
Implant Ø: 3.5
Platform: Mini
Code: 0004

Model: Gtz
Implant Ø: 4/4.5
Platform: Regular
Code: 0005

Model: Gtz
Implant Ø: 5
Platform: Wide
Code: 0005

Model: Gtz
Implant Ø: 5.5/6/6.5/7
Platform: Ultra-wide
Code: 0005

Model: Ntz
Implant Ø: 3.5
Platform: Mini
Code: 0004

Model: Ntz
Implant Ø: 4/4.5
Platform: Regular
Code: 0005

Model: Ntz
Implant Ø: 5
Platform: Wide
Code: 0005

Model: Ntz
Implant Ø: 5.5/6/6.5/7
Platform: Ultra-wide
Code: 0005

VULKAN IMPLANTS

Model: IN-Hex
Implant Ø: 3.3/3.75/4.2/5
Platform: 3.75
Code: 0040

WARANTEC (ONEPLANT)

Model: IU Implant System
Implant Ø: 3.3/3.6
Platform: Mini
Code: 0004

Model: UT Implant System
Implant Ø: 3.6
Platform: Mini
Code: 0004

WIN

Model: WIN
Implant Ø: 3.30/3.75/4.25/5
Platform: 3.75
Code: 0040_B

Model: WIN
Implant Ø: Universal
Platform: Universal
Code: 0025

XIVE

Model: Xive
Implant Ø: 3
Platform: 3
Code: 0084

Model: Xive
Implant Ø: 3.4
Platform: 3.4
Code: 0038

Model: Xive
Implant Ø: 3.8
Platform: 3.8
Code: 0039

Model: Xive
Implant Ø: 4.5
Platform: 4.5
Code: 0085

Model: Xive
Implant Ø: 5.5
Platform: 5.5
Code: 0086

YES IMPLANT

Model: S-SYSTEM
Implant Ø: 3.3/3.5
Platform: Narrow
Code: 0030

Model: S-SYSTEM
Implant Ø: 4/4.5
Platform: Regular
Code: 0030

Model: S-SYSTEM
Implant Ø: 5/5.5
Platform: Wide
Code: 0030

ZIACOM

Model: OEX
Implant Ø: 3.75/4.25
Platform: RP 4.1
Code: 0024

Model: ZINIC
Implant Ø: 3.7/4/4.3
Platform: RP 3.5
Code: 0040_B

Model: ZINIC Shorty
Implant Ø: 4.75
Platform: RP 3.5
Code: 0040_B

Model: GALAXY
Implant Ø: 3.4/3.7/4
Platform: RP 2.85
Code: 0004

ZIMMER

Model: Screw-Vent
Implant Ø: 3.7/4.1
Platform: 3.5
Code: 0040

Model: Screw-Vent
Implant Ø: 3.7/4.1
Platform: 3.5
Code: 0042

Model: Screw-Vent
Implant Ø: 4.7
Platform: 4.5
Code: 0041

Model: Screw-Vent
Implant Ø: 4.7
Platform: 4.5
Code: 0043

Model: Screw-Vent
Implant Ø: 6
Platform: 5.7
Code: 0080

Model: TSX Implant
Implant Ø: 3.1
Platform: 2.9
Code: 0178

Model: TSX Implant
Implant Ø: 3.7/4.1/4.7
Platform: 3.5
Code: 0040

Model: TSX Implant
Implant Ø: 3.7/4.1/4.7
Platform: 3.5
Code: 0042

Model: TSX Implant
Implant Ø: 5.4/6
Platform: 4.5
Code: 0041

Model: TSX Implant
Implant Ø: 5.4/6
Platform: 4.5
Code: 0043


COMPATIBILITIES AVAILABLE

Model: Swiss-Plus
Implant Ø: 3.7/4.1/4.8
Platform: 4.8
Code: 0074

Model: Swiss-Plus
Implant Ø: 3.7/4.1/4.8
Platform: 4.8
Code: 0037

Model: Eztetic
Implant Ø: 3.1
Platform: 2.9
Code: 0178

Model: Tapered Abutment
Implant Ø: Universal
Platform: Universal
Code: 0205



This printed version is probably
outdated and may not include all
products or compatibilities.

**Access and download the latest
version of this catalogue**

www.dynamicabutment.com

Some of the compatibilities listed may not appear in
the catalogue, as some of the associated products
are being developed and/or manufactured. If you
have any queries, please contact us.

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
0,3 mm				1,2 mm				2 mm			
R	31.322.001.01-2	43°	29°	31.322.001.02-2	25°	22°	25°	31.322.001.03-2	25°	31.322.001.04-2	
NR	31.322.001.01-2			31.312.001.02-2				31.312.001.03-2		31.312.001.04-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s
0,3 mm		CH-5mm	CH-7mm	CH-9mm	2 mm		CH-5mm	CH-7mm	CH-9mm
R	31.322.001.21-2	25°	20°	10°		25°	20°	10°	
NR	31.312.001.21-2				31.312.001.23-2				

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
		50.312.001.01-2	43.621.410.01-2	34.612.001.01-2	33.390.754.01-2	3	25°	23.412.001.01-2	54.315.001.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.103.01-2	10	50.312.001.04-2 (IG=3mm)	43.624.410.01-2		33.490.754.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
52.412.103.01-2	12				33.690.754.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
DYNAMIC SCREWS					STRAIGHT SCREWS					ANALOG	
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)			STRAIGHT SCREW	SCREWDRIVER Hex. 1.20		49.416.000.04-2 (13 mm)			
41.316.084.01-2		43.618.201.01-2	18		40.316.003.01-2	43.601.103.02-2		22.612.001.01-2	LAB SCANBODY	49.416.000.05-2 (13 mm)	E
		43.624.201.01-2	24				30.412.001.01-2	SCREWDRIVER		43.625.103.01-2	
		43.632.201.01-2	32								

MULTI-UNIT

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1 mm		2 mm		3 mm		4 mm	
R	42.302.001.01-2	42.302.001.02-2	42.302.001.03-2	42.302.001.04-2			

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c GINGIVAL HEIGHT α_s

1 mm 2 mm

R	31.323.002.01-2	45°	29°	31.323.002.02-2	25°
NR	31.313.002.01-2			31.313.002.02-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0.3 mm CH-5mm CH-7mm CH-9mm

R	31.323.002.21-2	25°	20°	10°
NR	31.313.002.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.408.101.01-2	8	50.313.002.01-2	43.621.410.01-2 43.624.410.01-2	34.613.002.01-2
52.410.101.01-2	10			
52.412.101.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.805.01-2	3	30°
33.490.805.01-2	4	
33.690.805.01-2	6	

23.413.002.01-2

SCANBODY	PEEK PINS	TYPE
54.315.004.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	SCREWDRIVER	43.601.103.02-2

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.084.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.316.003.01-2	43.601.103.02-2

22.613.002.01-2
LAB SCANBODY
30.413.002.01-2

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm 2 mm 3 mm 4 mm

R	42.303.002.01-2	42.303.002.02-2	42.303.002.03-2	42.303.002.04-2
---	-----------------	-----------------	-----------------	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
0,3 mm				0,5 mm			
R	31.322.003.01-2	45°	30°	31.322.003.02-2	25°	30°	
NR	31.312.003.01-2			31.312.003.02-2			

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s
1 mm		CH+5mm	CH+7mm	CH+9mm
R	31.322.003.23-2	30°	25°	15°
NR	31.312.003.23-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.003.01-2	43.621.410.01-2 43.624.410.01-2	34.612.003.01-2	33.390.716.01-2	3	25°
52.410.104.01-2	10				33.490.716.01-2	4	
52.412.104.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.065.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.320.003.02-2	43.601.103.02-2

ANALOG

22.612.003.01-2

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c
	1 mm			2 mm			3 mm			4 mm		
R	31.323.004.01-2	45°	29°	31.323.004.02-2	30°	20°	31.323.004.03-2	25°	25°	31.323.004.04-2	20°	25°
NR	31.313.004.01-2			31.313.004.02-2			31.313.004.03-2			31.313.004.04-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT	α_s	α_s	α_s
	1 mm	CH-5mm	CH-7mm	CH-9mm
R	31.323.004.21-2	25°	20°	10°
NR	31.313.004.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
		50.313.004.01-2	43.621.410.01-2	34.613.004.01-2	33.390.754.01-2	3	25°	23.413.004.02-2	54.315.004.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.103.01-2	10	50.313.004.03-2 (IG=3mm)	43.624.410.01-2	34.613.004.02-2	33.490.754.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
52.412.103.01-2	12				33.690.754.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
										49.416.000.04-2 (13 mm)	D
										49.416.000.05-2 (13 mm)	E

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27	ANALOG	LAB SCANBODY
41.316.076.01-2	41.316.118.01-2	43.618.201.01-2	18	40.316.005.02-2	43.601.105.01-2	22.613.004.01-2	30.413.002.01-2
		43.624.201.01-2	24				
		43.632.201.01-2	32				

MULTI-UNIT

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	1 mm	2 mm	3 mm	4 mm
R	42.303.004.01-2	42.303.004.02-2	42.303.004.03-2	42.303.004.04-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	2.5/3.9 mm	3.5/4.9 mm	4.5/5.9 mm
NR	48.312.004.02-2	48.312.004.03-2	48.312.004.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s
1 mm				2 mm				3 mm			4 mm		
R	31.324.005.01-2	38°	23°	31.324.005.02-2	25°	15°		31.324.005.03-2	20°		31.324.005.04-2	15°	
NR	31.314.005.01-2			31.314.005.02-2				31.314.005.03-2			31.314.005.04-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s
1 mm		CH-5mm	CH-7mm	CH-9mm
R	31.324.005.21-2	25°	20°	10°
NR	31.314.005.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHA- NK	α_{di}
		50.314.005.01-2	43.621.410.01-2 43.624.410.01-2	34.614.005.01-2	33.390.958.01-2	3	30°
52.410.102.01-2	10	50.314.005.03-2			33.490.958.01-2	4	
52.412.102.01-2	12	(IG-3mm)			33.690.958.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
41.320.090.01-2	41.320.137.01-2	43.618.201.01-2	18	40.320.005.03-2	43.601.105.01-2
		43.624.201.01-2	24		
		43.632.201.01-2	32		

ANALOG

22.614.005.01-2

LAB SCANBODY

30.413.002.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.005.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.105.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c
	1.2 mm			2 mm			3 mm			4 mm		
R	31.322.006.01-2	40°	20°	31.322.006.02-2	25°	12°	31.322.006.03-2	20°	30°	31.322.006.04-2	15°	30°
NR	31.312.006.01-2			31.312.006.02-2			31.312.006.03-2			31.312.006.04-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT	α_s	α_s	α_s
	1.2 mm	CH-5mm	CH-7mm	CH-9mm
R	31.322.006.21-2	30°	20°	15°
NR	31.312.006.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
					33.330.734.01-2	3	25°	23.412.006.01-2	54.315.006.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.105.01-2	10	50.312.006.03-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.612.006.01-2	33.430.734.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
52.412.105.01-2	12				33.630.734.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
										49.416.000.04-2 (13 mm) 49.416.000.05-2 (13 mm)	D E

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27	
41.316.072.01-2	41.316.115.01-2	43.618.201.01-2	18	40.316.005.01-2	43.625.105.01-2	22.612.006.01-2
		43.624.201.01-2	24			LAB SCANBODY
		43.632.201.01-2	32			30.412.001.01-2

MULTI-UNIT

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	1 mm	2 mm	3 mm	4 mm
R	42.302.006.01-2	42.302.006.02-2	42.302.006.03-2	42.302.006.04-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	2.5/3.9 mm	3.5/4.9 mm	4.5/5.9 mm
NR	48.312.006.02-2	48.312.006.03-2	48.312.006.04-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT 1.5 mm	α_s	α_c	GINGIVAL HEIGHT 2 mm	α_s	α_c	GINGIVAL HEIGHT 3 mm	α_s
R	31.323.007.01-2	38°	17°	31.323.007.02-2	25°	12°	31.323.007.03-2	25°
NR	31.313.007.01-2			31.313.007.02-2			31.313.007.03-2	

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT 1.5 mm	α_s	α_c	α_s
	CH-5mm	CH-7mm	CH-9mm	
R	31.323.007.21-2	25°	20°	10°
NR	31.313.007.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.408.101.01-2	8	50.313.007.01-2	43.621.410.01-2 43.624.410.01-2	34.613.007.01-2
52.410.101.01-2	10			
52.412.101.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.350.775.01-2	3	25°
33.450.775.01-2	4	
33.650.775.01-2	6	

23.413.007.01-2

SCANBODY	PEEK PINS	TYPE
54.315.007.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.105.01-2

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.074.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.318.005.02-2	43.601.105.01-2

22.613.007.01-2
LAB SCANBODY
30.413.002.01-2

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.303.007.01-2	42.303.007.02-2	42.303.007.03-2	42.303.007.04-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT 1.5/2.9 mm	GINGIVAL HEIGHT 2.5/3.9 mm	GINGIVAL HEIGHT 3.5/4.9 mm
NR	48.312.007.01-2	48.312.007.02-2	48.312.007.03-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0.5 mm

R	31.323.008.01-2	45°	30°
---	-----------------	-----	-----

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s
0.5 mm CH-5mm CH-7mm CH-9mm

R	31.323.008.21-2	25°	20°	10°
---	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY REFERENCE

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.408.113.01-2	8	50.313.008.01-2	43.621.410.01-2 43.624.410.01-2	34.613.008.01-2

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.370.716.01-2	3	30°
33.470.716.01-2	4	
33.670.716.01-2	6	

23.413.008.01-2

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.045.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.318.005.01-2	43.601.105.01-2

LAB SCANBODY 30.412.001.01-2

SCANBODY	PEEK PINS	TYPE
54.322.008.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
	SCREWDRIVER	43.625.105.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c			GINGIVAL HEIGHT α_s α_c			GINGIVAL HEIGHT α_s		
0.3 mm			0.5 mm			1 mm		
R	31.322.009.01-2	45°	25°	31.322.009.02-2	25°	25°	31.322.009.03-2	25°
NR	31.312.009.01-2			31.312.009.02-2			31.312.009.03-2	

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.009.01-2	43.621.410.01-2 43.624.410.01-2	34.612.009.01-2	33.390.716.01-2	3	25°
52.410.114.01-2	10				33.490.716.01-2	4	
52.412.114.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20	ANALOG
41.320.051.01-2	43.618.201.01-2	18	40.320.003.01-2	43.601.103.02-2	22.612.009.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0.3 mm

R	31.323.010.01-2	45°	29°
NR	31.313.010.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.313.010.01-2	43.621.410.01-2	34.613.010.01-2	33.390.716.01-2	3	30°
52.410.115.01-2	10	50.313.010.04-2 (IG=3mm)	43.624.410.01-2		33.490.716.01-2	4	
52.412.115.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20	ANALOG
41.318.065.01-2	43.618.201.01-2	18	40.318.003.01-2	43.601.103.02-2	22.613.010.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.413.002.01-2

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.303.010.01-2	42.303.010.02-2	42.303.010.03-2	42.303.010.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0,3 mm

R	31.322.011.01-2	25°	29°
NR	31.312.011.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s

0,3 mm CH-5mm CH-7mm CH-9mm 2 mm CH-5mm CH-7mm CH-9mm 3 mm CH-5mm CH-7mm CH-9mm

NR	31.312.011.21-2	25°	25°	15°	31.312.011.23-2	25°	20°	15°	31.312.011.24-2	25°	20°	10°
----	-----------------	-----	-----	-----	-----------------	-----	-----	-----	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.108.01-2	10	50.312.011.01-2	43.621.410.01-2	34.612.011.01-2	33.345.804.01-2	3	20°
52.412.108.01-2	12		43.624.410.01-2		33.445.804.01-2	4	
					33.645.804.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.094.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.04-2	43.601.105.01-2

LAB SCANBODY

30.412.001.01-2

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm

2 mm

3 mm

R	42.302.011.01-2	42.302.011.02-2	42.302.011.03-2
---	-----------------	-----------------	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,3 mm

R	31.323.012.01-2	25°	30°
NR	31.313.012.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s
0,3 mm CH-5mm CH-7mm CH-9mm 2 mm CH-5mm CH-7mm CH-9mm 3 mm CH-5mm CH-7mm CH-9mm

NR	31.313.012.21-2	25°	25°	15°	31.313.012.23-2	25°	20°	15°	31.313.012.24-2	25°	20°	10°
----	-----------------	-----	-----	-----	-----------------	-----	-----	-----	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.109.01-2	10	50.313.012.01-2	43.621.410.01-2 43.624.410.01-2	34.613.012.01-2	33.345.804.01-2	3	20°
52.412.109.01-2	12				33.445.804.01-2	4	
					33.645.804.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.094.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.04-2	43.601.105.01-2

LAB SCANBODY

30.413.002.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.012.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER 43.625.105.01-2		

MULTI-UNIT

GINGIVAL HEIGHT 1 mm GINGIVAL HEIGHT 2 mm GINGIVAL HEIGHT 3 mm

R	42.303.012.01-2	42.303.012.02-2	42.303.012.03-2
---	-----------------	-----------------	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0.7 mm			
R	31.323.013.01-2	43°	23°
NR	31.313.013.01-2		

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.074.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.320.007.02-2	43.601.107.01-2

LAB SCANBODY
30.413.002.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c
	1.2 mm			2 mm			3 mm		
R	31.322.014.01-2	41°	23°	31.322.014.02-2	25°	17°		20°	25°
NR	31.312.014.01-2			31.312.014.02-2			31.312.014.03-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT	α_s	α_s	α_s
	3 mm	CH-5mm	CH-7mm	CH-9mm
R	31.322.014.23-2	25°	20°	15°
NR	31.312.014.23-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
52.410.128.01-2	10	50.312.014.03-2 (IG-3mm)	43.621.415.01-2	34.612.014.01-2	33.345.804.01-2	3	25°	23.412.014.01-2	54.315.014.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
					33.445.804.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
					33.645.804.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
										49.416.000.04-2 (13 mm)	D
										49.416.000.05-2 (13 mm)	E

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
41.314.067.01-2	41.314.105.01-2	43.618.201.01-2	18	40.314.003.04-2	43.601.103.02-2
		43.624.201.01-2	24		
		43.632.201.01-2	32		

LAB SCANBODY

30.412.001.01-2

SCREWDRIVER

43.601.103.02-2

MULTI-UNIT

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	2 mm	3 mm	4 mm
R	42.302.014.02-2	42.302.014.03-2	42.302.014.04-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	3.5 / 4.9 mm	4.5 / 5.9 mm
NR	48.312.014.03-2	48.312.014.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s
1.7 mm				2.5 mm				3 mm			4 mm			5 mm		
R	31.323.015.01-2	43°	23°	31.323.015.02-2	25°	15°	15°	31.323.015.03-2	25°	25°	31.323.015.04-2	20°	20°	31.323.015.05-2	15°	15°
NR	31.313.015.01-2			31.313.015.02-2				31.313.015.03-2			31.313.015.04-2			31.313.015.05-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s
1.7 mm		CH-5mm	CH-7mm	CH-9mm	2.5/3.9 mm		CH-5mm	CH-7mm	CH-9mm	3.5 mm		CH-5mm	CH-7mm	CH-9mm	1 mm		CH-5mm	CH-7mm	CH-9mm
R	31.323.015.21-2	30°	25°	10°	31.323.015.22-2	25°	20°	10°		25°	20°	10°		30°	25°	20°			
NR	31.313.015.27-2				31.313.015.22-2				31.313.015.26-2				31.313.015.27-2						

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
52.410.104.01-2	10	50.313.015.01-2	43.621.410.01-2	34.613.015.01-2	33.390.805.01-2	3	25°	23.413.015.01-2	54.315.015.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.412.104.01-2	12	50.313.015.03-2 (IG-3mm)	43.624.410.01-2		33.490.805.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
					33.690.805.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
										49.416.000.04-2 (13 mm)	D
										49.416.000.05-2 (13 mm)	E

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.075.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.318.003.02-2	43.601.103.02-2

LAB SCANBODY

30.412.001.01-2

SCREWDRIVER

43.601.103.02-2

MULTI-UNIT

GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
1 mm	2 mm	3 mm	4 mm
R	42.303.015.01-2	42.303.015.02-2	42.303.015.03-2
			42.303.015.04-2

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
1.5/2.9 mm	2.5/3.9 mm	4.5/5.9 mm
NR	48.312.015.01-2	48.312.015.02-2
		48.312.015.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c GINGIVAL HEIGHT α_s α_s

0,8 mm 1,5 mm

R	31.322.016.01-2	45°	28°	31.322.016.02-2	25°	25°
NR	31.312.016.01-2			31.312.016.02-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

1,5 mm CH-5mm CH-7mm CH-9mm

NR	31.312.016.22-2	25°	25°	15°
----	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.408.106.01-2	8	50.312.016.01-2 50.312.016.04-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.610.016.01-2	33.360.754.01-2	3	25°
52.410.106.01-2	10				33.460.754.01-2	4	
52.412.106.01-2	12				33.660.754.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.071.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.05-2	43.601.105.01-2

LAB SCANBODY

30.412.001.01-2

MULTI-UNIT

GINGIVAL HEIGHT

2 mm

R	42.302.016.02-2
---	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
0.7 mm				1.5 mm			
R	31.323.017.01-2	45°	24°	31.323.017.02-2	25°	15°	
NR	31.313.017.01-2			31.313.017.02-2			

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s
0.7 mm		CH=5mm	CH=7mm	CH=9mm
R	31.323.017.21-2	30°	25°	15°
NR	31.313.017.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.408.101.01-2	8	50.313.017.04-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.613.017.01-2	33.360.756.01-2	3	30°
52.410.101.01-2	10				33.460.756.01-2	4	
52.412.101.01-2	12				33.660.756.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.317.073.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.317.005.01-2	43.601.105.01-2

LAB SCANBODY

30.413.002.01-2

MULTI-UNIT

GINGIVAL HEIGHT

2 mm

R	42.303.017.02-2
---	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1.2 mm

R	31.324.018.01-2	39°	18°
NR	31.314.018.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.314.018.01-2	43.621.410.01-2 43.624.410.01-2	34.614.018.01-2	33.360.756.01-2	3	30°
52.410.102.01-2	10				33.460.756.01-2	4	
52.412.102.01-2	12				33.660.756.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27	LAB SCANBODY
41.317.073.01-2	43.618.201.01-2	18	40.317.005.01-2	43.601.105.01-2	30.413.002.01-2
	43.624.201.01-2	24			
	43.632.201.01-2	32			

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0,8 mm			
R	31.322.019.01-2	45°	30°
NR	31.312.019.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.019.01-2	43.621.410.01-2 43.624.410.01-2	34.612.019.01-2	33.360.754.01-2	3	25°
52.410.105.01-2	10				33.460.754.01-2	4	
52.412.105.01-2	12				33.660.754.01-2	6	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.071.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.05-2	43.601.105.01-2

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0.6 mm

R	31.323.020.01-2	45°	30°
---	-----------------	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

REFERENCE SCANBODY

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.408.112.01-2	8	50.313.020.01-2	43.620.411.01-2	34.613.020.01-2	33.390.716.01-2	3	30°
					33.490.716.01-2	4	
					33.690.716.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.044.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.06-2	43.601.105.01-2

LAB SCANBODY

30.413.005.01-2

SCANBODY	PEEK PINS	TYPE
54.322.020.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
	SCREWDRIVER	43.625.105.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
1.5 mm				2 mm				3 mm				4 mm				5 mm			
R	31.322.021.01-2	43°	24°	31.322.021.02-2	25°	20°	20°	31.322.021.03-2	20°	30°	30°	31.322.021.04-2	15°	30°	30°	31.322.021.05-2	15°	20°	20°
NR	31.312.021.01-2			31.312.021.02-2				31.312.021.03-2				31.312.021.04-2				31.312.021.05-2			

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s
1.5 mm		CH-5mm	CH-7mm	CH-9mm	3 mm		CH-5mm	H-7mm	CH-9mm
R	31.322.021.21-2	25°	20°	10°	31.322.021.23-2	25°	20°	15°	15°
NR	31.312.021.21-2				31.312.021.23-2				

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL		SHANK	α_{di}		SCANBODY	PEEK PINS		TYPE
		50.312.021.03-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.612.021.01-2	33.335.754.01-2	3	25°	23.412.021.01-2	54.315.015.21-2		49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A	
52.410.103.01-2	10				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B							
52.412.103.01-2	12				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C							
DYNAMIC SCREWS				STRAIGHT SCREWS				ANALOG			SCREWDRIVER		
DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	LAB SCANBODY					43.625.108.01-2		
41.316.073.01-2	41.316.108.01-2	43.618.201.01-2	18	40.316.008.02-2	43.601.108.01-2	22.612.021.01-2							
		43.624.201.01-2	24			30.412.001.01-2							
		43.632.201.01-2	32										

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	ANALOG
41.316.073.01-2	41.316.108.01-2	43.618.201.01-2	18	40.316.008.02-2	43.601.108.01-2	22.612.021.01-2
		43.624.201.01-2	24			LAB SCANBODY
		43.632.201.01-2	32			30.412.001.01-2

MULTI-UNIT

ANGULATED MULTI-UNIT 28°

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1 mm		2 mm		3 mm		4 mm		3.5/5.4 mm		4.5/6.4 mm	
R	42.302.021.01-2	42.302.021.02-2	42.302.021.03-2	42.302.021.04-2	NR	48.612.021.03-2	48.612.021.04-2				

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1.5/2.9 mm		2.5/3.9 mm		4.5/5.9 mm	
NR	48.312.021.01-2	48.312.021.02-2	48.312.021.04-2		

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
1.3 mm				2 mm				3 mm				4 mm				5 mm			
R	31.323.022.01-2	40°	19°	31.323.022.02-2	25°	14°	31.323.022.03-2	20°	30°	31.323.022.04-2	15	30	31.323.022.05-2	15°	20°				
NR	31.313.022.01-2			31.313.022.02-2			31.313.022.03-2			31.313.022.04-2			31.313.022.05-2						

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_c	α_s	α_c	α_s	α_c
1.3 mm		CH-5mm	CH-7mm	CH-9mm	3 mm		CH-5mm	CH-5mm	CH-7mm	CH-7mm	CH-9mm	CH-9mm
R	31.323.022.21-2	30°	25°	10°	31.323.022.23-2	20°	25°	20°	25°	10°	25°	
NR	31.313.022.21-2				31.313.022.23-2							

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY		HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL		SHANK	α_{di}	SCANALOG		SCANBODY	PEEK PINS		TYPE
52.408.101.01-2	8	50.313.022.01-2 50.313.022.03-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.613.022.01-2	33.335.758.01-2	3	30°	23.413.022.01-2		54.315.022.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A			
52.410.101.01-2	10				33.435.758.01-2	4						49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B		
52.412.101.01-2	12				33.635.758.01-2	6							49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C	
DYNAMIC SCREWS					STRAIGHT SCREWS						22.613.022.01-2	49.416.000.04-2 (13 mm)		D	
DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	ANALOG		30.413.002.01-2	LAB SCANBODY					49.416.000.05-2 (13 mm)	E
41.320.075.01-2	41.320.117.01-2	43.618.201.01-2	18	40.320.008.02-2	43.625.108.01-2	22.613.022.01-2				30.413.002.01-2			LAB SCANBODY		49.416.000.05-2 (13 mm)
		43.624.201.01-2	24												
		43.632.201.01-2	32												
												SCREWDRIVER		43.625.108.01-2	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	ANALOG	
41.320.075.01-2	41.320.117.01-2	43.618.201.01-2	18	40.320.008.02-2	43.625.108.01-2	22.613.022.01-2	
		43.624.201.01-2	24			LAB SCANBODY	
		43.632.201.01-2	32			30.413.002.01-2	

MULTI-UNIT

ANGULATED MULTI-UNIT 28°

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1 mm		2 mm		3 mm		4 mm		3.5/5.4 mm		4.5/6.4 mm	
R	42.303.022.01-2	42.303.022.02-2	42.303.022.03-2	42.303.022.04-2					NR	48.612.022.03-2	48.612.022.04-2

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1.5/2.9 mm		2.5/3.9 mm	
NR	48.312.022.01-2	48.312.022.02-2	

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0.3 mm

R	31.322.023.01-2	45°	30°
NR	31.312.023.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}	
		50.312.023.01-2	43.621.410.01-2 43.624.410.01-2	34.612.023.01-2	33.390.805.01-2	3	25°	23.412.023.01-2
52.410.103.01-2	10				33.490.805.01-2	4		
52.412.103.01-2	12				33.690.805.01-2	6		

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	ANALOG
41.316.059.01-2	43.618.201.01-2	18	40.316.008.01-2	43.625.108.01-2	22.612.023.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.412.001.01-2

Diagram illustrating the measurement of gingival height. The diagram shows a cross-section of a tooth with a gingival margin. The distance from the gingival margin to the gingival height is labeled as 0.5 mm. The angles are labeled as α_s and α_c .

DYNAMIC 3TIBASE

DYNAMIC SCANBODY (LAB/CLIN)

DYNAMIC MILLING TOOL

SCANBODY OP

DYNAMIC SCREWS

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER UNIGRIP
-------------------	------------------------

ANALOG

SCANBODY		PEEK PINS		TYPE
54.315.015.21-2	49.414.000.01-2 (6 mm)	49.415.000.01-2 (9 mm)	49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm)	49.415.000.02-2 (9 mm)	49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm)	49.415.000.03-2 (9 mm)	49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)			D
	49.416.000.05-2 (13 mm)			E
	SCREWDRIVER		43.625.108.01-2	

MULTI-UNIT

LIBRARY OPTIONS: GH = Gingival Height CH = Cement Height IG = Adaptor 3mm α_s = Standard maximum angulation α_c = Captive maximum angulation α_d = Direct to implant maximum angulation R = Rotational / Non-Engaging NR = Non Rotational / Engaging

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,3 mm

R	31.323.025.01-2	45°	30°
---	-----------------	-----	-----

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s
0,3 mm CH-5mm CH-7mm CH-9mm

R	31.323.025.21-2	30°	25°	10°
---	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

REFERENCE SCANBODY

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.408.112.01-2	8	50.313.025.02-2	43.620.411.01-2	34.613.025.01-2
			43.621.410.01-2	
52.410.111.01-2	10	50.313.025.01-2	43.624.410.01-2	

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.716.01-2	3	30°
33.490.716.01-2	4	
33.690.716.01-2	6	

23.413.025.01-2

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREW Ø2,6	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.039.01-2	41.314.050.31-2 (Temporary or Zirconio Direct MU)	43.618.201.01-2	18
		43.624.201.01-2	24
		43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER UNIGRIP
40.314.008.01-2	43.625.108.01-2

ANALOG

22.613.025.01-2
LAB SCANBODY
30.413.005.01-2

SCANBODY	PEEK PINS	TYPE
54.322.025.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
	SCREWDRIVER	43.625.108.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c GINGIVAL HEIGHT α_s α_c

0.5 mm 1.2 mm

R	31.322.026.01-2	45°	29°	31.322.026.02-2	25°	22°
NR	31.312.026.01-2			31.312.026.02-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0.5 mm CH-5mm CH-7mm CH-9mm

R	31.322.026.21-2	25°	20°	10°
NR	31.312.026.21-2			

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.108.01-2	10	50.312.026.04-2	43.621.410.01-2 43.624.410.01-2	34.612.026.01-2
52.412.108.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.805.01-2	3	25°
33.490.805.01-2	4	
33.690.805.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.075.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER UNIGRIP
40.318.008.01-2	43.625.108.01-2

ANALOG

22.612.026.01-2
LAB SCANBODY
30.412.001.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.026.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.108.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
0.3 mm				1.2 mm			
R	31.323.027.01-2	35°	29°	31.323.027.02-2	25°	22°	
NR	31.313.027.01-2			31.313.027.02-2			

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s
0.3 mm		CH-5mm	CH-7mm	CH-9mm
R	31.323.027.21-2	25°	20°	10°
NR	31.313.027.21-2			

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.109.01-2	10	50.314.027.01-2	43.621.410.01-2 43.624.410.01-2	34.614.027.01-2
52.412.109.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.958.01-2	3	30°
33.490.958.01-2	4	
33.690.958.01-2	6	

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.090.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER UNIGRIP
40.320.008.03-2	43.625.108.01-2

ANALOG

22.613.027.01-2

LAB SCANBODY

30.413.002.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.027.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.108.01-2

STANDARD DYNAMIC TIBASEGINGIVAL HEIGHT α_s α_c

0.3 mm

R	31.324.028.01-2	35°	30°
NR	31.314.028.01-2		

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.109.01-2	10	50.314.028.01-2	43.621.410.01-2 43.624.410.01-2	34.614.028.01-2
52.412.109.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.958.01-2	3	30°
33.490.958.01-2	4	
33.690.958.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.090.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER UNIGRIP
40.320.008.03-2	43.625.108.01-2

ANALOG

22.614.028.01-2
LAB SCANBODY
30.413.002.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.028.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.108.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c
	1.2 mm			2 mm			3 mm			4 mm		
R	31.322.029.01-2	30°	23°	31.322.029.02-2	25°	15°	31.322.029.03-2	20°	25°	31.322.029.04-2	15°	25°
NR	31.312.029.01-2			31.312.029.02-2			31.312.029.03-2			31.312.029.04-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT	α_s	α_s	α_s	GINGIVAL HEIGHT	α_s	α_s	α_s
	1.2 mm	CH-5mm	CH-7mm	CH-9mm	2 mm	CH-5mm	H-7mm	CH-9mm
R	31.322.029.21-2	25°	20°	20°		25°	20°	15°
NR	31.312.029.21-2				31.312.029.22-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
		50.312.029.01-2	43.621.410.01-2	34.613.029.01-2	33.345.804.01-2	3	20°	23.412.029.01-2	54.315.029.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.103.01-2	10	50.312.029.03-2 (IG=3mm)	43.624.410.01-2		33.445.804.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
52.412.103.01-2	12				33.645.804.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
										49.416.000.04-2 (13 mm)	D
										49.416.000.05-2 (13 mm)	E

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
41.316.094.01-2	41.316.132.01-2	43.618.201.01-2	18	40.316.003.02-2	43.601.103.02-2
		43.624.201.01-2	24		
		43.632.201.01-2	32		

LAB SCANBODY

30.412.001.01-2

SCREWDRIVER

43.601.103.02-2

MULTI-UNIT

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	1 mm	2 mm	3 mm	4 mm
R	42.302.029.01-2	42.302.029.02-2	42.302.029.03-2	42.302.029.04-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	2.5/3.9 mm	3.5/4.9 mm	4.5/5.9 mm
NR	48.312.029.02-2	48.312.029.03-2	48.312.029.04-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c
	1.1 mm			2 mm			3 mm			4 mm		
R	31.323.030.01-2	42°	25°	31.323.030.02-2	25°	15°	31.323.030.03-2	20°	30°	31.323.030.04-2	15°	30°
NR	31.313.030.01-2			31.313.030.02-2			31.313.030.03-2			31.313.030.04-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT	α_s	α_s	α_s	GINGIVAL HEIGHT	α_s	α_s	α_s	GINGIVAL HEIGHT	α_s	α_s	α_s
	1.1 mm	CH-5mm	CH-7mm	CH-9mm	2 mm	CH-5mm	CH-7mm	CH-9mm	3 mm	CH-5mm	CH-7mm	CH-9mm
R	31.323.030.21-2	25°	20°	10°	31.323.030.22-2	25°	20°	15°	31.323.030.23-2	25°	20°	10°
NR	31.313.030.21-2				31.313.030.22-2				31.313.030.23-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
52.408.101.01-2	8	50.313.030.01-2 50.313.030.03-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.613.030.01-2	33.345.808.01-2	3	30°	23.413.030.01-2	54.315.030.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.101.01-2	10				33.445.808.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
52.412.101.01-2	12				33.645.808.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
DYNAMIC SCREWS				STRAIGHT SCREWS				LAB SCANBODY 30.413.002.01-2		49.416.000.04-2 (13 mm)	D
DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20	49.416.000.05-2 (13 mm)	E				
41.320.079.01-2	41.320.125.01-2	43.618.201.01-2	18	40.320.003.04-2	43.601.103.02-2	SCREWDRIVER 1.20	43.601.103.02-2				
		43.624.201.01-2	24								
		43.632.201.01-2	32								

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20	LAB SCANBODY
41.320.079.01-2	41.320.125.01-2	43.618.201.01-2	18	40.320.003.04-2	43.601.103.02-2	30.413.002.01-2
		43.624.201.01-2	24			
		43.632.201.01-2	32			

MULTI-UNIT

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	1 mm	2 mm	3 mm	4 mm	5 mm
R	42.303.030.01-2	42.303.030.02-2	42.303.030.03-2	42.303.030.04-2	42.303.030.05-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	1.5/2.9 mm	2.5/3.9 mm	3.5/4.9 mm	4.5/5.9 mm
NR	48.312.030.01-2	48.312.030.02-2	48.312.030.03-2	48.312.030.04-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c
	1.3 mm			2 mm			3 mm			4 mm		
R	31.322.033.01-2	38°	18°	31.322.033.02-2	20°	14°	31.322.033.03-2	15°	25°	31.322.033.04-2	15°	25°
NR	31.312.033.01-2			31.312.033.02-2			31.312.033.03-2			31.312.033.04-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT	α_s	α_s	α_s	GINGIVAL HEIGHT	α_s	α_s	α_s
	1.3 mm	CH-5mm	CH-7mm	CH-9mm	3 mm	CH-5mm	CH-7mm	CH-9mm
R	31.322.033.21-2	25°	20°	10°	31.322.033.23-2	20°	15°	10°
NR	31.312.033.21-2				31.312.033.23-2			

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL
SCANALOG
SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
52.408.106.01-2	8	50.312.033.03-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.612.033.01-2	33.315.804.01-2	3	25°	23.412.033.01-2	54.315.033.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.106.01-2	10				33.415.804.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
52.412.106.01-2	12				33.615.804.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
DYNAMIC SCREWS				STRAIGHT SCREWS							
DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREW/DRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6	ANALOG					
41.316.078.01-2	41.316.124.01-2	43.618.201.01-2	18	40.316.007.01-2	43.601.107.01-2	22.612.033.01-2					
		43.624.201.01-2	24			LAB SCANBODY					
		43.632.201.01-2	32			30.412.001.01-2					
								SCREWDRIVER		43.601.107.01-2	

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6
41.316.078.01-2	41.316.124.01-2	43.618.201.01-2	18	40.316.007.01-2	43.601.107.01-2
		43.624.201.01-2	24		
		43.632.201.01-2	32		

ANALOG

22.612.033.01-2

LAB SCANBODY

30.412.001.01-2

SCREWDRIVER

43.601.107.01-2

MULTI-UNIT

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	1 mm	2 mm	3 mm	4 mm
R	42.302.033.01-2	42.302.033.02-2	42.302.033.03-2	42.302.033.04-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	2.5/3.9 mm	3.5/4.9 mm	4.5/5.9 mm
NR	48.312.033.02-2	48.312.033.03-2	48.312.033.04-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c	GINGIVAL HEIGHT	α_s	α_c
	1.1 mm			2 mm			3 mm			4 mm		
R	31.323.035.01-2	39°	18°	31.323.035.02-2	20°	14°	31.323.035.03-2	15°	30°	31.323.035.04-2	15°	30°
NR	31.313.035.01-2			31.313.035.02-2			31.313.035.03-2			31.313.035.04-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT	α_s	α_s	α_s	GINGIVAL HEIGHT	α_s	α_s	α_s
	1.1 mm	CH-5mm	CH-7mm	CH-9mm	3 mm	CH-5mm	CH-7mm	CH-9mm
R	31.323.035.21-2	25°	20°	10°	31.323.035.23-2	20°	15°	10°
NR	31.313.035.21-2				31.313.035.23-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
		50.313.035.01-2	43.621.410.01-2	34.613.035.01-2	33.315.804.01-2	3	25°	23.413.035.01-2	54.315.035.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.107.01-2	10	50.313.035.03-2	43.624.410.01-2		33.415.804.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
52.412.107.01-2	12	(IG-3mm)			33.615.804.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
										49.416.000.04-2 (13 mm)	D
										49.416.000.05-2 (13 mm)	E

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6	ANALOG	LAB SCANBODY
41.316.078.01-2	41.316.124.01-2	43.618.201.01-2	18	40.316.007.01-2	43.601.107.01-2	22.613.035.01-2	30.413.002.01-2
		43.624.201.01-2	24				
		43.632.201.01-2	32				

MULTI-UNIT

	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	1 mm	2 mm	3 mm	4 mm
R	42.303.035.01-2	42.303.035.02-2	42.303.035.03-2	42.303.035.04-2

ANGULATED MULTI-UNIT 20°

	GINGIVAL HEIGHT	GINGIVAL HEIGHT
	2.5/3.9 mm	3.5/4.9 mm
NR	48.312.035.02-2	48.312.035.03-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0,6 mm

R	31.323.037.01-2	45°	25°
NR	31.313.037.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0,6 mm

CH=5mm

CH=7mm

CH=9mm

R	31.323.037.21-2	30°	25°	15°
NR	31.313.037.21-2			

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL
SCANALOG

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}	
52.408.110.01-2	8	50.313.037.04-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.613.037.01-2	33.315.708.01-2	3	30°	23.413.037.01-2
52.410.110.01-2	10				33.415.708.01-2	4		
52.412.110.01-2	12				33.615.708.01-2	6		

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.067.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER TORX T6	ANALOG
40.320.007.01-2	43.601.107.01-2	22.613.037.01-2
		LAB SCANBODY
		30.413.004.01-2

MULTI-UNIT

GINGIVAL HEIGHT

1 mm

R	42.303.037.01-2
---	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,7 mm

R	31.322.038.01-2	45°	29°
NR	31.312.038.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s
0,7 mm CH-5mm CH-7mm CH-9mm

R	31.322.038.21-2	30°	25°	10°
NR	31.312.038.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.038.01-2	43.621.410.01-2 43.624.410.01-2	34.612.038.01-2	33.345.804.01-2	3	25°
52.410.103.01-2	10				33.445.804.01-2	4	
52.412.103.01-2	12				33.645.804.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.081.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.316.004.02-2	43.601.104.01-2

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s
0,7 mm				2 mm		
R	31.323.039.01-2	45°	29°	31.323.039.03-2	25°	
NR	31.313.039.01-2			31.313.039.03-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s
0,7 mm		CH-5mm	CH-7mm	CH-9mm
R	31.323.039.21-2	30°	25°	10°
NR	31.313.039.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.114.01-2	10	50.313.039.01-2	43.621.410.01-2 43.624.410.01-2	34.613.039.01-2	33.345.856.01-2	3	25°
52.412.114.01-2	12				33.445.856.01-2	4	
					33.645.856.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.081.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.316.004.02-2	43.601.104.01-2

LAB SCANBODY

30.413.002.01-2

MULTI-UNIT

GINGIVAL HEIGHT

1 mm	
R	42.303.039.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
0,6 mm				1,5 mm				3 mm				4 mm			
R	31.322.040.01-2	45°	30°	31.322.040.02-2	25°	25°	31.322.040.03-2	20°	30°	31.322.040.04-2	15°	30°	31.322.040.05-2	10°	23°
NR	31.312.040.01-2			31.312.040.02-2			31.312.040.03-2			31.312.040.04-2			31.312.040.05-2		
NR (Friction Fit)	31.312.042.01-2														

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT	GINGIVAL HEIGHT	α_s	α_s	α_s		
0.6 mm		CH+5mm	CH+7mm	CH+9mm	1 mm		CH+5mm	CH+7mm	CH+9mm	2 mm		3 mm		CH+5mm	CH+7mm	CH+9mm
R	31.322.040.21-2	25°	20°	10°	31.322.040.29-2	30°	25°	20°			31.322.040.23-2	25°	20°	15°		
NR	31.312.040.21-2				31.312.040.29-2				31.312.040.28-2	31.312.040.23-2						

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL
SCANALOG
SCANBODY OP

SCANBODY		HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}			SCANBODY	PEEK PINS	TYPE					
52.408.101.01-2	8	50.312.040.03-2 (IG+3mm)	43.621.410.01-2 43.624.410.01-2	34.612.040.01-2		33.370.716.01-2	3	25°	23.412.040.01-2	54.315.040.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A						
52.410.101.01-2	10					33.470.716.01-2	4						49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B				
52.412.101.01-2	12					33.670.716.01-2	6											
DYNAMIC SCREWS						STRAIGHT SCREWS												
DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER		SCREWDRIVER LENGTH (mm)		STRAIGHT SCREW	SCREWDRIVER Hex. 1.27	ANALOG										
41.317.071.01-2	41.317.106.01-2	43.618.201.01-2	18	40.317.004.01-2	43.601.104.01-2	22.612.040.01-2												
		43.624.201.01-2	24			LAB SCANBODY												
		43.632.201.01-2	32			30.412.001.01-2												
								SCREWDRIVER		43.625.105.01-2								

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.317.071.01-2	41.317.106.01-2	43.618.201.01-2	18
		43.624.201.01-2	24
		43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.317.004.01-2	43.601.104.01-2

ANALOG

22.612.040.01-2
LAB SCANBODY
30.412.001.01-2

MULTI-UNIT
ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
1 mm	2 mm	3 mm	4 mm
R	42.302.040.01-2	42.302.040.02-2	42.302.040.03-2
			42.302.040.04-2

GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
2,5 mm	3,5 mm	4,5 mm
NR	48.312.040.02-2	48.312.040.03-2
		48.312.040.04-2

STANDARD DYNAMIC TIBASE											
GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
0.6 mm				1.5 mm				3 mm			
R	31.322.040.01-2	45°	30°	31.322.040.02-2	25°	25°		31.322.040.03-2	20°	31.322.040.04-2	15°
NR	31.312.040.01-2			31.312.040.02-2				31.312.040.03-2		31.312.040.04-2	
NR (Friction Fit)	31.312.042.01-2										

DYNAMIC 3TIBASE													
GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s
0.6 mm		CH-5mm	CH-7mm	CH-9mm	1 mm		CH-5mm	CH-7mm	CH-9mm	2 mm		3 mm	
R	31.322.040.21-2	25°	20°	10°	31.322.040.29-2	30°	25°	20°			31.322.040.23-2	25°	20°
NR	31.312.040.21-2				31.312.040.29-2					31.312.040.28-2			

DYNAMIC MILLING TOOL		
DYNAMIC MILLING TOOL	SHANK	α_{di}
33.370.716.01-2	3	25°
33.470.716.01-2	4	
33.670.716.01-2	6	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.071.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
0.4 mm				1.5 mm			
R	31.323.041.01-2	45°	30°		31.323.041.02-2	30°	25°
NR	31.313.041.01-2				31.313.041.02-2		
NR (Friction Fit)	31.313.043.01-2						

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s
0.4 mm		CH-5mm	CH-7mm	CH-9mm	1.5 mm		CH-5mm	CH-7mm	CH-9mm
R	31.323.041.21-2	30°	20°	10°		31.323.041.22-2	30°	25°	15°
NR	31.313.041.21-2					31.313.041.22-2			

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL****SCANALOG****SCANBODY OP**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.313.041.01-2			33.370.716.01-2	3	30°
52.410.102.01-2	10	50.313.041.03-2 (IG=3mm)	43.621.410.01-2	34.613.041.01-2	33.470.716.01-2	4	
52.412.102.01-2	12		43.624.410.01-2		33.670.716.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
41.317.071.01-2	41.317.106.01-2	43.618.201.01-2	18	40.317.004.01-2	43.601.104.01-2
		43.624.201.01-2	24		
		43.632.201.01-2	32		

ANALOG

22.613.041.01-2

LAB SCANBODY

30.413.002.01-2

SCANBODY**PEEK PINS**

TYPE

54.315.041.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.105.01-2

STANDARD DYNAMIC TIBASE						
GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s α_c
0.4 mm				1.5 mm		
R	31.323.041.01-2	45°	30°	31.323.041.02-2	30°	25°
NR	31.313.041.01-2			31.313.041.02-2		

DYNAMIC 3TIBASE				
GINGIVAL HEIGHT		α_s	α_s	α_s
0.4 mm		CH-5mm	CH-7mm	CH-9mm
R	31.323.041.21-2	30°	20°	10°
NR	31.313.041.21-2			

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.071.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL		
DYNAMIC MILLING TOOL	SHANK	α_{di}
33.370.716.01-2	3	30°
33.470.716.01-2	4	
33.670.716.01-2	6	

LAB SCANBODY
30.413.002.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT 1 mm	α_s	α_c	GINGIVAL HEIGHT 2 mm	α_s
R	31.322.044.01-2	42°	23°	31.322.044.02-2	25°
NR	31.312.044.01-2			31.312.044.02-2	

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT 1 mm	α_s	CH-5mm α_s	CH-7mm α_s	CH-9mm α_s
R	31.322.044.21-2	25°	20°	10°	
NR	31.312.044.21-2				

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.044.01-2	43.621.410.01-2 43.624.410.01-2	34.612.044.01-2	33.390.716.01-2	3	25°
52.410.105.01-2	10				33.490.716.01-2	4	
52.412.105.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.065.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.318.003.01-2	43.601.103.02-2

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s
1 mm				2 mm		
R	31.323.045.01-2	43°	22°	31.323.045.02-2	25°	
NR	31.313.045.01-2			31.313.045.02-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s
1 mm		CH-5mm	CH-7mm	CH-9mm
R	31.323.045.21-2	25°	20°	10°
NR	31.313.045.21-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.313.045.01-2	43.621.410.01-2 43.624.410.01-2	34.613.045.01-2	33.390.716.01-2	3	30°
52.410.118.01-2	10				33.490.716.01-2	4	
52.412.118.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.065.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.318.003.01-2	43.601.103.02-2

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASEGINGIVAL HEIGHT α_s α_c

1 mm

R	31.324.046.01-2	42°	21°
NR	31.314.046.01-2		

DYNAMIC 3TIBASEGINGIVAL HEIGHT α_s α_s α_s

1 mm

CH-5mm

CH-7mm

CH-9mm

R	31.324.046.21-2	30°	20°	10°
NR	31.314.046.21-2			

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.314.046.01-2	43.621.410.01-2 43.624.410.01-2	34.614.046.01-2	33.390.716.01-2	3	30°
52.410.125.01-2	10				33.490.716.01-2	4	
					33.690.716.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
41.318.065.01-2	43.618.201.01-2	18	40.318.003.01-2	43.601.103.02-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0,6 mm

R	31.322.047.01-2	45°	30°
NR	31.312.047.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0,6 mm

CH-5mm

CH-7mm

CH-9mm

NR	31.312.047.21-2	30°	25°	20°
----	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.409.123.01-2	9	50.312.047.01-2	43.621.410.01-2 43.624.410.01-2	34.612.047.01-2	33.390.716.01-2	3	25°
52.410.123.01-2	10				33.490.716.01-2	4	
52.412.123.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6
41.320.074.01-2	43.618.201.01-2	18	40.320.007.02-2	43.601.107.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASEGINGIVAL HEIGHT α_s α_c

0,6 mm

R	31.323.048.01-2	45°	30°
NR	31.313.048.01-2		

DYNAMIC 3TIBASEGINGIVAL HEIGHT α_s α_s α_s

0,6 mm CH-5mm CH-7mm CH-9mm

R	31.323.048.21-2	30°	25°	20°
NR	31.313.048.21-2			

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.409.123.01-2	9	50.312.047.01-2	43.621.410.01-2 43.624.410.01-2	34.612.047.01-2	33.390.716.01-2	3	30°
52.410.123.01-2	10				33.490.716.01-2	4	
52.412.123.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6
41.320.074.01-2	43.618.201.01-2	18	40.320.007.02-2	43.601.107.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0.6 mm			
R	31.321.049.01-2	45°	30°
NR	31.311.049.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.311.049.01-2	43.621.410.01-2 43.624.410.01-2	34.611.049.01-2	33.325.472.01-2'	3	25°
52.410.116.01-2	10				33.425.472.01-2'	4	
52.412.116.01-2	12				33.625.472.01-2'	6	

* Only for titanium and soft materials

DYNAMIC SCREWS			STRAIGHT SCREWS	
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
41.314.064.01-2	43.618.201.01-2	18	40.314.004.01-2	43.601.104.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c GINGIVAL HEIGHT α_s

0.5 mm 2 mm

R	31.323.051.01-2	45°	27°	31.323.051.03-2	25°
NR	31.313.051.01-2			31.313.051.03-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

2 mm CH-5mm CH-7mm CH-9mm

R	31.323.051.23-2	25°	20°	15°
NR	31.313.051.23-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.050.04-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.612.050.01-2	33.335.676.01-2	3	25°
52.410.117.01-2	10				33.435.676.01-2	4	
52.412.117.01-2	12				33.635.676.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.064.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.318.004.03-2	43.601.104.01-2

LAB SCANBODY

30.412.001.01-2

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm 2 mm

R	42.302.050.01-2	42.302.050.02-2
---	-----------------	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		α_c	GINGIVAL HEIGHT α_s		
0.5 mm			2 mm		
R	31.323.051.01-2	45°	25°	31.323.051.03-2	25°
NR	31.313.051.01-2			31.313.051.03-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s
2 mm		CH-5mm	CH-7mm	CH-9mm
R	31.323.051.23-2	25°	20°	15°
NR	31.313.051.23-2			

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.313.051.01-2	43.621.410.01-2 43.624.410.01-2	34.613.051.01-2	33.335.676.01-2	3	25°
52.410.118.01-2	10	50.313.051.04-2			33.435.676.01-2	4	
52.412.118.01-2	12	(IG=3mm)			33.635.676.01-2	6	

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.064.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.318.004.03-2	43.601.104.01-2

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0.3 mm

R	31.324.052.01-2	45°	27°
NR	31.314.052.01-2		

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.314.052.04-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.614.052.01-2	33.335.676.01-2	3	30°
52.410.102.01-2	10				33.435.676.01-2	4	
52.412.102.01-2	12				33.635.676.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
41.318.064.01-2	43.618.201.01-2	18	40.318.004.03-2	43.601.104.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0,6 mm			
R	31.323.054.01-2	45°	25°
NR	31.313.054.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.314.054.01-2	43.621.410.01-2 43.624.410.01-2	34.614.054.01-2	33.345.856.01-2	3	30°
52.410.119.01-2	10				33.445.856.01-2	4	
52.412.119.01-2	12				33.645.856.01-2	6	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.067.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS
40.318.012.01-2

LAB SCANBODY
30.413.002.01-2

STANDARD DYNAMIC TIBASEGINGIVAL
HEIGHT α_s α_c

0.5 mm

R	31.324.057.01-2	45°	27°
NR	31.314.057.01-2		

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.408.101.01-2	8	50.314.057.01-2	43.621.410.01-2 43.624.410.01-2	34.614.057.01-2	33.390.805.01-2	3	30°
52.410.101.01-2	10				33.490.805.01-2	4	
52.412.101.01-2	12				33.690.805.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20	ANALOG
41.316.084.01-2	43.618.201.01-2	18	40.316.003.01-2	43.601.103.02-2	22.614.057.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.414.003.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0.5 mm			
R	31.324.058.01-2	45°	30°
NR	31.314.058.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.314.058.01-2	43.621.410.01-2 43.624.410.01-2	34.614.058.01-2	33.390.716.01-2	3	30°
52.410.118.01-2	10				33.490.716.01-2	4	
52.412.118.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.047.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.320.003.01-2	43.601.103.02-2

ANALOG
22.614.058.01-2
LAB SCANBODY
30.414.003.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,7 mm

R	31.324.059.01-2	45°	27°
NR	31.314.059.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR
		50.313.010.01-2	43.621.410.01-2 43.624.410.01-2
52.410.115.01-2	10	50.313.010.04-2 (IG=3mm)	
52.412.115.01-2	12		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.065.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.716.01-2	3	25
33.490.716.01-2	4	
33.690.716.01-2	6	

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.318.003.01-2	43.601.103.02-2

LAB SCANBODY

30.414.003.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0.5 mm

R	31.324.060.01-2	45°	30°
NR	31.314.060.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.122.01-2	10	50.314.060.01-2	43.621.410.01-2 43.624.410.01-2	34.614.060.01-2	33.390.716.01-2	3	30°
					33.490.716.01-2	4	
					33.690.716.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.20	ANALOG
41.320.060.01-2	43.618.201.01-2	18	40.320.003.02-2	43.601.103.02-2	22.614.060.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.415.007.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0.6 mm

R	31.324.061.01-2	45°	30°
NR	31.314.061.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.314.061.01-2	43.621.410.01-2 43.624.410.01-2	34.614.061.01-2	33.390.958.01-2	3	30°
52.410.125.01-2	10				33.490.958.01-2	4	
					33.690.958.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	ANALOG
41.325.067.01-2	43.618.201.01-2	18	40.325.008.01-2	43.625.108.01-2	22.614.060.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.415.007.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT				α_s	α_c
0.6 mm					
R	31.323.066.01-2	45°	30°		

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.039.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER UNIGRIP
40.314.008.01-2	43.625.108.01-2

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,6 mm

R	31.323.074.01-2	45°	30°
NR	31.313.074.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}	
52.408.110.01-2	8	50.313.074.01-2	43.621.410.01-2 43.624.410.01-2	34.613.074.01-2	33.330.708.01-2	3	30°	23.413.074.01-2
52.410.110.01-2	10				33.430.708.01-2	4		
52.412.110.01-2	12				33.630.708.01-2	6		

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Sq. 1.30	ANALOG
41.320.050.01-2	43.618.201.01-2	18	40.320.007.04-2	43.601.107.01-2	22.613.074.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.415.007.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s
1 mm				2 mm				3 mm			4 mm		
R	31.322.075.01-2	42°	24°	31.322.075.02-2	20°	15°	31.322.075.03-2	20°	31.322.075.04-2	15°			
NR				31.312.075.02-2			31.312.075.03-2		31.312.075.04-2				

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s					
1 mm		CH-5mm	CH-7mm	CH-9mm	
R	31.322.075.21-2	30°	20°	15°	

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.105.01-2	10	50.312.075.03-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.612.075.01-2	33.330.734.01-2	3	25°
52.412.105.01-2	12				33.430.734.01-2	4	
					33.630.734.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.077.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW

40.318.013.01-2

ANALOG

22.612.075.01-2

LAB SCANBODY

30.412.001.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.015.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E

MULTI-UNIT

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
2 mm		3 mm		4 mm	
R	42.302.075.02-2	42.302.075.03-2	42.302.075.04-2		

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,4 mm

R	31.324.080.01-2	45°	30°
NR	31.314.080.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.124.01-2	10	50.314.080.01-2	43.621.410.01-2 43.624.410.01-2	34.614.080.01-2
52.412.124.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.370.716.01-2	3	30°
33.470.716.01-2	4	
33.670.716.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.317.071.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.317.004.01-2	43.601.104.01-2

ANALOG

22.614.080.01-2
LAB SCANBODY
30.414.003.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.080.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.104.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0,4 mm			
R	31.325.081.01-2	41°	18°
NR	31.315.081.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.315.081.01-2	43.621.410.01-2 43.624.410.01-2	34.615.081.01-2	33.335.676.01-2	3	30°
52.410.126.01-2	10				33.435.676.01-2	4	
52.412.126.01-2	12				33.635.676.01-2	6	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.064.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.318.004.03-2	43.601.104.01-2

LAB SCANBODY
30.414.003.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1,2 mm

R	31.322.082.01-2	45°	25°
NR	31.312.082.01-2		

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.082.01-2	43.621.410.01-2 43.624.410.01-2	34.612.082.01-2	33.345.804.01-2	3	25°
52.410.105.01-2	10				33.445.804.01-2	4	
52.412.105.01-2	12				33.645.804.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.074.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW
40.316.012.01-2

LAB SCANBODY

30.412.001.01-2

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm 2 mm 3 mm 4 mm

R	42.302.082.01-2	42.302.082.02-2	42.302.082.03-2	42.302.082.04-2
---	-----------------	-----------------	-----------------	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
1,2 mm

R	31.323.083.01-2	45°	25°
NR	31.313.083.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.103.01-2	10	50.313.083.01-2	43.621.410.01-2 43.624.410.01-2	34.613.083.01-2	33.345.856.01-2	3	30°
52.412.103.01-2	12				33.445.856.01-2	4	
					33.645.856.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.076.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT
SCREW

40.318.012.02-2

LAB SCANBODY

30.413.002.01-2

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm 2 mm 3 mm 4 mm

R	42.302.083.01-2	42.302.083.02-2	42.302.083.03-2	42.302.083.04-2
---	-----------------	-----------------	-----------------	-----------------

STANDARD DYNAMIC TIBASEGINGIVAL HEIGHT α_s

1 mm

R	31.321.084.01-2	40°
NR	31.311.084.01-2	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.076.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER Star 1.50
40.314.003.03-2	43.601.103.02-2

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0.3 mm

R	31.324.085.01-2	45°	30°
NR	31.314.085.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0.3 mm

CH-5mm

CH-7mm

CH-9mm

NR	31.314.085.21-2	30°	25°	20°
----	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.314.085.01-2	43.621.410.01-2 43.624.410.01-2	34.614.085.01-2	33.345.856.01-2	3	25°
52.410.117.01-2	10				33.445.856.01-2	4	
52.412.117.01-2	12				33.645.856.01-2	6	

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
41.316.081.01-2	43.618.201.01-2	18	40.316.004.02-2	43.601.104.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s
1 mm

R	31.325.086.01-2	45°
NR	31.315.086.01-2	

DYNAMIC SCANBODY (LAB/CLIN)

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR
52.410.117.01-2	10	50.314.085.01-2	43.621.410.01-2
52.412.117.01-2	12		43.624.410.01-2

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.081.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.345.856.01-2	3	25°
33.445.856.01-2	4	
33.645.856.01-2	6	

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.316.004.02-2	43.601.104.01-2

LAB SCANBODY

30.415.007.01-2

STANDARD DYNAMIC TIBASE
GINGIVAL HEIGHT α_s

0,3 mm

R	31.321.087.01-2	25°
NR	31.311.087.01-2	

DYNAMIC 3TIBASE
GINGIVAL HEIGHT
GINGIVAL HEIGHT
 α_s
 α_s
 α_s

0,3 mm

2 mm

CH-5mm

CH-7mm

CH-9mm

NR	31.311.087.21-2	31.311.087.23-2	25°	20°	15°
----	-----------------	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR
52.410.132.01-2	10	50.311.087.04-2 (IG-3mm)	43.621.415.01-2
52.412.132.01-2	12		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.094.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.04-2	43.601.105.01-2

LAB SCANBODY

30.410.006.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.087.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.105.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0,3 mm

R	31.324.088.01-2	25°	25°
NR	31.314.088.01-2		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.094.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.320.005.04-2	43.601.105.01-2

LAB SCANBODY

30.414.003.01-2

SCANBODY OP

SCANBODY

PEEK PINS

TYPE

54.315.088.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		54.315.088.21-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		α_c	GINGIVAL HEIGHT α_s		α_c	GINGIVAL HEIGHT α_s		α_c
1 mm			2 mm			3 mm		
R	31.321.091.01-2	45°	24°	31.311.090.02-2	25°	31.321.091.03-2	20°	
NR	31.321.091.01-2			31.311.090.02-2		31.314.091.03-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s		α_s	α_s	α_s
1 mm		CH-5mm	CH-7mm	CH-9mm
NR	31.311.090.21-2	30°	25°	15°

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.102.01-2	10	50.311.090.03-2 (IG=3mm)	43.621.415.01-2	34.611.090.01-2	33.390.958.01-2	3	25°
52.412.102.01-2	12				33.490.958.01-2	4	
					33.690.958.01-2	6	

*Only for R

*Only for titanium and soft materials

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.074.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.314.005.01-2	43.601.105.01-2

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT 1 mm	α_s	α_c	GINGIVAL HEIGHT 2 mm	α_s	GINGIVAL HEIGHT 3 mm	α_s	α_c	GINGIVAL HEIGHT 4 mm	α_s	α_c
R	31.324.091.01-2	38°	18°	31.324.091.02-2	25°	31.324.091.03-2	20°	25°	31.324.091.04-2	15°	25°
NR	31.324.091.01-2			31.324.091.02-2		31.324.091.03-2			31.324.091.04-2		

DYNAMIC 3TIBASE

	GINGIVAL HEIGHT 1 mm	α_s	α_s	α_s
	1 mm	CH-5mm	CH-7mm	CH-9mm
NR	31.314.091.21-2	30°	25°	15°

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.102.01-2	10	50.314.091.03-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.614.09.01-2	33.390.958.01-2	3	30°
52.410.102.01-2	12				33.490.958.01-2	4	
52.410.102.01-2	12				33.690.958.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
41.320.082.01-2	43.618.201.01-2	18	40.320.005.01-2	43.601.105.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.303.091.01-2	42.303.091.01-2	42.303.091.01-2	42.303.091.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		α_c	GINGIVAL HEIGHT α_s		α_s	GINGIVAL HEIGHT α_s	
1 mm			2 mm			3 mm	
R	31.315.092.01-2	45°	25°	31.315.092.02-2	25°	31.325.092.03-2	20°
NR	31.315.092.01-2			31.315.092.02-2		31.325.092.03-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s		α_s	α_s	α_s
1 mm		CH-5mm	CH-7mm	CH-9mm
NR	31.315.092.21-2	30°	25°	15°

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.315.092.01-2	43.621.410.01-2	34.615.092.01-2	33.390.958.01-2	3	30°
52.410.129.01-2	10	50.315.092.03-2 (IG=3mm)	43.624.410.01-2		33.490.958.01-2	4	
					33.690.958.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.082.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.320.005.01-2	43.601.105.01-2

LAB SCANBODY

30.415.007.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,6 mm

R	31.324.096.01-2	45°	30°
NR	31.324.096.01-2		

D. S.(LAB/CLIN)

DIGITAL ANALOG

SCREWDRIVER
ADAPTORDIGITAL
ANALOG

43.621.410.01-2
43.624.410.01-2

34.614.096.01-2

DYNAMIC SCREWS

DYNAMIC
SCREWDYNAMIC
SCREWDRIVERSCREWDRIVER
LENGTH (mm)

41.320.067.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL

DYNAMIC
MILLING TOOL

SHANK

 α_{di}

33.315.708.01-2	3	30°
33.415.708.01-2	4	
33.615.708.01-2	6	

STRAIGHT SCREWS

STRAIGHT
SCREWSCREWDRIVER
Hex. 1.27

40.320.007.01-2	43.601.107.01-2
-----------------	-----------------

SCANALOG

23.414.096.01-2

ANALOG

22.614.096.01-2
LAB SCANBODY 30.414.008.01-2

SCANBODY OP

SCANBODY

PEEK PINS

TYPE

54.315.096.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.107.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0.5 mm

R	31.323.101.01-2	45°	30°
---	-----------------	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL			SCANALOG	REFERENCE SCANBODY			
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS		TYPE
52.409.133.01-2	9	50.313.101.01-2	50.313.101.01-2 43.624.410.01-2	34.613.101.01-2	33.335.676.01-2	3	30°	23.413.101.01-2	54.322.101.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A	
					33.435.676.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B	
					33.635.676.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C	
DYNAMIC SCREWS				STRAIGHT SCREWS				LAB SCANBODY 30.413.005.01-2		49.416.000.04-2 (13 mm) 49.416.000.05-2 (13 mm)	D E	
DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27	CAPS	mm					
41.314.043.01-2	41.314.050.31-2	43.618.201.01-2	18	40.314.007.01-2	43.601.107.01-2	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8					
		43.624.201.01-2	24									
		43.632.201.01-2	32					49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)		6		
						49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8					
										SCREWDRIVER	43.601.107.01-2	

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1.8 mm

R	31.322.102.01-2	38°	18°
NR	31.312.102.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s

1 mm CH-5mm CH-7mm CH-9mm 1.8 mm CH-5mm CH-7mm CH-9mm 3 mm CH-5mm CH-7mm CH-9mm

R	31.322.102.29-2	30°	25°	20°	31.322.102.21-2	25°	15°	10°	31.322.102.23-2	20°	20°	15°
NR	31.312.102.29-2				31.312.102.21-2				31.312.102.23-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}		SCANBODY	PEEK PINS	TYPE
		50.312.102.03-2 (IG-3mm)	43.621.415.01-2	34.612.102.01-2	33.345.856.01-2	3	25°	23.412.102.01-2	54.315.102.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
52.410.128.01-2	10				33.445.856.01-2	4				49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
					33.635.856.01-2	6				49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
										49.416.000.04-2 (13 mm)	D
										49.416.000.05-2 (13 mm)	E

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.317.065.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.317.005.02-2	43.601.105.01-2

LAB SCANBODY
30.412.001.01-2

SCREWDRIVER	43.625.105.01-2
-------------	-----------------

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm 2 mm 3 mm 4 mm

R	42.302.102.01-2	42.302.102.02-2	42.302.102.03-2	42.302.102.04-2
---	-----------------	-----------------	-----------------	-----------------

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT

2.5/3.9 mm

NR	48.302.102.02-2
----	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
1,2 mm			
R	31.322.109.01-2	45°	29°
NR	31.312.109.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.109.01-2	43.621.415.01-2	34.612.109.01-2	33.360.754.01-2'	3	25°
52.410.128.01-2	10				33.460.754.01-2'	4	
					33.660.754.01-2'	6	

*Only for R

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.070.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.314.005.02-2	43.601.105.01-2

ANALOG
22.612.109.01-2
LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,4 mm

R	31.320.110.01-2	45°	30°
NR	31.310.110.01-2		

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.110.04-2 IG-3mm	43.621.410.01-2 43.624.410.01-2	34.610.110.01-2	33.360.756.01-2*	3	20°
52.410.117.01-2	10				33.460.756.01-2*	4	
52.412.117.01-2	12				33.660.756.01-2*	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6	LAB SCANBODY
41.318.083.01-2	43.618.201.01-2	18	40.318.007.01-2	43.601.107.01-2	30.410.006.01-2
	43.624.201.01-2	24			
	43.632.201.01-2	32			

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s
0,4 mm				2,5 mm		
R	31.323.111.01-2	45°	30°			25°
NR	31.313.111.01-2			31.313.111.03-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.110.04-2 IG-3mm	43.621.410.01-2 43.624.410.01-2	34.610.110.01-2	33.360.756.01-2	3	20°
52.410.117.01-2	10				33.460.756.01-2	4	
52.412.117.01-2	12				33.660.756.01-2	6	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.083.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER TORX T6
40.318.007.01-2	43.601.107.01-2

LAB SCANBODY
30.413.002.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT 1 mm	α_s	GINGIVAL HEIGHT 2 mm	α_s	α_s
R	31.321.119.01-2	35°	31.321.119.02-2	30°	30°
NR	31.311.119.01-2		31.311.119.02-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.311.119.03-2 IG-3mm	43.621.415.01-2	34.611.119.01-2	33.360.756.01-2	3	25
52.410.132.01-2	10				33.460.756.01-2	4	
52.412.132.01-2	12				33.660.756.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1,27	LAB SCANBODY
41.316.080.01-2	43.618.201.01-2	18	40.316.005.07-2	43.601.105.01-2	30.410.006.01-2
	43.624.201.01-2	24			
	43.632.201.01-2	32			

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.301.119.01-2	42.301.119.02-2	42.301.119.03-2	42.301.119.04-2

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT

3.5/4.9 mm

NR	48.312.119.03-2
----	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_s

1 mm

R	31.323.121.01-2	45°	25°
NR	31.313.121.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s

2 mm CH-5mm CH-7mm CH-9mm 3 mm CH-5mm CH-7mm CH-9mm

NR	31.313.121.22-2	25°	20°	15°	31.313.121.23-2	25°	20°	10°
----	-----------------	-----	-----	-----	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.108.01-2	10	50.312.120.03-2 IG-3mm	43.621.410.01-2 43.624.410.01-2	34.612.120.01-2	33.360.754.01-2	3	20°
52.412.108.01-2	12				33.460.754.01-2	4	
					33.660.754.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
41.316.080.01-2	43.618.201.01-2	18	40.316.005.07-2	43.601.105.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm 2 mm 3 mm 4 mm

R	42.303.121.01-2	42.303.121.02-2	42.303.121.03-2	42.303.121.04-2
---	-----------------	-----------------	-----------------	-----------------

ANGULATED MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT

2.5 mm 3.5 mm

NR	48.302.120.02-2	48.312.120.03-2
----	-----------------	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1 mm

R	31.323.121.01-2	45°	25°
NR	31.313.121.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s GINGIVAL HEIGHT α_s α_s α_s

2 mm CH-5mm CH-7mm CH-9mm 3 mm CH-5mm CH-7mm CH-9mm

NR	31.313.121.22-2	25°	20°	15°	31.313.121.23-2	25°	20°	10°
----	-----------------	-----	-----	-----	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.313.121.01-2	43.621.410.01-2	34.613.121.01-2	33.360.754.01-2	3	20°
52.410.109.01-2	10	50.313.121.03-2 IG=3mm	43.624.410.01-2		33.460.754.01-2	4	
52.412.109.01-2	12				33.660.754.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
41.316.080.01-2	43.618.201.01-2	18	40.316.005.07-2	43.601.105.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

MULTI-UNIT

GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT GINGIVAL HEIGHT

1 mm 2 mm 3 mm 4 mm

R	42.303.121.01-2	42.303.121.02-2	42.303.121.03-2	42.303.121.04-2
---	-----------------	-----------------	-----------------	-----------------

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
1.4 mm			
R	31.324.124.01-2	42°	19°
NR	31.314.124.01-2		

DYNAMIC SCANBODY (LAB/CLIN)			DIGITAL ANALOG
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR
		50.314.124.01-2	43.621.410.01-2 43.624.410.01-2
52.410.125.01-2	10		
			34.614.124.01-2

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.075.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL		
DYNAMIC MILLING TOOL	SHANK	α_{di}
33.335.758.01-2	3	30°
33.435.758.01-2	4	
33.635.758.01-2	6	

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER UNIGRIP
40.320.008.02-2	43.625.108.01-2

LAB SCANBODY
30.412.003.01-2

SCANBODY OP		
SCANBODY	PEEK PINS	TYPE
54.315.124.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.625.108.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		α_c	GINGIVAL HEIGHT α_s		α_c
1.1 mm			2 mm		
R	31.323.125.01-2	42°	20°	31.323.125.02-2	25°
NR	31.313.125.01-2			31.313.125.02-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s		α_s	α_s	GINGIVAL HEIGHT α_s		α_s	α_s	α_s
1.1 mm		CH-5mm	CH-7mm	CH-9mm	3 mm		CH-5mm	CH-7mm
R	31.323.125.21-2	30°	25°	15°		20°	15°	10°
NR	31.313.125.21-2				31.313.125.23-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREW/DRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.313.125.01-2	43.621.410.01-2 43.624.410.01-2	34.613.125.01-2	33.315.804.01-2	3	30°
52.410.117.01-2	10	50.313.125.03-2			33.415.804.01-2	4	
52.412.117.01-2	12	IG-3mm			33.615.804.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6
41.316.078.01-2	43.618.201.01-2	18	40.316.007.01-2	43.601.107.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.413.002.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.125.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.103.02-2

MULTI-UNIT

GINGIVAL HEIGHT		GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
1 mm		2 mm	3 mm	4 mm
R	42.303.125.01-2	42.303.125.02-2	42.303.125.03-2	42.303.125.04-2

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT		GINGIVAL HEIGHT	GINGIVAL HEIGHT
1.5 mm		2.5 mm	3.5 mm
NR	48.312.125.01-2	48.312.125.02-2	48.312.125.03-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c			
2.5 mm			
NR	31.322.128.01-2	45°	30°

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.044.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.320.003.05-2	43.601.103.01-2

LAB SCANBODY
30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,3 mm

R	31.325.129.01-2	35°	30°
NR	31.315.129.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.315.129.01-2	43.621.410.01-2 43.624.410.01-2	34.615.129.01-2	33.390.958.01-2	3	30°
52.410.130.01-2	10				33.490.958.01-2	4	
					33.690.958.01-2	6	

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	ANALOG
41.320.090.01-2	43.618.201.01-2	18	40.320.008.03-2	43.625.108.01-2	22.615.129.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.415.007.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0.5 mm			
R	31.322.130.01-2	30°	29°
NR	31.312.130.01-2		

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.081.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.08-2	43.601.105.01-2

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0.5 mm

R	31.323.131.01-2	45°	29°
NR	31.313.131.01-2		

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0.5 mm CH-5mm CH-7mm CH-9mm

R	31.323.131.21-2	30°	20°	15°
NR	31.313.131.21-2			

DYNAMIC SCREWS

DYNAMIC SCREW DYNAMIC SCREWDRIVER SCREWDRIVER LENGTH (mm)

41.316.081.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS

STRAIGHT SCREW SCREWDRIVER Hex. 1.27

40.316.005.08-2	43.601.105.01-2
-----------------	-----------------

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0.5 mm			
R	31.324.132.01-2	45°	28°
NR	31.314.132.01-2		

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.081.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.316.005.08-2	43.601.105.01-2

LAB SCANBODY
30.414.003.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1 mm

R	31.320.135.01-2	45°	30°
NR	31.320.135.01-2		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.080.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.315.804.01-2	3	25°
33.415.804.01-2	4	
33.415.804.01-2	6	

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.314.007.02-2	43.601.107.01-2

LAB SCANBODY

30.410.006.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.135.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.107.01-2

STANDARD DYNAMIC TIBASE

	GINGIVAL HEIGHT 0.7 mm	α_s	α_c	GINGIVAL HEIGHT 1.5 mm	α_s	α_c	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm	α_s	GINGIVAL HEIGHT 5 mm	α_s
R	31.320.136.01-2	45°	30°	31.320.136.02-2	25°	20°	31.320.136.04-2	31.320.136.05-2	20°	31.320.136.06-2	15°
NR	31.310.136.01-2			31.310.136.02-2			31.310.136.04-2	31.310.136.05-2		31.310.136.06-2	

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.136.01-2	43.621.415.01-2	34.610.136.01-2	33.360.754.01-2	3	25°
52.410.128.01-2	10	50.310.136.04-2			33.460.754.01-2	4	
		IG=3mm			33.660.754.01-2	6	

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.25	LAB SCANBODY
41.316.071.01-2	43.618.201.01-2	18	40.316.004.03-2	43.601.104.01-2	30.410.006.01-2
	43.624.201.01-2	24			
	43.632.201.01-2	32			

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.300.136.01-2	42.300.136.02-2	42.300.136.3-2	42.300.136.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
 0,6 mm

R	31.324.137.01-2	45°	30°
NR	31.314.137.01-2		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.044.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER TORX T6	
40.320.007.04-2	43.601.107.01-2	
		LAB SCANBODY
		30.414.008.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s
1.5 mm		45°
R	31.320.145.01-2	
NR	31.310.145.01-2	

DYNAMIC SCANBODY
(LAB/CLIN)

SCANBODY	HEIGHT mm
52.410.128.01-2	10

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.315.078.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS

STRAIGHT SCREW	SCREWDRIVER
40.314.005.01-2	43.601.105.01-2

LAB SCANBODY
30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1.5 mm

R	31.323.149.01-2	45°	29°
NR	31.313.149.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.161.01-2	43.621.415.01-2	34.610.161.01-2	33.320.704.01-2*	3	25°
52.410.132.01-2	10				33.420.704.01-2*	4	
52.412.132.01-2	12				33.620.704.01-2*	6	

*Only for R

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.079.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT
SCREW

40.316.014.01-2

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		
0.5 mm		
R	31.323.150.01-2	45°

D.S. (L./C.)

SCREWDRIVER
ADAPTOR

43.620.411.01-2

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.046.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.314.004.04-2	43.601.104.01-2

LAB SCANBODY

30.413.005.01-2

STANDARD DYNAMIC TIBASE**GINGIVAL
HEIGHT** α_s α_c

0.3 mm

R	31.323.151.01-2	45°	30°
---	-----------------	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.409.123.01-2	9	50.313.151.01-2	43.621.410.01-2 43.624.410.01-2	34.613.151.01-2	33.390.716.01-2	3	30°
					33.490.716.01-2	4	
					33.690.716.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP
41.314.039.01-2	43.618.201.01-2	18	40.314.008.01-2	43.625.108.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

DYNAMIC SCREWS

DYNAMIC SCREW

41.317.040.01-2

REFERENCE SCANBODY		
SCANBODY	PEEK PINS	TYPE
54.322.020.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
SCREWDRIVER		43.625.105.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1.5 mm

R	31.320.159.01-2	41°	17°
NR	31.310.159.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.159.01-2	43.621.415.01-2	34.610.159.01-2	33.335.754.01-2*	3	25°
52.410.128.01-2	10				33.435.754.01-2*	4	
					33.635.754.01-2*	6	

*Only for R

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER UNIGRIP	ANALOG
41.314.067.02-2	43.618.201.01-2	18	40.314.008.02-2	43.625.108.01-2	22.610.159.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.410.006.01-2

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.300.159.01-2	42.300.159.02-2	42.300.159.03-2	42.300.159.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0.5 mm

R	31.320.160.01-2	45°	30°
NR	31.310.160.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}	
52.410.131.01-2	10	50.310.160.01-2	43.621.415.01-2	34.610.160.01-2	33.315.804.01-2	3	25°	23.410.160.01-2
52.412.131.01-2	12				33.415.804.01-2	4		
					33.615.804.01-2	6		

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER TORX T6	ANALOG
41.316.078.01-2	43.618.201.01-2	18	40.316.007.01-2	43.601.107.01-2	22.610.160.01-2
	43.624.201.01-2	24			LAB SCANBODY
	43.632.201.01-2	32			30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1.5 mm

R	31.320.161.01-2	45°	25°
NR	31.310.161.01-2		

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.161.01-2	43.621.415.01-2	34.610.161.01-2	33.320.704.01-2*	3	25°
52.410.132.01-2	10				33.420.704.01-2*	4	
52.412.132.01-2	12				33.620.704.01-2*	6	

*Only for R

DYNAMIC SCREWS

STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.079.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT
SCREW

30.410.006.01-2

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
1.5 mm			
R	31.324.162.01-2	45°	24°
NR	31.314.162.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.161.01-2	43.621.415.01-2	34.610.161.01-2	33.320.704.01-2*	3	25°
52.410.132.01-2	10				33.420.704.01-2*	4	
52.412.132.01-2	12				33.620.704.01-2*	6	

*Only for R

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.079.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS
40.316.014.01-2

LAB SCANBODY
30.414.003.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,5 mm

R	31.323.163.01-2	45°	30°
---	-----------------	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.408.112.01-2	8	50.313.163.01-2	43.620.411.01-2	34.613.163.01-2	33.390.716.01-2	3	30°
					33.490.716.01-2	4	
					33.690.716.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.039.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

**STRAIGHT
SCREW**

40.314.014.01-2

LAB SCANBODY

30.413.005.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
0,9 mm			
R	31.320.166.01-2	45°	30°
NR	31.310.166.01-2		

DYNAMIC SCANBODY (LAB/CLIN)				DIGITAL ANALOG	DYNAMIC MILLING TOOL		
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.310.166.03-2 IG-3mm	43.621.415.01-2	34.610.166.01-2	33.330.734.01-2	3	20°
52.410.128.01-2	10				33.430.734.01-2	4	
					33.630.734.01-2	6	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.084.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
40.314.004.02-2	43.601.104.01-2

LAB SCANBODY
30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
0,9 mm

R	31.322.167.01-2	43°	30°
NR	31.312.167.01-2		

DYNAMIC SCANBODY (LAB/CLIN)**DIGITAL ANALOG****DYNAMIC MILLING TOOL**

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.313.167.03-2 (IG= 3mm)	43.621.410.01-2 43.624.410.01-2	34.613.167.01-2	33.330.734.01-2	3	20°
52.410.117.01-2	10				33.430.734.01-2	4	
52.412.117.01-2	12				33.630.734.01-2	6	

DYNAMIC SCREWS**STRAIGHT SCREWS**

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)	STRAIGHT SCREW	SCREWDRIVER Hex. 1.25
41.316.084.02-2	43.618.201.01-2	18	40.316.004.01-2	43.601.104.01-2
	43.624.201.01-2	24		
	43.632.201.01-2	32		

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c			GINGIVAL HEIGHT α_s α_c			GINGIVAL HEIGHT α_s		
0,6 mm			1,5 mm			3 mm		
R	31.322.169.01-2	45°	29°	31.322.169.02-2	25°	15°	31.322.169.04-2	20°
NR	31.312.169.01-2			31.312.169.02-2			31.312.169.04-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s		α_s α_s		GINGIVAL HEIGHT α_s α_s		α_s α_s		α_s	
1,5 mm		CH-5mm	CH-7mm	CH-9mm	3 mm		CH-5mm	CH-7mm	CH-9mm
R	31.322.169.22-2	30°	25°	15°	31.322.169.24-2	25°	20°	15°	
NR	31.312.169.22-2				31.312.169.24-2				

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.169.01-2	43.621.410.01-2	34.612.169.01-2	33.330.734.01-2	3	25°
52.410.117.01-2	10	50.312.169.04-2 IG=3mm	43.624.410.01-2		33.430.734.01-2	4	
52.412.117.01-2	12				33.630.734.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.317.070.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.412.001.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.169.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.104.01-2

MULTI-UNIT

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1 mm		2 mm		3 mm		4 mm	
R	42.302.169.01-2	42.302.169.02-2	42.302.169.03-2	42.302.169.04-2			

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1,5/2,9 mm		2,5/3,9 mm		4,5/5,9 mm	
NR	48.312.169.01-2	48.312.169.02-2	48.312.169.03-2		

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s
0.3 mm

R	31.322.170.01-2	38°
NR	31.312.170.01-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s
0.3 mm CH-5mm CH-7mm CH-9mm

R	31.322.170.21-2	30°	20°	15°
NR	31.312.170.21-2			

DYNAMIC SCANBODY (L/C)

SCANBODY HEIGHT
mm

52.410.104.01-2	10
52.412.104.01-2	12

DYNAMIC SCREWS

DYNAMIC SCREW DYNAMIC SCREWDRIVER SCREWDRIVER LENGTH (mm)

41.316.079.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL

DYNAMIC MILLING TOOL SHANK α_{di}

33.390.754.01-2	3	25°
33.490.754.01-2	4	
33.690.754.01-2	6	

SCANBODY OP

SCANBODY PEEK PINS TYPE

54.315.170.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.103.02-2

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s
0.3 mm		
R	31.323.171.01-2	35°
NR	31.313.171.01-2	

DYNAMIC SCANBODY (L/C)

SCANBODY	HEIGHT mm
52.410.117.01-2	10
52.412.117.01-2	12

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.079.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.754.01-2	3	25°
33.490.754.01-2	4	
33.690.754.01-2	6	

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c
 0.5 mm

R	31.323.173.01-2	45°	30°
---	-----------------	-----	-----

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.040.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.413.005.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		
0.5 mm		
NR	31.310.176.01-2	35°

DYNAMIC SCANBODY (LAB/CLIN)			DIGITAL ANALOG
SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR
52.408.138.01-2	8	50.310.176.01-2	43.621.415.01-2
			34.610.176.01-2

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.044.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL		
DYNAMIC MILLING TOOL	SHANK	α_{di}
33.360.756.01-2	3	35°
33.460.756.01-2	4	
33.660.756.01-2	6	

LAB SCANBODY
30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

1.5 mm

R	31.320.178.01-2	45°
NR	31.310.178.01-2	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.080.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		
0,4 mm		
R	31.322.181.01-2	45°

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s				
0,4 mm CH-5mm CH-7mm CH-9mm				
R	31.322.181.21-2	30°	30°	20°

DYNAMIC SCANBODY (LAB/CLIN)

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR
52.408.112.01-2	8	50.312.181.01-2	43.620.411.01-2

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.043.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DYNAMIC MILLING TOOL

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.360.756.01-2	3	30°
33.460.756.01-2	4	
33.660.756.01-2	6	

SCANALOG

23.412.181.01-2

LAB SCANBODY
30.413.005.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0,5 mm

R	31.322.183.01-2	45°
---	-----------------	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

REFERENCE SCANBODY

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.408.136.01-2	8	50.312.183.01-2	43.620.411.01-2	34.612.183.01-2

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.330.734.01-2	3	30°
33.430.734.01-2	4	
33.630.734.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.048.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

SCANBODY	PEEK PINS	TYPE
54.322.183.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c			GINGIVAL HEIGHT α_s α_c			GINGIVAL HEIGHT α_s		
1.2 mm			2.5 mm			3.5 mm		
R	31.323.186.01-2	40°	30°	31.323.186.02-2	20°	18°	31.323.186.03-2	15°
NR	31.313.186.01-2			31.313.186.02-2			31.313.186.03-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s			GINGIVAL HEIGHT α_s α_s α_s					
1.2 mm CH-5mm CH-7mm CH-9mm			2.5/3.9 mm CH-5mm CH-7mm CH-9mm					
R	31.323.186.21-2	30°	25°	15°	31.323.186.22-2	30°	25°	15°
NR	31.313.186.21-2				31.313.186.22-2			

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.408.101.01-2	8	50.313.186.04-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.613.186.01-2	33.330.734.01-2	3	25°
52.410.101.01-2	10				33.430.734.01-2	4	
52.412.101.01-2	12				33.630.734.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.084.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

MULTI-UNIT

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1 mm		2 mm		3 mm		4 mm		5 mm	
R	42.303.186.01-2	42.303.186.02-2	42.303.186.03-2	42.303.186.04-2	42.303.186.05-2				

ANGULATED MULTI-UNIT 20°

GINGIVAL HEIGHT		GINGIVAL HEIGHT	
2.5/3.9 mm		3.5/4.9 mm	
NR	48.312.186.02-2	48.312.186.03-2	

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.186.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s
0.3 mm				0.5 mm			1 mm		
R	31.322.009.01-2	45°	25°	31.322.009.02-2	25°	31.322.009.03-2	25°		25°
NR	31.312.009.01-2			31.312.009.02-2		31.312.009.03-2			

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
		50.312.187.01-2	43.621.410.01-2 43.624.410.01-2	34.612.187.01-2	33.390.716.01-2	3	25°
52.410.114.01-2	10				33.490.716.01-2	4	
52.412.114.01-2	12				33.690.716.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.059.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s
1.5 mm		45°
R	31.320.188.01-2	
NR	31.310.188.01-2	

DYNAMIC SCANBODY
(LAB/CLIN)

SCANBODY	HEIGHT mm
52.410.128.01-2	10

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.315.078.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY
30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

1.8 mm

R	31.320.190.01-2	45°
NR	31.310.190.01-2	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.084.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s
0,8 mm		
R	31.322.191.01-2	45°
NR	31.312.191.01-2	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.084.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER
40.316.003.01-2	43.601.103.02-2

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0,3 mm

R	31.323.192.01-2	45°
---	-----------------	-----

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.048.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.413.005.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		
0.3 mm		
R	31.323.193.01-2	45°

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.051.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER
40.318.005.01-2	43.601.105.01-2

LAB SCANBODY
30.413.005.01-2

REFERENCE SCANBODY		
SCANBODY	PEEK PINS	TYPE
54.322.193.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
	SCREWDRIVER	43.625.105.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0.5 mm

R	31.323.195.01-2	45°
---	-----------------	-----

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0.5 mm

CH=5mm

CH=7mm

CH=9mm

R	31.323.195.21-2	30°	25°	20°
---	-----------------	-----	-----	-----

DYNAMIC MILLING TOOL

DYNAMIC MILLING TOOL SHANK α_{di}

33.390.716.01-2	3	25
33.490.716.01-2	4	
33.690.716.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW

DYNAMIC SCREW Ø2,6

DYNAMIC SCREWDRIVER

SCREWDRIVER LENGTH (mm)

41.317.041.01-2	41.317.052.36-2	43.618.201.01-2	18
		43.624.201.01-2	24
		43.632.201.01-2	32

LAB SCANBODY

30.413.005.01-2

REFERENCE SCANBODY

SCANBODY

PEEK PINS

TYPE

54.322.195.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
	SCREWDRIVER	
	43.601.104.01-2	

STANDARD DYNAMIC TIBASE								
GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s
1.2 mm			2 mm			3 mm		
R	31.320.196.01-2	40°	31.320.196.02-2	25°		31.320.196.03-2	25°	
NR	31.310.196.01-2		31.310.196.02-2			31.310.196.03-2		

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.086.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY
30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s
1,2 mm				2 mm				3 mm		
R	31.322.197.01-2	35°	25°	31.322.197.02-2	20°	25°		31.322.197.03-2	20°	
NR	31.312.197.01-2			31.312.197.02-2				31.312.197.03-2		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.086.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s
1.2 mm		40°
R	31.324.198.01-2	
NR	31.314.198.01-2	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.086.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY
30.414.003.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0,3 mm

R	31.322.205.01-2	45°
---	-----------------	-----

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.408.112.01-2	8	50.312.205.01-2	43.620.411.01-2	34.612.205.01-2	33.390.716.01-2	3	30°
					33.490.716.01-2	4	
					33.690.716.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.317.040.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		α_c	GINGIVAL HEIGHT α_s		α_c	GINGIVAL HEIGHT α_s	
1.5 mm			3 mm			4 mm	
R	31.320.207.01-2	35°	15°	31.320.207.03-2	20°	31.320.207.04-2	15°
NR	31.310.207.01-2			31.310.207.03-2		31.310.207.04-2	

DYNAMIC 3TIBASE

GINGIVAL HEIGHT		α_s	α_s	α_s	GINGIVAL HEIGHT		α_s	α_s	α_s
1.5 mm		CH-5mm	CH-7mm	CH-9mm	3 mm		CH-5mm	CH-7mm	CH-9mm
R	31.320.207.21-2	25°	20°	15°	31.320.207.23-2	20°	15°	10°	
NR	31.310.207.21-2				31.310.207.23-2				

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.103.01-2	10	50.310.207.03-2 IG-3mm	43.621.410.01-2 43.624.410.01-2	34.610.207.01-2	33.345.856.01-2*	3	30°
52.412.103.01-2	12				33.445.856.01-2*	4	
					33.645.856.01-2*	6	

*Only for R

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.066.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.410.006.01-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.207.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.107.01-2

MULTI-UNIT

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
1 mm		2 mm		3 mm		4 mm	
R	42.302.207.01-2	42.302.207.02-2	42.302.207.03-2	42.302.207.04-2			

ANGULATED MULTI-UNIT

GINGIVAL HEIGHT		GINGIVAL HEIGHT		GINGIVAL HEIGHT	
2.5/3.9mm		3.5/4.9 mm		4.5/5.9 mm	
NR	48.312.207.02-2	48.312.207.03-2	48.312.207.04-2		

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

1.5 mm

R	31.324.208.01-2	45°	20°
NR	31.314.208.01-2		

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL
SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.103.01-2	10	50.310.207.03-2 IG-3mm	43.621.410.01-2 43.624.410.01-2	34.614.208.01-2
52.412.103.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.345.856.01-2*	3	30°
33.445.856.01-2*	4	
33.645.856.01-2*	6	

*Only for **R**

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.066.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.414.003.01-2

SCANBODY	PEEK PINS	TYPE
54.315.207.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.107.01-2

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.302.207.01-2	42.302.207.02-2	42.302.207.03-2	42.302.207.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
1 mm			2 mm				3 mm			
R	31.323.223.01-2	30°	31.323.223.02-2		25°	15°	31.323.223.03-2		20°	10°
NR	31.323.223.01-2		31.323.223.02-2				31.313.223.03-2			

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.078.02-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

DIGITAL ANALOG

DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
34.613.223.01-2	33.315.604.01-2	3	25°
	33.415.604.01-2	4	
	33.615.604.01-2	6	

DYNAMIC MILLING TOOL

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.223.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.103.02-2

LAB SCANBODY

30.413.002.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0.5 mm

R	31.320.229.01-2	40°
NR	31.310.229.01-2	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.064.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c
1.8 mm			
R	31.322.236.01-2	20°	25°
NR	31.312.236.01-2		

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.075.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.236.21-2	49.414.000.01-2 (6 mm)	A
	49.415.000.01-2 (9 mm)	
	49.416.000.01-2 (13 mm)	
	49.414.000.02-2 (6 mm)	B
	49.415.000.02-2 (9 mm)	
	49.416.000.02-2 (13 mm)	
	49.414.000.03-2 (6 mm)	C
	49.415.000.03-2 (9 mm)	
	49.416.000.03-2 (13 mm)	
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.103.02-2

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0,5 mm

R	31.323.242.01-2	40°
---	-----------------	-----

DIGITAL ANALOG

DIGITAL
ANALOG

34.613.242.01-2

DYNAMIC MILLING TOOL

DYNAMIC
MILLING TOOL SHANK α_{di}

33.370.716.01-2

3

33.470.716.01-2

4

33.670.716.01-2

6

25°

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREW Ø2,6	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.318.045.01-2	41.318.055.01-2 (Direct MU)	43.618.201.01-2	18
		43.624.201.01-2	24
		43.632.201.01-2	32

LAB SCANBODY

30.413.005.01-2

REFERENCE SCANBODY

SCANBODY	PEEK PINS	TYPE
54.322.242.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3,8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
	SCREWDRIVER	43.601.104.01-2

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.084.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

MULTI-UNIT		
GINGIVAL HEIGHT		GINGIVAL HEIGHT
2 mm		3 mm
R	42.303.243.02-2	42.303.243.03-2

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.316.003.01-2	43.601.103.02-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s
0.5 mm

R	31.323.245.01-2	40°
---	-----------------	-----

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.039.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.413.005.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s
1 mm			2 mm		
R	31.322.246.01-2	25°	31.322.246.02-2	25°	
NR	31.312.246.01-2		31.312.246.02-2		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.080.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s
0.3 mm

R	31.323.247.01-2	40°
---	-----------------	-----

DYNAMIC SCREWS

SCREWDRIVER LENGTH (mm)	SCREWDRIVER LENGTH (mm)
43.618.201.01-2	18
43.624.201.01-2	24
43.632.201.01-2	32

LAB SCANBODY

30.413.005.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s
1 mm

R	31.320.249.01-2	40°
NR	31.310.249.01-2	

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG	DYNAMIC MILLING TOOL	SHANK	α_{di}
52.410.128.01-2	10	50.310.249.03-2 IG=3mm	43.621.415.01-2	34.610.249.01-2	33.315.804.01-2	3	25°
					33.415.804.01-2	4	
					33.615.804.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.080.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c	GINGIVAL HEIGHT		α_s	α_c
1mm				1.5 mm				3 mm			
R	31.322.251.01-2	40°	25°	31.322.251.02-2	40°	20°	31.322.251.04-2	25°	10°		
NR	31.312.251.01-2			31.312.251.02-2			31.312.251.04-2				

DYNAMIC SCREWS

DYNAMIC SCREW	HIGH DYNAMIC SCREW	DYNAMIC SCREW/DRIVER	SCREW/DRIVER LENGTH (mm)
41.316.064.02-2	40.316.007.01-2	43.618.201.01-2	18
		43.624.201.01-2	24
		43.632.201.01-2	32

LAB SCANBODY

30.412.001.01-2

DYNAMIC SCREWS	
DYNAMIC SCREW	HIGH DYNAMIC SCREW
41.318.065.01-2	40.318.003.01-2

STRAIGHT SCREWS	
STRAIGHT SCREW	SCREWDRIVER Hex. 1.20
40.318.003.01-2	43.601.103.02-2

MULTI-UNIT				
	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.302.257.01-2	42.302.257.02-2	42.302.257.03-2	42.302.257.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s α_c

0,3 mm

R	31.323.025.01-2	40°	30°
---	-----------------	-----	-----

DYNAMIC 3TIBASE

GINGIVAL HEIGHT α_s α_s α_s

0,3 mm

CH-5mm

CH-7mm

CH-9mm

R	31.323.025.21-2	30°	25°	10°
---	-----------------	-----	-----	-----

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL
SCANALOG
REFERENCE SCANBODY

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.408.112.01-2	8	50.313.025.02-2	43.621.410.01-2	34.613.025.01-2
52.410.111.01-2	10		43.624.410.01-2	
			43.620.411.01-2	

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.390.716.01-2	3	30
33.490.716.01-2	4	
33.690.716.01-2	6	

23.413.025.01-2

DYNAMIC SCREWS
STRAIGHT SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.314.039.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

STRAIGHT SCREW	SCREWDRIVER Hex. 1.27
40.314.005.04-2	43.601.105.01-2

ANALOG

22.613.025.01-2
LAB SCANBODY
30.413.005.01-2

SCANBODY	PEEK PINS	TYPE
54.322.025.31-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
	CAPS	mm
	49.418.000.01-2 (Regular) 49.418.000.02-2 (Wide)	3.8
	49.419.000.01-2 (Regular) 49.419.000.02-2 (Wide)	6
	49.420.000.01-2 (Regular) 49.420.000.02-2 (Wide)	8
	SCREWDRIVER	43.625.108.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0.5 mm

R	31.320.260.01-2	40°
NR	31.310.260.01-2	

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.103.01-2	10	50.310.207.03-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.610.260.01-2
52.412.103.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.345.856.01-2	3	30°
33.445.856.01-2	4	
33.645.856.01-2	6	

23.410.260.01-2

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.066.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.410.006.01-2

SCANBODY	PEEK PINS	TYPE
54.315.207.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.107.01-2

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.302.207.01-2	42.302.207.02-2	42.302.207.03-2	42.302.207.04-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0.5 mm

R	31.323.261.01-2	40°
NR	31.313.261.01-2	

DYNAMIC SCANBODY (LAB/CLIN)

DIGITAL ANALOG

DYNAMIC MILLING TOOL

SCANALOG

SCANBODY OP

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.103.01-2	10	50.310.207.03-2 (IG-3mm)	43.621.410.01-2 43.624.410.01-2	34.613.261.01-2
52.412.103.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.345.856.01-2	3	25°
33.445.856.01-2	4	
33.645.856.01-2	6	

23.413.261.01-2

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.066.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.413.004.01-2

SCANBODY	PEEK PINS	TYPE
54.315.207.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.107.01-2

MULTI-UNIT

GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT	GINGIVAL HEIGHT
1 mm	2 mm	3 mm	4 mm
R	42.302.207.01-2	42.302.207.02-2	42.302.207.03-2
			42.302.207.04-2

STANDARD DYNAMIC TIBASE
GINGIVAL HEIGHT α_s

0.5 mm

R	31.324.262.01-2	40°
NR	31.314.262.01-2	

DYNAMIC SCANBODY (LAB/CLIN)
DIGITAL ANALOG
DYNAMIC MILLING TOOL

SCANBODY	HEIGHT mm	ADAPTOR	SCREWDRIVER ADAPTOR	DIGITAL ANALOG
52.410.103.01-2	10	50.310.207.03-2 (IG=3mm)	43.621.410.01-2 43.624.410.01-2	34.614.262.01-2
52.412.103.01-2	12			

DYNAMIC MILLING TOOL	SHANK	α_{di}
33.345.856.01-2	3	25°
33.445.856.01-2	4	
33.445.856.01-2	6	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.316.066.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

MULTI-UNIT

	GINGIVAL HEIGHT 1 mm	GINGIVAL HEIGHT 2 mm	GINGIVAL HEIGHT 3 mm	GINGIVAL HEIGHT 4 mm
R	42.302.207.01-2	42.302.207.02-2	42.302.207.03-2	42.302.207.04-2

SCANBODY OP

SCANBODY	PEEK PINS	TYPE
54.315.207.21-2	49.414.000.01-2 (6 mm) 49.415.000.01-2 (9 mm) 49.416.000.01-2 (13 mm)	A
	49.414.000.02-2 (6 mm) 49.415.000.02-2 (9 mm) 49.416.000.02-2 (13 mm)	B
	49.414.000.03-2 (6 mm) 49.415.000.03-2 (9 mm) 49.416.000.03-2 (13 mm)	C
	49.416.000.04-2 (13 mm)	D
	49.416.000.05-2 (13 mm)	E
SCREWDRIVER		43.601.107.01-2

LAB SCANBODY

30.414.008.01-2

DYNAMIC 3TIBASE			
GINGIVAL HEIGHT		α_s	α_s
0.3 mm		CH= 7mm	CH= 9mm
R	31.323.264.21-2	25°	25°

MINI SCANBODY

53.413.264.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s		GINGIVAL HEIGHT α_s		
1.5 mm		2 mm		
R	31.322.265.03-2	35°	31.322.265.04-2	25°
NR	31.312.265.03-2		31.312.265.04-2	

DYNAMIC SCREWS

DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
43.618.201.01-2	18
43.624.201.01-2	24
43.632.201.01-2	32

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s	GINGIVAL HEIGHT		α_s
0.5 mm			2 mm		
R	31.320.266.01-2	35°	31.310.266.04-2	25°	
NR	31.310.266.01-2		31.320.266.04-2		

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.068.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.410.006.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT		α_s
0.5 mm		
R	31.322.267.01-2	40°
NR	31.312.267.01-2	

DYNAMIC SCREWS		
DYNAMIC SCREW	DYNAMIC SCREWDRIVER	SCREWDRIVER LENGTH (mm)
41.320.068.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY
30.412.001.01-2

STANDARD DYNAMIC TIBASE

GINGIVAL HEIGHT α_s

0,5 mm


R	31.323.268.01-2	45°
NR	31.313.268.01-2	

DYNAMIC SCREWS

DYNAMIC SCREW	DYNAMIC SCREW/DRIVER	SCREW/DRIVER LENGTH (mm)
41.320.068.01-2	43.618.201.01-2	18
	43.624.201.01-2	24
	43.632.201.01-2	32

LAB SCANBODY

30.413.002.01-2





This printed version is probably
outdated and may not include all
products or compatibilities.


**Access and download the latest
version of this catalogue**

www.dynamicabutment.com

DAS LIBRARIES









DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS



MARCA	MODELO	PLATAFORMA	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
-------	--------	------------	-----------	--------	--------	-----------------	------------------	-------------	-------------	-------------	------------------	------	-------	-----------------	------------

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
AB	MULTI-UNIT D1-P14	CONICAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
AB	MULTI-UNIT D1-P64	CONICAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
AB	l22	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
AB	l5	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
AB	l55	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
AB	l10	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
AB	l15	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
AB	l2	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ACE	MULTI-UNIT	CONICAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
ACE	INFINITY EXTERNAL HEX	NP 3.5	41.316.059.01-2	1.60	20 Ncm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
ACE	INFINITY INTERNAL HEX	4.50	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ACE	INFINITY EXTERNAL HEX	WP 5	41.320.047.01-2	2.00	25 Ncm	4.70	3.00	3.30	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
ACE	INFINITY EXTERNAL HEX	RP 4.1	41.320.060.01-2	2.00	25 Ncm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
ACE	INFINITY OCTAGON	WP 6.5	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
ACE	INFINITY TRI-CAM	4.30	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
ACE	INFINITY TRI-CAM	5.00	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
ACE	INFINITY INTERNAL HEX	3.50	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ACE	INFINITY TRI-CAM	3.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
ACE	INFINITY OCTAGON	RP 4.8	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
ADIN	MULTI-UNIT	CONICAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
ADIN	TOUAREG CLOSEFIT	UNP	41.315.078.01-2	No-80	15 Ncm	7.80	2.45	5.77	6.00	1.80	2.30	conical	65°	45° chamfer	Hex. 1.70
ADIN	TOUAREG CLOSEFIT	NP	41.315.078.01-2	No-80	15 Ncm	7.80	2.45	5.77	6.00	1.80	2.30	conical	65°	45° chamfer	Hex. 1.70
ADIN	TOUAREG CLOSEFIT	RP	41.316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70
ADIN	SWELL INTERNAL HEX	3.45	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	SWELL INTERNAL HEX	3.60	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	SWELL INTERNAL HEX	4.00	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	SWELL INTERNAL HEX	4.60	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	TOUAREG-S/TOUAREG-OS	3.45	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	TOUAREG-S/TOUAREG-OS	3.60	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	TOUAREG-S/TOUAREG-OS	4.00	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	TOUAREG-S/TOUAREG-OS	5.00	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ADIN	TOUAREG CLOSEFIT	WP	41.320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
ADVAN	GFA	UNIVERSAL GTB	41.316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70
ALFAGATE	CONICAL INTERNAL HEXAGON	3.50	41.316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70
ALFAGATE	MULTI-UNIT	UNIVERSAL	41.317.040.01-2	N1-72	25 Ncm	4.00	2.10	2.50	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
ALFAGATE	Bioactive S-line/SLA	3.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ALFAGATE	CONICAL INTERNAL HEXAGON	3.50	41.316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
ALFAGATE	MULTI-UNIT	UNIVERSAL	41.317.040.01-2	N1-72	25 Ncm	4.00	2.10	2.50	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
ALFAGATE	Bioactive S-line/SLA	3.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ALFAGATE	Bioactive S-line/SLA	3.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ALPHA BIO	NEO/NICE	3.50	41.316.071.01-2	1.60	20 Ncm	7.10	2.80	5.20	5.54	1.56	2.30	conical	60°	45° chamfer	Hex. 1.70
ALPHA BIO	MULTI UNIT	UNIVERSAL	41.317.041.01-2	N1-72	25 Ncm	4.10	1.90	2.30	2.47	1.63	2.30	conical	55°	45° chamfer	Hex. 1.70
ALPHA BIO	CONICAL STANDARD (CS)	STANDARD	41.317.070.01-2	N1-72	25 Ncm	7.00	2.20	4.96	5.80	1.20	2.30	conical	30°	45° chamfer	Hex. 1.70
ALPHA BIO	DFI	UNIVERSAL	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ALPHA BIO	ATID	UNIVERSAL	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ALPHA BIO	SPI	UNIVERSAL	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ALPHA BIO	ICE	UNIVERSAL	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ANCLADEN	ANCLALOCK	3.50	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ANKYLOS	BALANCE BASE	MULTI UNIT	41.316.048.02-2	1.60	20 Ncm	4.80	2.85	3.00	3.58	1.22	2.30	conical	31°	45° chamfer	Hex. 1.70
ANKYLOS	C/X NP/RP/WP	3.5/4.5/5.5/7	41.318.077.01-2	1.80	25 Ncm	7.70	2.50	5.59	6.50	1.20	2.30	conical	30°	45° chamfer	Hex. 1.70
ANTHOGR	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
ANTHOGR	AXIOM REG/PX	3.40	41.316.079.01-2	1.60	20 Ncm	7.90	2.30	5.42	6.60	1.30	2.30	conical	20°	45° chamfer	Hex. 1.70
ANTHOGR	AXIOM REG/PX	4.00	41.316.079.01-2	1.60	20 Ncm	7.90	2.30	5.42	6.60	1.30	2.30	conical	20°	45° chamfer	Hex. 1.70
ANTHOGR	AXIOM REG/PX	4.60	41.316.079.01-2	1.60	20 Ncm	7.90	2.30	5.42	6.60	1.30	2.30	conical	20°	45° chamfer	Hex. 1.70
ANTHOGR	AXIOM REG/PX	5.20	41.316.079.01-2	1.60	20 Ncm	7.90	2.30	5.42	6.60	1.30	2.30	conical	20°	45° chamfer	Hex. 1.70
ANTHOGR	OSSFIT	4.80	41.320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
ANTHOGR	OSSFIT	6.50	41.320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
ANTHOGR	ANTHOFIT HE	R (4.1)	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
ANTHOGR	ANTHOFIT HE	L (5)	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
ANTHOGR	ANTHOFIT OI	3.50	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
ANTHOGR	ANTHOFIT OI	4.00	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
ANTHOGR	ANTHOFIT OI	5.00	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
ARDS	SMART	3.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ARDS	CLASSIC	3.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ARDS	PREMIUM	3.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
ASTRA	CONE 20°	REGULAR/WIDE	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
ASTRA	CONE 45°	REGULAR/WIDE	41.314.052.01-2	1.40	15 Ncm	5.20	2.90	1.80	-	3.40	2.30	straight	-	45° chamfer	Hex. 1.70
ASTRA	EVOLUTION	3.0	41.314.074.01-2	1.40	15 Ncm	7.40	3.55	5.00	5.99	1.41	2.30	conical	25°	45° chamfer	Hex. 1.70
ASTRA	EVOLUTION	3.60	41.316.072.01-2	1.60	20 Ncm	7.20	3.50	5.20	5.82	1.38	2.30	conical	30°	45° chamfer	Hex. 1.70
ASTRA	EVOLUTION	RP	41.316.076.01-2	1.60	20 Ncm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
ASTRA	EVOLUTION	4.20	41.318.074.01-2	1.80	25 Ncm	7.40	3.80	5.80	6.04	1.36	2.30	conical	50°	45° chamfer	Hex. 1.70
ASTRA	EVOLUTION	4.80	41.320.082.01-2	2.00	25 Ncm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
ASTRA	EVOLUTION	5.40	41.320.082.01-2	2.00	25 Ncm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
ASTRA	LILAC	WP	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
ASTRA	EVOLUTION EXTERNAL CONE	UNIVERSAL	41.318.045.01-2	1.80	25 Ncm	4.50	2.30	2.81	2.90	1.60	2.30	conical	70°	45° chamfer	Hex. 1.70
AVINENT	TRANSEPITELIAL	REGULAR	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
AVINENT	HI	3.50	41.317.073.01-2	N1-72	25 Ncm	7.30	2.50	5.50	5.77	1.53	2.30	conical	60°	45° chamfer	Hex. 1.70
AVINENT	HI	4.10	41.317.073.01-2	N1-72	25 Ncm	7.30	2.50	5.50	5.77	1.53	2.30	conical	60°	45° chamfer	Hex. 1.70
AVINENT	HE	3.50	41.318.051.01-2	1.80	25 Ncm	5.10	2.70	3.55	3.70	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
AVINENT	HE	4.10	41.320.068.01-2	2.00	25 Ncm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
AVINENT	HE	5.10	41.325.067.01-2	2.50	25 Ncm	6.70	4.60	5.10	-	1.60	2.85	straight	-	semi-sphere	Hex. 1.70
BEGO	MULTI PLUS	UNIVERSAL	41.314.046.01-2	1.40	15 Ncm	4.60	2.50	2.50	3.17	1.43	2.30	conical	35°	45° chamfer	Hex. 1.70
BEGO	RS/RSX	3.0	41.314.064.01-2	1.40	15 Ncm	6.40	2.20	4.21	5.15	1.25	2.30	conical	25°	45° chamfer	Hex. 1.70
BEGO	EXTERNAL HEXAGON	NARROW MINI (2.7/2.9/3.1)	41.316.059.01-2	1.60	20 Ncm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
BEGO	S/RI/RS/RSX	3.25/3.75	41.318.064.01-2	1.80	25 Ncm	6.40	3.45	4.70	5.10	1.30	2.30	conical	35°	45° chamfer	Hex. 1.70
BEGO	S/RI/RS/RSX	4.10	41.318.064.01-2	1.80	25 Ncm	6.40	3.45	4.70	5.10	1.30	2.30	conical	35°	45° chamfer	Hex. 1.70
BEGO	S/RI/RS/RSX	4.50	41.318.064.01-2	1.80	25 Ncm	6.40	3.45	4.70	5.10	1.30	2.30	conical	35°	45° chamfer	Hex. 1.70
BEGO	S/RI/RS/RSX	5.50	41.318.064.01-2	1.80	25 Ncm	6.40	3.45	4.70	5.10	1.30	2.30	conical	35°	45° chamfer	Hex. 1.70
BIO3 IMPLANTS	PROGRESSIVE / ADVANCED	WIDE	41.318.083.01-2	1.80	25 Ncm	8.30	4.25	6.79	6.95	1.35	2.30	conical	60°	45° chamfer	Hex. 1.70
BIOCONCEPT	BV TAPERED BONE LEVEL	MULTI UNIT	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
BIOCONCEPT	BC BONE LEVEL	NARROW	41.316.078.01-2	1.60	20 Ncm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
BIOCONCEPT	BC BONE LEVEL	REGULAR	41.316.078.01-2	1.60	20 Ncm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
BIOCONCEPT	BV TAPERED BONE LEVEL	NARROW	41.316.094.01-2	1.60	20 Ncm	9.40	2.90	7.62	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
BIOCONCEPT	BC TISSUE LEVEL STANDARD	REGULAR	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
BIOCONCEPT	BC TISSUE LEVEL STANDARD PLUS	REGULAR	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
BIOCONCEPT	BC TISSUE LEVEL TAPERED EFFECT	REGULAR	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
BIOCONCEPT	BV TAPERED BONE LEVEL	REGULAR	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
BIOGENESIS	ATICON (CONO 20°/45°)	REGULAR/WIDE	41.314.040.01-2	1.40	15 Ncm	4.00	1.85	2.00	2.76	1.24	2.30	conical	31°	45° chamfer	Hex. 1.70
BIOGENESIS	3ICON	MINI (PINK)	41.316.059.01-2	1.60	20 Ncm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
BIOGENESIS	3ICON	WIDE (YELLOW)	41.320.047.01-2	2.00	25 Ncm	4.70	3.00	3.30	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
BIOGENESIS	3ICON	REGULAR (BLUE)	41.320.060.01-2	2.00	25 Ncm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
BIOGENESIS	ITICON	4.80	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
BIOGENESIS	ATICON	PINK/BLUE	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
BIOHORIZONS	MULTI-UNIT		41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
BIOHORIZONS	INTERNAL	3 (GREY)	41.317.065.01-2	N1-72	25 Ncm	6.50	2.40	4.70	5.18	1.33	2.30	conical	45°	45° chamfer	Hex. 1.70
BIOHORIZONS	INTERNAL	3.5 (YELLOW)	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
BIOHORIZONS	INTERNAL	4.5 (GREEN)	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
BIOHORIZONS	INTERNAL	5.7 (BLUE)	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
BIOHORIZONS	TAPERED TISSUE LEVEL	3.5 (YELLOW)	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
BIOHORIZONS	TAPERED TISSUE LEVEL	4.5 (GREEN)	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
BIOHORIZONS	TAPERED TISSUE LEVEL	5,7 (BLUE)	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIOHORIZONS	EXTERNAL	3,7 (YELLOW)	41.318.065.01-2	1,80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOHORIZONS	EXTERNAL	5,2 (BLUE)	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOLOK	EXTERNAL HEXAGON	3,45	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
BIOMET 3i	LOW PROFILE	UNIVERSAL	41.314.039.01-2	1,40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
BIOMET 3i	CERTAIN	3,40	41.316.084.01-2	1,60	20 N·cm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
BIOMET 3i	CERTAIN	4,10	41.316.084.01-2	1,60	20 N·cm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
BIOMET 3i	CERTAIN	5,00	41.316.084.01-2	1,60	20 N·cm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
BIOMET 3i	OSSEOTITE	4,10	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOMET 3i	OSSEOTITE	5,00	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOMET 3i	OSSEOTITE	3,40	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
BIONER	HIKELT	3,95	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIONER	HIKELT	4,90	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIONER	IKELT	5,00	41.320.047.01-2	2,00	25 N·cm	4,70	3,00	3,30	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BIONER	IKELT/BIKELT	4,10	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIONER	TRANSEPITELIAL	REGULAR	41.314.039.01-2	1,40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
BIONER	TRANSEPITELIAL TOP DM	3,50	41.314.039.01-2	1,40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
BIONER	TRANSEPITELIAL A-5M	REGULAR	41.314.039.01-2	1,40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
BIONER	TOP-DM	3,50	41.316.074.01-2	1,60	20 N·cm	7,40	2,70	5,43	6,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
BIONER	TOP-DM	4,00	41.316.074.01-2	1,60	20 N·cm	7,40	2,70	5,43	6,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
BIONER	TOP-DM	5,00	41.316.074.01-2	1,60	20 N·cm	7,40	2,70	5,43	6,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
BIONNOVATION	INTERNAL HEXAGON	SP	41.316.084.01-2	1,60	20 N·cm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
BIONNOVATION	EXTERNAL HEXAGON	SP	41.318.065.01-2	1,80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIONNOVATION	INTERNAL HEXAGON	RP	41.318.075.01-2	1,80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BIONNOVATION	INTERNAL HEXAGON	WP	41.318.075.01-2	1,80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BIONNOVATION	EXTERNAL HEXAGON	4,10	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
BIONNOVATION	EXTERNAL HEXAGON	5,00	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
BIONNOVATION	EXTERNAL HEXAGON	RP	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
BIONNOVATION	EXTERNAL HEXAGON	WP	41.325.067.01-2	2,50	25 N·cm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
BIOTEC	SPR/CIM	3,30	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIOTEC	SPR/CIM	3,75	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIOTEC	SPR/SPTT/CIM	4,20	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIOTEC	SPR/SPTT/CIM	5,00	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIOTEC	SPR/SPTT/CIM	6,00	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BIOTECH	KONTACT	3,00	41.312.078.01-2	1,20	15 N·cm	7,80	2,65	5,97	6,55	1,25	2,30	conical	45°	45° chamfer	Hex. 1,70
BIOTECH	MULTI UNIT	UNIVERSAL	41.314.040.02-2	1,40	15 N·cm	4,00	1,70	2,25	2,70	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
BIOTECH	KONTACT	3,6/4,2/4,8/5,4	41.314.076.01-2	1,40	15 N·cm	7,60	2,40	5,88	6,35	1,25	2,30	conical	45°	45° chamfer	Hex. 1.213

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
BIOTECH	SMILEA /SMILEA CONIC EXTERNO	NARROW (3,6)	41320.065.01-2	2.00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
BIOTEM	BR FIXTURE	MINI	41316.059.01-2	1.60	20 N·cm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOTEM	AR FIXTURE R-TYPE/T-TYPE	MINI	41316.076.01-2	1.60	20 N·cm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOTEM	BR FIXTURE	REGULAR	41320.060.01-2	2.00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOTEM	AR FIXTURE R-TYPE/T-TYPE	REGULAR	41320.090.01-2	2.00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOTEM	AR FIXTURE R-TYPE/T-TYPE	WIDE	41320.090.01-2	2.00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BIOTEM	BR FIXTURE	WIDE	41325.067.01-2	2.50	25 N·cm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
BLUE SKY BIO	MAX	NP	41316.076.01-2	1.60	20 N·cm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BLUE SKY BIO	INTERNAL HEX	NARROW	41317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BLUE SKY BIO	INTERNAL HEX	REGULAR	41317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BLUE SKY BIO	INTERNAL HEX	WIDE	41317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BLUE SKY BIO	TRILOBE	3,50	41318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BLUE SKY BIO	ONE STAGE	REGULAR	41320.074.01-2	2.00	25 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BLUE SKY BIO	TRILOBE	4,30	41320.090.01-2	2.00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BLUE SKY BIO	TRILOBE	5,00	41320.090.01-2	2.00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BREDENT MEDICAL	COPA SKY	4,00	41316.064.02-2	1.60	20 N·cm	6,40	1,60	4,08	5,48	0,92	2,30	conical	15°	45° chamfer	Hex. 1,70
BREDENT MEDICAL	NARROW SKY	NP 3,5	41318.083.01-2	1.80	25 N·cm	8,30	4,25	6,79	6,95	1,35	2,30	conical	60°	45° chamfer	Hex. 1,70
BREDENT MEDICAL	BLUE SKY / BLUE SKY CLASSIC	4,00	41318.083.01-2	1.80	25 N·cm	8,30	4,25	6,79	6,95	1,35	2,30	conical	60°	45° chamfer	Hex. 1,70
BTI	MULTI-IM CONNECTION	UNIVERSAL 4,1	41314.039.01-2	1.40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
BTI	MULTI-IM CONNECTION	WIDE 5,5	41314.039.01-2	1.40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
BTI	MULTI-IM	UNIVERSAL 4,1	41314.039.01-2	1.40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
BTI	INTERNAL CONNECTION	ESTRECHA 3,5	41318.065.01-2	1.80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTI	INTERNAL CONNECTION	UNIVERSAL 4,1	41318.065.01-2	1.80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTI	INTERNAL CONNECTION	UNIVERSAL 4,1	41318.065.01-2	1.80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTI	INTERNAL CONNECTION	WIDE 5,5	41318.065.01-2	1.80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTI	EXTERNAL CONNECTION	TINY 3,5	41320.051.01-2	2.00	25 N·cm	5,10	3,10	3,60	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTI	EXTERNAL CONNECTION	UNIVERSAL 4,1	41320.060.01-2	2.00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTI	EXTERNAL CONNECTION	WIDE 5,5	41320.060.01-2	2.00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTK	KLASSIC INTERNAL	3,5 IR	41317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BTK	KONIC INTERNAL	3,5 IR	41317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
BTK	KLASSIC/KONIC EXTERNAL	3,4 EN	41318.065.01-2	1.80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
BTK	KLASSIC INTERNAL	4,5 IM	41318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BTK	KLASSIC INTERNAL	5 IW	41318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BTK	KONIC INTERNAL	4,5 IM	41318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BTK	KONIC INTERNAL	5 IW	41318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
BTK	KLASSIC/KONIC EXTERNAL	4,1 EN	41320.060.01-2	2.00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
BTK	KLASSIC/KONIC EXTERNAL	5 EW	41.320.060.01-2	2.00	25 Ncm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
BTK	SAFE TL	REGULAR	41.320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
BW	EXTERNAL HEXAGON	3.30	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
BW	CONICAL INTERNAL HEXAGON	4.00	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
BW	EXTERNAL HEXAGON	5.00	41.320.047.01-2	2.00	25 Ncm	4.70	3.00	3.30	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
BW	EXTERNAL HEXAGON	4.10	41.320.060.01-2	2.00	25 Ncm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
CAMLOG	COMFOUR	3.3/3.8/4.3	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
CAMLOG	CONELOG SCREW LINE	3.30	41.316.080.01-2	1.60	20 Ncm	8.00	3.00	6.30	6.51	1.49	2.30	conical	60°	45° chamfer	Hex. 1.70
CAMLOG	CONELOG SCREW LINE	3.80	41.316.080.01-2	1.60	20 Ncm	8.00	3.00	6.30	6.51	1.49	2.30	conical	60°	45° chamfer	Hex. 1.70
CAMLOG	CONELOG SCREW LINE	4.30	41.316.080.01-2	1.60	20 Ncm	8.00	3.00	6.30	6.51	1.49	2.30	conical	60°	45° chamfer	Hex. 1.70
CAMLOG	CAMLOG SCREW LINE	3.30	41.316.094.01-2	1.60	20 Ncm	9.40	2.90	7.62	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CAMLOG	CAMLOG SCREW LINE	3.80	41.316.094.01-2	1.60	20 Ncm	9.40	2.90	7.62	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CAMLOG	CAMLOG SCREW LINE	4.30	41.316.094.01-2	1.60	20 Ncm	9.40	2.90	7.62	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CAMLOG	ISY		41.316.094.01-2	1.60	20 Ncm	9.40	2.90	7.62	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CAMLOG	COMFOUR	5/6	41.320.047.01-2	2.00	25 Ncm	4.70	3.00	3.30	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
CAMLOG	CAMLOG SCREW LINE	5.00	41.320.094.01-2	2.00	25 Ncm	9.40	3.00	7.84	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CAMLOG	CAMLOG SCREW LINE	6.00	41.320.094.01-2	2.00	25 Ncm	9.40	3.00	7.84	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CONEXAO	AR MORSE	UNIVERSAL	41.316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CONEXAO	AR MORSE	UNIVERSAL	41.316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CONEXAO	FLASH	UNIVERSAL	41.316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CONEXAO	TORQ	UNIVERSAL	41.316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CONEXAO	EXPAND	UNIVERSAL	41.316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
CONEXAO	GRIP ACTIVES	UNIVERSAL	41.320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
CONEXAO	GRIP ACTIVES	UNIVERSAL	41.320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
CONEXAO	GRIP POROUS	UNIVERSAL	41.320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
CONEXAO	GRIP POROUS	UNIVERSAL	41.320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
CORTEX	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
CORTEX	CONICAL PLATFORM	NP	41.314.070.01-2	1.40	15 Ncm	7.00	2.30	5.39	5.65	1.35	2.30	conical	60°	45° chamfer	Hex. 1.70
CORTEX	CONICAL PLATFORM	RP	41.316.076.01-2	1.60	20 Ncm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
CORTEX	MAGIX	RP	41.316.076.01-2	1.60	20 Ncm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
CORTEX	SATURN	3.50	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
CORTEX	CLASSIX/DYNAMIX	3.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
CORTEX	CONICAL PLATFORM	WP	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
C-TECH	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
C-TECH	EL ESTHETIC LINE	HEX INT. 3.8/4.3/5.1	41.318.080.01-2	1.80	25 Ncm	8.00	4.00	6.50	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
DENTAL IMPLANT KA	BONE LEVEL	NP	41.316.076.01-2	1.60	20 Ncm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
DENTAL IMPLANT KA	TISSUE LEVEL	RP	41.316.076.01-2	1.60	20 Ncm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAL IMPLANT KA	TISSUE LEVEL	WP	41.316.076.01-2	1.60	20 Ncm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAL TECH	IMPLOGIC	HEX. INT. 4.5	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTAURUM	MULTI-UNIT BAR	CONO	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAURUM	MULTI-UNIT ANGLEDFIX	CONO	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAURUM	MULTI-UNIT BRIDGE	CONO	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAURUM	TIOLOX	3.50	41.316.055.01-2	1.60	20 Ncm	5.50	2.43	2.85	4.20	1.30	2.30	conical	23°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOX	4.50	41.316.055.01-2	1.60	20 Ncm	5.50	2.43	2.85	4.20	1.30	2.30	conical	23°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOGIC	S	41.316.081.01-2	1.60	20 Ncm	8.10	3.00	6.35	6.73	1.37	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOGIC	M	41.316.081.01-2	1.60	20 Ncm	8.10	3.00	6.35	6.73	1.37	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOGIC	L	41.316.081.01-2	1.60	20 Ncm	8.10	3.00	6.35	6.73	1.37	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTEGRIS	INTERNAL HEXAGON	3.3	41.318.071.01-2	1.80	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTAURUM	INTERNAL HEXAGON	3.8/4.1	41.318.071.01-2	1.80	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTAURUM	INTERNAL HEXAGON	4.5	41.318.071.01-2	1.80	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTAURUM	INTERNAL HEXAGON	5.5	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTIS	S-CLEAN	REGULAR	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIS	S-CLEAN	MINI	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIS	S-CLEAN	WIDE	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIUM	SIMPLE LINE II	HEX. EXT. 4.8	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTIUM	SIMPLE LINE II	HEX. EXT. 6.5	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTIUM	MULTI-UNIT NR LINE	4.30	41.316.048.01-2	1.60	20 Ncm	4.80	2.60	2.93	3.50	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIUM	MULTI-UNIT NR LINE	5.50	41.316.048.01-2	1.60	20 Ncm	4.80	2.60	2.93	3.50	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIUM	NR LINE	3.20	41.316.084.02-2	1.60	20 Ncm	8.40	2.70	5.85	6.85	1.55	2.30	conical	30°	45° chamfer	Hex. 1.70
DENTIUM	NR LINE	3.60	41.316.084.02-2	1.60	20 Ncm	8.40	2.70	5.85	6.85	1.55	2.30	conical	30°	45° chamfer	Hex. 1.70
DENTIUM	NR LINE	4.30	41.316.084.02-2	1.60	20 Ncm	8.40	2.70	5.85	6.85	1.55	2.30	conical	30°	45° chamfer	Hex. 1.70
DENTIUM	NR LINE	5.00	41.316.084.02-2	1.60	20 Ncm	8.40	2.70	5.85	6.85	1.55	2.30	conical	30°	45° chamfer	Hex. 1.70
DENTIUM	NR LINE	6.00	41.316.084.02-2	1.60	20 Ncm	8.40	2.70	5.85	6.85	1.55	2.30	conical	30°	45° chamfer	Hex. 1.70
DENTIUM	MULTI-UNIT SUPER LINE AND IMPLANTUM	3.6/4/4.5/5/6	41.318.051.02-2	1.80	25 Ncm	5.10	2.70	3.38	3.65	1.45	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIUM	SIMPLE LINE II	OCT. INT. 4.8	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
DENTIUM	SIMPLE LINE II	OCT. INT. 6.5	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
DENTIUM	SUPER LINE AND IMPLANTUM	3.6/4/4.5/5/6	41.320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTOFLEX	SUPREME MORSE CONE	3.30	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
DENTOFLEX	SUPREME MORSE CONE	4.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
DIO	UF	NARROW	41.314.067.01-2	1.40	15 Ncm	6.70	2.31	5.00	5.45	1.25	2.30	conical	45°	45° chamfer	Hex. 1.70
DIO	SCREW MULTI UNIT	UNIVERSAL	41.316.040.01-2	1.60	20 Ncm	4.00	2.07	2.30	2.47	1.53	2.30	conical	60°	45° chamfer	Hex. 1.70
DIO	EXTERNAL	NARROW 3.5	41.316.059.01-2	1.60	20 Ncm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
DIO	SM SYSTEM	NARROW	41.318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	INTERNAL OCTA	4.8	41.320.044.01-2	2,00	25 N·cm	4,40	2,45	2,95	3,10	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	INTERNAL OCTA	6.0	41.320.044.01-2	2,00	25 N·cm	4,40	2,45	2,95	3,10	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	INTERNAL OCTA	6.5	41.320.044.01-2	2,00	25 N·cm	4,40	2,45	2,95	3,10	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	EXTERNAL	REGULAR 4.1	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	SM SYSTEM	REGULAR/WIDE	41.320.074.01-2	2,00	25 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	UF	REGULAR	41.320.079.01-2	2,00	25 N·cm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	EXTERNAL	WIDE 5.1	41.325.067.01-2	2,50	25 N·cm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	EXTERNAL HEXAGON	3.30	41.318.065.01-2	1.80	25 N·cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	3.30	41.318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	EXTERNAL HEXAGON	4.10	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
DSP BIOMEDICAL	EXTERNAL HEXAGON	5.00	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	3.75	41.320.074.01-2	2,00	25 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	4.00	41.320.074.01-2	2,00	25 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	5.00	41.320.074.01-2	2,00	25 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	CRESTAL LEVEL	3.50	41.320.074.01-2	2,00	25 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MULTI-UNIT CONICAL ABUTMENT	UNIVERSAL	41.314.039.01-2	1.40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
EASY IMPLANT	HEXCEL-S	4.10	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	HEXCEL-S	5.00	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	HEXCEL-S	3.30	41.320.065.01-2	2,00	25 N·cm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EASY IMPLANT	MASTER C	4.50	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER C	5.00	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER L	3.30	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER L	3.75	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER L	4.25	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER L	4.75	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER S	3.75	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER S	4.25	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER S	4.75	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	ALL-SPIRAL	REGULAR	41.318.048.01-2	1.80	25 N·cm	4,80	2,80	3,22	3,65	1,15	2,30	conical	31°	semi-sphere	Hex. 1,70
ECKERMANN	WINNER	3/3.5/4	41.318.071.01-2	1.80	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ECKERMANN	WINNER	4/4.5/5	41.318.071.01-2	1.80	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ECKERMANN	TRIPLO	4.00	41.318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	TRIPLO	4.50	41.318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	TRIPLO	5.00	41.318.075.01-2	1.80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	HEXAGON	4.10	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ELITE MEDICAL	EXTERNAL HEXAGON	NARROW	41.316.059.01-2	1.60	20 N·cm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,358

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
ELITE MEDICAL	EXTERNAL HEXAGON	REGULAR	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ELITE MEDICAL	EXTERNAL HEXAGON	WIDE	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO MORSE	NARROW	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO MORSE	REGULAR	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO MORSE	WIDE	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	COLOSSO ECI/CRC EVOLUTION	NARROW	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	COLOSSO ECI/CRC EVOLUTION	REGULAR	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	COLOSSO ECI/CRC EVOLUTION	WIDE	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	NARROW	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	REGULAR	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	WIDE	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EQUINOX	UNITI	3.30	41.314.076.01-2	1,40	15 Ncm	7,60	2,40	5,88	6,35	1,25	2,30	conical	45°	45° chamfer	Hex. 1,70
EMFILS	UNITI	3.70	41.316.081.01-2	1,60	20 Ncm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
EMFILS	UNITI	4.30	41.316.081.01-2	1,60	20 Ncm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
EMFILS	UNITI	5.60	41.316.081.01-2	1,60	20 Ncm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
EUROTEKNIKA	MULTI-UNIT TETRA	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
EUROTEKNIKA	MULTI-UNIT ANGULATED CONICAL	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
EUROTEKNIKA	NATURACTIS/NATURALL	3.00	41.314.074.01-2	1,40	15 Ncm	7,40	3,55	5,00	5,99	1,41	2,30	conical	25°	45° chamfer	Hex. 1,70
EUROTEKNIKA	NATURACTIS	3.40	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATURACTIS	3.80	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATURACTIS	4.30	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATURACTIS	4.80	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATEA/NATURALL	NARROW	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATEA/NATURALL	REGULAR	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	UNEVA	4.10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	AESTHETICA	4.20	41.320.070.01-2	2,00	25 Ncm	7,00		5,60	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	AESTHETICA	4.80	41.320.070.01-2	2,00	25 Ncm	7,00	3,00	5,60	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	AESTHETICA	6.50	41.320.070.01-2	2,00	25 Ncm	7,00	3,00	5,60	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATEA/NATURALL	WIDE	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EVO IMPLANTS	EVOSOFT / EOPLUS	HEX. INT. 3,8/4,2/4,8/5,5	41.316.094.01-2	1,60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
EUROTEKNIKA	EVOSHORT	HEX. INT. 4,8/5,5/6,2	41.316.094.01-2	1,60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
FRIALIT	MULTI-UNIT MP		41.316.044.01-2	1,60	20 Ncm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
GALIMPLANT	STRAIGHT MULTI-POSITION ABUTMENT	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
GALIMPLANT	MULTI-POSITION ANGLED ABUTMENT	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
GALIMPLANT	INTERNAL CONNECTION	3.50	41.316.071.01-2	1,60	20 Ncm	7,10	2,80	5,20	5,54	1,56	2,30	conical	60°	45° chamfer	Hex. 1,70
GALIMPLANT	INTERNAL CONNECTION	4.00	41.316.071.01-2	1,60	20 Ncm	7,10	2,80	5,20	5,54	1,56	2,30	conical	60°	45° chamfer	Hex. 1,70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
GALIMPLANT	INTERNAL CONNECTION	5.00	41.316.071.01-2	1.60	20 N·cm	7.10	2.80	5.20	5.54	1.56	2.30	conical	60°	45° chamfer	Hex. 1.70
GALIMPLANT	STRAIGHT-AESTHETIC MULTI-POSITION ABUTMENT	UNIVERSAL	41.320.047.01-2	2.00	25 N·cm	4.70	3.00	3.30	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
GALIMPLANT	EXTERNAL CONNECTION	4.00	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
CORTEX	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1.40	15 N·cm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
CORTEX	CONICAL PLATFORM	NP	41.314.070.01-2	1.40	15 N·cm	7.00	2.30	5.39	5.65	1.35	2.30	conical	60°	45° chamfer	Hex. 1.70
CORTEX	CONICAL PLATFORM	RP	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
CORTEX	MAGIX	RP	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
CORTEX	SATURN	3.50	41.317.071.01-2	N1-72	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
CORTEX	CLASSIX/DYNAMIX	3.75	41.317.071.01-2	N1-72	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
CORTEX	CONICAL PLATFORM	WP	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
C-TECH	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1.40	15 N·cm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
C-TECH	EL ESTHETIC LINE	HEX. INT. 3.8/4.3/5.1	41.318.080.01-2	1.80	25 N·cm	8.00	4.00	6.50	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
DENTAL IMPLANT KA	BONE LEVEL	NP	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAL IMPLANT KA	TISSUE LEVEL	RP	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAL IMPLANT KA	TISSUE LEVEL	WP	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAL TECH	IMPLOGIC	HEX. INT. 4.5	41.317.071.01-2	N1-72	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTAURUM	MULTI-UNIT BAR	CONO	41.316.044.01-2	1.60	20 N·cm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAURUM	MULTI-UNIT ANGLEDFIX	CONO	41.316.044.01-2	1.60	20 N·cm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAURUM	MULTI-UNIT BRIDGE	CONO	41.316.044.01-2	1.60	20 N·cm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTAURUM	TIOLOX	3.50	41.316.055.01-2	1.60	20 N·cm	5.50	2.43	2.85	4.20	1.30	2.30	conical	23°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOX	4.50	41.316.055.01-2	1.60	20 N·cm	5.50	2.43	2.85	4.20	1.30	2.30	conical	23°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOGIC	S	41.316.081.01-2	1.60	20 N·cm	8.10	3.00	6.35	6.73	1.37	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOGIC	M	41.316.081.01-2	1.60	20 N·cm	8.10	3.00	6.35	6.73	1.37	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTAURUM	TIOLOGIC	L	41.316.081.01-2	1.60	20 N·cm	8.10	3.00	6.35	6.73	1.37	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTEGRIS	INTERNAL HEXAGON	3.3	41.318.071.01-2	1.80	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTEGRIS	INTERNAL HEXAGON	3.8/4.1	41.318.071.01-2	1.80	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTEGRIS	INTERNAL HEXAGON	4.5	41.318.071.01-2	1.80	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
DENTEGRIS	INTERNAL HEXAGON	5.5	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTIS	S-CLEAN	REGULAR	41.320.079.01-2	2.00	25 N·cm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIS	S-CLEAN	MINI	41.320.079.01-2	2.00	25 N·cm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIS	S-CLEAN	WIDE	41.320.079.01-2	2.00	25 N·cm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
DENTIUM	SIMPLE LINE II	HEX. EXT. 4.8	41.316.044.01-2	1.60	20 N·cm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTIUM	SIMPLE LINE II	HEX. EXT. 6.5	41.316.044.01-2	1.60	20 N·cm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
DENTIUM	MULTI-UNIT NR LINE	4.30	41.316.048.01-2	1.60	20 N·cm	4.80	2.60	2.93	3.50	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
DENTIUM	MULTI-UNIT NR LINE	5,50	41.316.048.01-2	1,60	20 Ncm	4,80	2,60	2,93	3,50	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
DENTIUM	NR LINE	3,20	41.316.084.02-2	1,60	20 Ncm	8,40	2,70	5,85	6,85	1,55	2,30	conical	30°	45° chamfer	Hex. 1,70
DENTIUM	NR LINE	3,60	41.316.084.02-2	1,60	20 Ncm	8,40	2,70	5,85	6,85	1,55	2,30	conical	30°	45° chamfer	Hex. 1,70
DENTIUM	NR LINE	4,30	41.316.084.02-2	1,60	20 Ncm	8,40	2,70	5,85	6,85	1,55	2,30	conical	30°	45° chamfer	Hex. 1,70
DENTIUM	NR LINE	5,00	41.316.084.02-2	1,60	20 Ncm	8,40	2,70	5,85	6,85	1,55	2,30	conical	30°	45° chamfer	Hex. 1,70
DENTIUM	NR LINE	6,00	41.316.084.02-2	1,60	20 Ncm	8,40	2,70	5,85	6,85	1,55	2,30	conical	30°	45° chamfer	Hex. 1,70
DENTIUM	MULTI-UNIT SUPER LINE AND IMPLANTUM	3,6/4/4,5/5/6	41.318.051.02-2	1,80	25 Ncm	5,10	2,70	3,38	3,65	1,45	2,30	conical	45°	45° chamfer	Hex. 1,70
DENTIUM	SIMPLE LINE II	OCT. INT. 4,8	41.320.067.01-2	2,00	25 Ncm	6,70	2,30	3,59	5,70	1,00	2,58	conical	15°	45° chamfer	Hex. 1,70
DENTIUM	SIMPLE LINE II	OCT. INT. 6,5	41.320.067.01-2	2,00	25 Ncm	6,70	2,30	3,59	5,70	1,00	2,58	conical	15°	45° chamfer	Hex. 1,70
DENTIUM	SUPER LINE AND IMPLANTUM	3,6/4/4,5/5/6	41.320.079.01-2	2,00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
DENTOFLEX	SUPREME MORSE CONE	3,30	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DENTIUM	SUPREME MORSE CONE	4,50	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	UF	NARROW	41.314.067.01-2	1,40	15 Ncm	6,70	2,31	5,00	5,45	1,25	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	SCREW MULTI UNIT	UNIVERSAL	41.316.040.01-2	1,60	20 Ncm	4,00	2,07	2,30	2,47	1,53	2,30	conical	60°	45° chamfer	Hex. 1,70
DIO	EXTERNAL	NARROW 3,5	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	SM SYSTEM	NARROW	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	INTERNAL OCTA	4,8	41.320.044.01-2	2,00	25 Ncm	4,40	2,45	2,95	3,10	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	INTERNAL OCTA	6,0	41.320.044.01-2	2,00	25 Ncm	4,40	2,45	2,95	3,10	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	INTERNAL OCTA	6,5	41.320.044.01-2	2,00	25 Ncm	4,40	2,45	2,95	3,10	1,30	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	EXTERNAL	REGULAR 4,1	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	SM SYSTEM	REGULAR/WIDE	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DIO	UF	REGULAR	41.320.079.01-2	2,00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
DIO	EXTERNAL	WIDE 5,1	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	EXTERNAL HEXAGON	3,30	41.318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	3,30	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	EXTERNAL HEXAGON	4,10	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
DSP BIOMEDICAL	EXTERNAL HEXAGON	5,00	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	3,75	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	4,00	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	HEXAGON CONE MORSE	5,00	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
DSP BIOMEDICAL	CRESTAL LEVEL	3,50	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MULTI-UNIT CONICAL ABUTMENT	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
EASY IMPLANT	HEXCEL-S	4,10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	HEXCEL-S	5,00	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	HEXCEL-S	3,30	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EASY IMPLANT	MASTER C	4,50	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER C	5,00	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,343

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
EASY IMPLANT	MASTER L	3,30	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER L	3,75	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER L	4,25	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER L	4,75	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER S	3,75	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER S	4,25	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EASY IMPLANT	MASTER S	4,75	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	ALL-SPIRAL	REGULAR	41.318.048.01-2	1,80	25Ncm	4,80	2,80	3,22	3,65	1,15	2,30	conical	31°	semi-sphere	Hex. 1,70
ECKERMANN	WINNER	3/3,5/4	41.318.071.01-2	1,80	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ECKERMANN	WINNER	4/4,5/5	41.318.071.01-2	1,80	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ECKERMANN	TRIPLO	4,00	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	TRIPLO	4,50	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	TRIPLO	5,00	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ECKERMANN	HEXAGON	4,10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ELITE MEDICAL	EXTERNAL HEXAGON	NARROW	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ELITE MEDICAL	EXTERNAL HEXAGON	REGULAR	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ELITE MEDICAL	EXTERNAL HEXAGON	WIDE	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO MORSE	NARROW	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO MORSE	REGULAR	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO MORSE	WIDE	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	COLOSSO ECI/CRC EVOLUTION	NARROW	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	COLOSSO ECI/CRC EVOLUTION	REGULAR	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	COLOSSO ECI/CRC EVOLUTION	WIDE	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	NARROW	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	REGULAR	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	WIDE	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
EQUINOX	UNITI	3,30	41.314.076.01-2	1,40	15 Ncm	7,60	2,40	5,88	6,35	1,25	2,30	conical	45°	45° chamfer	Hex. 1,70
EQUINOX	UNITI	3,70	41.316.081.01-2	1,60	20 Ncm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
EQUINOX	UNITI	4,30	41.316.081.01-2	1,60	20 Ncm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
EQUINOX	UNITI	5,60	41.316.081.01-2	1,60	20 Ncm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
EUROTEKNIKA	MULTI-UNIT TETRA	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
EUROTEKNIKA	MULTI-UNIT ANGULATED CONICAL	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
EUROTEKNIKA	NATURACTIS/NATURALL	3,00	41.314.074.01-2	1,40	15 Ncm	7,40	3,55	5,00	5,99	1,41	2,30	conical	25°	45° chamfer	Hex. 1,70
EUROTEKNIKA	NATURACTIS	3,40	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATURACTIS	3,80	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATURACTIS	4,30	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATURACTIS	4,80	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,380

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
EUROTEKNIKA	NATEA/NATURALL	NARROW	41316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATEA/NATURALL	REGULAR	41316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	UNEVA	4,10	41320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	AESTHETICA	4,20	41320.070.01-2	2,00	25 Ncm	7,00		5,60	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	AESTHETICA	4,80	41320.070.01-2	2,00	25 Ncm	7,00	3,00	5,60	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	AESTHETICA	6,50	41320.070.01-2	2,00	25 Ncm	7,00	3,00	5,60	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
EUROTEKNIKA	NATEA/NATURALL	WIDE	41320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
EVO IMPLANTS	EVOSOFT / EVOPLUS	HEX. INT. 3,8/4,2/4,8/5,5	41316.094.01-2	1,60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
EVO IMPLANTS	EVOSHORT	HEX. INT. 4,8/5,5/6,2	41316.094.01-2	1,60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
FRIALIT	MULTI-UNIT MP		41316.044.01-2	1,60	20 Ncm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
GALIMPLANT	STRAIGHT MULTI-POSITION ABUTMENT	UNIVERSAL	41314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
GALIMPLANT	MULTI-POSITION ANGLED ABUTMENT	UNIVERSAL	41314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
GALIMPLANT	INTERNAL CONNECTION	3,50	41316.071.01-2	1,60	20 Ncm	7,10	2,80	5,20	5,54	1,56	2,30	conical	60°	45° chamfer	Hex. 1,70
GALIMPLANT	INTERNAL CONNECTION	4,00	41316.071.01-2	1,60	20 Ncm	7,10	2,80	5,20	5,54	1,56	2,30	conical	60°	45° chamfer	Hex. 1,70
GALIMPLANT	INTERNAL CONNECTION	5,00	41316.071.01-2	1,60	20 Ncm	7,10	2,80	5,20	5,54	1,56	2,30	conical	60°	45° chamfer	Hex. 1,70
GALIMPLANT	STRAIGHT-AESTHETIC MULTI-POSITION ABUTMENT	UNIVERSAL	41320.047.01-2	2,00	25 Ncm	4,70	3,00	3,30	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
GALIMPLANT	EXTERNAL CONNECTION	4,00	41320.068.01-2	2,00	25 Ncm	6,80	4,40	5,30	5,40	1,40	2,30	conical	60°	45° chamfer	Hex. 1,70
GC TECH	AADVA STANDARD / TAPERED IMPLANTS	NARROW	41316.086.01-2	1,60	20 Ncm	8,60	3,00	7,20	-	1,40	2,30	straight	-	45° chamfer	Hex. 1,70
GC TECH	AADVA STANDARD / TAPERED IMPLANTS	REGULAR	41316.086.01-2	1,60	20 Ncm	8,60	3,00	7,20	-	1,40	2,30	straight	-	45° chamfer	Hex. 1,70
GC TECH	AADVA STANDARD / TAPERED IMPLANTS	WIDE	41316.086.01-2	1,60	20 Ncm	8,60	3,00	7,20	-	1,40	2,30	straight	-	45° chamfer	Hex. 1,70
GLOBAL D (TEKKA)	IN-KONE MULTI UNIT	4,30	41318.044.01-2	1,80	25 Ncm	4,40	2,75	3,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
GLOBAL D (TEKKA)	IN-KONE UNIVERSAL	5,00	41318.077.02-2	1,80	25 Ncm	7,70	2,20	6,08	6,35	1,35	2,30	conical	60°	45° chamfer	Hex. 1,70
GLOBAL D (TEKKA)	IN-KONE PRIMO	5,00	41318.077.02-2	1,80	25 Ncm	7,70	2,20	6,08	6,35	1,35	2,30	conical	60°	45° chamfer	Hex. 1,70
GMI	UNIVERSAL	PS-RP 4,8	41314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
GMI	AVANTGARD	RP	41316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
GMI	PHOENIX	NARROW 3,3	41318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
GMI	FRONTIER	RP 3,3	41318.071.01-2	1,80	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
GMI	FRONTIER	WP 4,3	41318.071.01-2	1,80	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
GMI	PHOENIX	STANDARD 4,1	41320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
GMI	PHOENIX	WIDE 5,1 (E/C 2,7)	41320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
GMI	PHOENIX	WIDE 5,1 (E/C 3,4)	41325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
GMI	PHOENIX	EXTRA WIDE 6	41325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
GT MEDICAL	BEST-FIT EXTERNAL HEXAGON	NARROW	41316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
GT MEDICAL	BEST-FIT EXTERNAL HEXAGON	REGULAR	41320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
GT MEDICAL	BEST-FIT INTERNAL OCTAGON	REGULAR	41320.067.01-2	2,00	25 Ncm	6,70	2,30	3,59	5,70	1,00	2,58	conical	15°	45° chamfer	Hex. 1,70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
GT MEDICAL	BEST-FIT INTERNAL HEXAGON	WIDE	41320.094.01-2	2.00	25 Ncm	9.40	3.00	7.84	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
GT MEDICAL	BEST-FIT EXTERNAL HEXAGON	WIDE	41325.067.01-2	2.50	25 Ncm	6.70	4.60	5.10	-	1.60	2.85	straight	-	semi-sphere	Hex. 1.70
HAHN IMPLANT	TAPERED IMPLANT	3 (Grey)	41314.070.01-2	1.40	15 Ncm	7.00	2.30	5.39	5.65	1.35	2.30	conical	60°	45° chamfer	Hex. 1.70
HAHN IMPLANT	TAPERED IMPLANT	3.5 (Purple)	41316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70
HAHN IMPLANT	TAPERED IMPLANT	5 (Blue)	41320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
HAHN IMPLANT	TAPERED IMPLANT	7 (Green)	41320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
HAHN IMPLANT	MULTI UNIT ABUTMENT SYSTEM	UNIVERSAL	41314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
HIOSSSEN	ETII SA / ETIII SA	MINI	41316.094.01-2	1.60	20 Ncm	9.40	2.90	7.62	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
HIOSSSEN	ETIII BA	MINI	41316.094.01-2	1.60	20 Ncm	9.40	2.90	7.62	8.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
HIOSSSEN	ETII SA / ETIII SA	REGULAR	41320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
HIOSSSEN	ETIII BA	REGULAR	41320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
HI-TEC	TAPERED SELF THREAD	3.50	41317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
HI-TEC	TAPERED SELF THREAD	4.50	41317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
HI-TEC	LOGIC PLUS	3.70	41317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
HI-TEC	LOGIC PLUS	3.90	41317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
IBS	MAGICORE	4/4.5	41316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
IBS	MAGICORE	5.50	41316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
IBS	MAGIC FC	3.80	41320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
IBS	N.R. FIX	3.80	41320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
IDO IMPLANT	I DO IMPLANT	UNIVERSAL	41320.079.01-2	2.00	25 Ncm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
IHDE DENTAL (IMBIODENT)	BONE LEVEL PLUS	3.30	41316.078.01-2	1.60	20 Ncm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
IHDE DENTAL (IMBIODENT)	BONE LEVEL PLUS	4.10	41316.078.01-2	1.60	20 Ncm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
IHDE DENTAL (IMBIODENT)	BONE LEVEL PLUS	4.80	41316.078.01-2	1.60	20 Ncm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
IMPLACIL DE BORTOLI	EXTERNAL HEXAGON	3.50	41316.059.01-2	1.60	20 Ncm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLACIL DE BORTOLI	INTERNAL HEXAGON	3.50	41318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLACIL DE BORTOLI	CONE MORSE	3.50	41318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLACIL DE BORTOLI	CONE MORSE	4.00	41318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLACIL DE BORTOLI	CONE MORSE	4.50	41318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLACIL DE BORTOLI	CONE MORSE	5.50	41318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLACIL DE BORTOLI	EXTERNAL HEXAGON	4.00	41320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
IMPLACIL DE BORTOLI	EXTERNAL HEXAGON	5.00	41320.065.01-2	2.00	25 Ncm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
IMPLACIL DE BORTOLI	INTERNAL HEXAGON	4.00	41.320.074.01-2	2.00	25 Ncm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLACIL DE BORTOLI	INTERNAL HEXAGON	5.00	41.320.074.01-2	2.00	25 Ncm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	SWISH ACTIVE	3.00	41.316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70
IMPLANT DIRECT	INTERACTIVE	3.00	41.316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70
IMPLANT DIRECT	LEGACY	3.50	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
IMPLANT DIRECT	LEGACY	4.50	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
IMPLANT DIRECT	LEGACY	5.70	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
IMPLANT DIRECT	REACTIVE	3.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REPLANT	3.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REPLUS	3.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	SWISHPLANT	4.80	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
IMPLANT DIRECT	SWISHPLANT	6.50	41.320.067.01-2	2.00	25 Ncm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
IMPLANT DIRECT	SWISH ACTIVE	3.40	41.320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
IMPLANT DIRECT	INTERACTIVE	3.40	41.320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REACTIVE	4.30	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REACTIVE	5.00	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REPLANT	4.30	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REPLANT	5.00	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REPLANT	6.00	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REPLUS	4.30	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT DIRECT	REPLUS	5.00	41.320.090.01-2	2.00	25 Ncm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
IMPLANT GENESIS	AKTIVE SYSTEM	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
IMPLASA HÖSCHT	DOUBLEGOTIC/GOTIC/NORDIC	3.50	41.318.074.01-2	1.80	25 Ncm	7.40	3.80	5.80	6.04	1.36	2.30	conical	50°	45° chamfer	Hex. 1.70
INTRA LOCK	INTRAHEX MULTI UNIT 17°	UNIVERSAL	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
INTRA LOCK	INTRAHEX MULTI UNIT	UNIVERSAL	41.316.059.01-2	1.60	20 Ncm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
INTRA LOCK	INTRAHEX	3.75/4	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
INTRA LOCK	INTRAHEX	4.75	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
INTRA LOCK	UNIHEX	REGULAR	41.320.060.01-2	2.00	25 Ncm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
INTRA LOCK	UNIHEX	WIDE	41.320.070.01-2	2.00	25 Ncm	7.00	3.00	5.60	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
INTRA LOCK	INDEX	NARROW	41.320.082.01-2	2.00	25 Ncm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
INTRA LOCK	INDEX	STANDARD	41.320.082.01-2	2.00	25 Ncm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
INTRA LOCK	INDEX	WIDE	41.320.082.01-2	2.00	25 Ncm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
JDENTALCARE	CONICAL ABUTMENT	UNIVERSAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
JDENTALCARE	JDEVOLUTION S	3.20	41.316.073.01-2	1.60	20 Ncm	7.30	2.20	4.71	5.56	1.74	2.30	conical	35°	45° chamfer	Hex. 1.70
JDENTALCARE	JDEVOLUTION/JDEVOLUTION PLUS	3.70	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
JDENTALCARE	JDEVOLUTION/JDEVOLUTION PLUS	4.00	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
JDENTALCARE	JDEVOLUTION/JDEVOLUTION PLUS	5.00	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
JDENTALCARE	JD ICON PLUS	3.70	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON PLUS	4.00	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON PLUS	4.80	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON PLUS T	3.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON PLUS T	3.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON PLUS T	3.50	41.318.075.01-2	1.80	25 Ncm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON	3.90	41.320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON	4.00	41.320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
JDENTALCARE	JD ICON	4.70	41.320.075.01-2	2.00	25 Ncm	7.50	2.75	5.90	6.19	1.31	2.30	conical	35°	semi-sphere	Hex. 1.70
KENTRON	INTERNAL HEXAGON	NP 3.8	41.316.084.01-2	1.60	20 Ncm	8.40	3.50	6.80	-	1.60	2.30	straight	-	semi-sphere	Hex. 1.70
KENTRON	INTERNAL HEXAGON	RP 4.5	41.316.084.01-2	1.60	20 Ncm	8.40	3.50	6.80	-	1.60	2.30	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	PRIMACONNEX MULTI UNIT	SD/RD	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	PRIMACONNEX MULTI UNIT	WD	41.316.044.01-2	1.60	20 Ncm	4.40	2.50	2.90	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	INTERNAL TILOBE PRIMACONNEX	3.50	41.318.065.01-2	1.80	25 Ncm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	INTERNAL TILOBE PRIMACONNEX	4.10	41.318.065.01-2	1.80	25 Ncm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	INTERNAL TILOBE PRIMACONNEX	5.00	41.318.065.01-2	1.80	25 Ncm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	RESTORE	RD 4.1	41.320.060.01-2	2.00	25 Ncm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	RESTORE	WD 5	41.325.054.01-2	2.50	25 Ncm	5.40	4.60	4.10	-	1.30	2.87	straight	-	semi-sphere	Hex. 1.70
KEYSTONE	RESTORE	SD 3.3	41.318.044.01-2	1.80	25 Ncm	4.40	2.75	3.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
KLOCKNER	KL	NARROW	41.316.059.01-2	1.60	20 Ncm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KLOCKNER	VEGA	NV	41.316.074.01-2	1.60	20 Ncm	7.40	2.70	5.43	6.00	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
KLOCKNER	NK2-SK2	4.20	41.318.044.01-2	1.80	25 Ncm	4.40	2.75	3.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
KLOCKNER	ESSENTIAL CONE	4.50	41.318.067.01-2	1.80	25 Ncm	6.70	2.35	4.93	5.40	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70
KLOCKNER	ESSENTIAL CONE	6.00	41.318.067.01-2	1.80	25 Ncm	6.70	2.35	4.93	5.40	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70
KLOCKNER	KL	REGULAR	41.320.060.01-2	2.00	25 Ncm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
KLOCKNER	KL	WIDE	41.325.067.01-2	2.50	25 Ncm	6.70	4.60	5.10	-	1.60	2.85	straight	-	semi-sphere	Hex. 1.70
KLOCKNER	MULTI UNIT PERMANENT	UNIVERSAL	41.314.040.02-2	1.40	15 Ncm	4.00	1.70	2.25	2.70	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70
KLOCKNER	VEGA	RV	41.318.076.01-2	1.80	25 Ncm	7.60	2.52	5.73	6.20	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
KLOCKNER	ESSENTIAL CONE OCTACONE 12°	4.50	41.318.051.02-2	1.80	25 Ncm	5.10	2.70	3.38	3.65	1.45	2.30	conical	45°	45° chamfer	Hex. 1.70
KLOCKNER	ESSENTIAL CONE OCTACONE 25°	4.50	41.318.051.02-2	1.80	25 Ncm	5.10	2.70	3.38	3.65	1.45	2.30	conical	45°	45° chamfer	Hex. 1.70
LASAK	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1.40	15 Ncm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
LASAK	BIONIQ	2.90	41.314.084.01-2	1.40	15 Ncm	8.40	2.50	5.92	6.85	1.55	2.30	conical	35°	45° chamfer	Hex. 1.70
LASAK	BIONIQ	3.5/4/5	41.316.084.02-2	1.60	20 Ncm	8.40	2.70	5.85	6.85	1.55	2.30	conical	30°	45° chamfer	Hex. 1.70
LEADER	TIXOS MORSE CONEXION	3.60	41.316.084.01-2	1.60	20 Ncm	8.40	3.50	6.80	-	1.60	2.30	straight	-	semi-sphere	Hex. 1.70
LEADER	TIXOS INTERNAL HEXAGON	3.50	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
LEADER	TIXOS INTERNAL HEXAGON	4.00	41.317.071.01-2	N1-72	25 Ncm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.519

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
LEADER	TIXOS INTERNAL HEXAGON	5.00	41.317.071.01-2	N1-72	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
LEADER	TIXOS INTERNAL HEXAGON	6.00	41.317.071.01-2	N1-72	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
LEADER	TIXOS EXTERNAL HEXAGON	4.10	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
LEADER	TIXOS EXTERNAL HEXAGON	5.00	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
LIFECORE	EXTERNAL HEXAGON	NARROW	41.318.065.01-2	1.80	25 N·cm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
LIFECORE	EXTERNAL HEXAGON	REGULAR	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
LIFECORE	EXTERNAL HEXAGON	WIDE	41.325.067.01-2	2.50	25 N·cm	6.70	4.60	5.10	-	1.60	2.85	straight	-	semi-sphere	Hex. 1.70
MEDENTIKA	MULTI-UNIT	UNIVERSAL	41.314.039.01-2	1.40	15 N·cm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
MEDENTIS	ICX TEMPLANT	3.45	41.316.078.01-2	1.60	20 N·cm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
MEDENTIS	ICX TEMPLANT	3.75	41.316.078.01-2	1.60	20 N·cm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
MEDENTIS	ICX TEMPLANT	4.10	41.316.078.01-2	1.60	20 N·cm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
MEDENTIS	ICX TEMPLANT	4.80	41.316.078.01-2	1.60	20 N·cm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
MEGAGEN	MULTI UNIT N TYPE	UNIVERSAL	41.314.039.01-2	1.40	15 N·cm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
MEGAGEN	MULTI UNIT S TYPE	UNIVERSAL	41.314.045.01-2	1.40	15 N·cm	4.50	2.30	2.50	3.25	1.25	2.30	conical	31°	45° chamfer	Hex. 1.70
MEGAGEN	MINI	3/3.25	41.314.067.01-2	1.40	15 N·cm	6.70	2.31	5.00	5.45	1.25	2.30	conical	45°	45° chamfer	Hex. 1.70
MEGAGEN	ANYONE EXTERNAL	SMALL 3.5	41.316.059.01-2	1.60	20 N·cm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	EZ PLUS	SMALL	41.316.084.01-2	1.60	20 N·cm	8.40	3.50	6.80	-	1.60	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	ANYRIDGE	SMALL	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	ANYRIDGE	REGULAR	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	ANYRIDGE	WIDE	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	ANYRIDGE	SUPERWIDE	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	ANYONE EXTERNAL	REGULAR 4.1	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	ANYONE EXTERNAL	REGULAR 4.5	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	ANYONE EXTERNAL	WIDE 5	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
MEGAGEN	ANYONE EXTERNAL	SUPERWIDE 5.5	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
MEGAGEN	ANYONE INTERNAL	GENERAL	41.320.079.01-2	2.00	25 N·cm	7.90	3.30	6.32	6.50	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
MEGAGEN	EZ PLUS	REGULAR	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	EZ PLUS	WIDE	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	RESCUE	4.50	41.325.054.01-2	2.50	25 N·cm	5.40	4.60	4.10	-	1.30	2.87	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	RESCUE	5.00	41.325.054.01-2	2.50	25 N·cm	5.40	4.60	4.10	-	1.30	2.87	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	RESCUE	5.50	41.325.054.01-2	2.50	25 N·cm	5.40	4.60	4.10	-	1.30	2.87	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	RESCUE	6.00	41.325.054.01-2	2.50	25 N·cm	5.40	4.60	4.10	-	1.30	2.87	straight	-	semi-sphere	Hex. 1.70
MEGAGEN	CONE ABUTMENT	5.80	41.320.044.01-2	2.00	25 N·cm	4.40	2.45	2.95	3.10	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70
MEGAGEN	CONE ABUTMENT	3.80	41.320.050.01-2	2.00	25 N·cm	5.00	2.80	3.35	3.60	1.40	2.30	conical	31°	semi-sphere	Hex. 1.70
MEGAGEN	CONE ABUTMENT	4.80	41.320.050.01-2	2.00	25 N·cm	5.00	2.80	3.35	3.60	1.40	2.30	conical	31°	semi-sphere	Hex. 1.70
MICRODENT	SYSTEM	3.50	41.318.044.01-2	1.80	25 N·cm	4.40	2.75	3.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
MICRODENT	UNIVERSAL	3.50	41.318.044.01-2	1.80	25 N·cm	4.40	2.75	3.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.556

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
MICRODENT	EKTOS	3,50	41.318.071.01-2	1,80	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
MICRODENT	GENIUS	3,50	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	GENIUS	4,00	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	GENIUS	4,50	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	TRYLOGIC	3,80	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	TRYLOGIC	4,30	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	SYSTEM	4,20	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	SYSTEM	5,10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	UNIVERSAL	4,10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	UNIVERSAL	5,10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	TRYLOGIC	5,10	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
MICRODENT	MK	5,60	41.325.054.01-2	2,50	25 Ncm	5,40	4,60	4,10	-	1,30	2,87	straight	-	semi-sphere	Hex. 1,70
MIS	SEVEN/M4	NARROW	41.316.071.01-2	1,60	20 Ncm	7,10	2,80	5,20	5,54	1,56	2,30	conical	60°	45° chamfer	Hex. 1,70
MIS	C1	NARROW	41.316.071.01-2	1,60	20 Ncm	7,10	2,80	5,20	5,54	1,56	2,30	conical	60°	45° chamfer	Hex. 1,70
MIS	SEVEN/M4	STANDARD	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
MIS	SEVEN/M4	WIDE	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
MIS	V3	STANDARD	41.317.073.01-2	N1-72	25 Ncm	7,30	2,50	5,50	5,77	1,53	2,30	conical	60°	45° chamfer	Hex. 1,70
MIS	C1	STANDARD	41.317.073.01-2	N1-72	25 Ncm	7,30	2,50	5,50	5,77	1,53	2,30	conical	60°	45° chamfer	Hex. 1,70
MIS	C1	WIDE	41.317.073.01-2	N1-72	25 Ncm	7,30	2,50	5,50	5,77	1,53	2,30	conical	60°	45° chamfer	Hex. 1,70
MIS	LANCE	STANDARD	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MIS	LANCE	WIDE	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MIS	MULTI-UNIT	GENERAL	41.316.044.01-2	1,60	20 Ncm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MOZO GRAU	MG INHEX	2,3 MINI	41.314.070.01-2	1,40	15 Ncm	7,00	2,30	5,39	5,65	1,35	2,30	conical	60°	45° chamfer	Hex. 1,70
MOZO GRAU	MG INHEX	2,8 STANDARD	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MOZO GRAU	MG OSSEOUS	4,1 STANDARD	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MOZO GRAU	MG OSSEOUS	3,4 MINI	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
MOZO GRAU	MG INHEX	3,8 MAXI	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MOZO GRAU	MG OSSEOUS	5 MAXI	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
MPI	CM INTERNAL HEXAGON	REGULAR	41.316.076.01-2	1,60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MPI	HE EXTERNAL HEXAGON	5,00	41.320.047.01-2	2,00	25 Ncm	4,70	3,00	3,30	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
MPI	HE EXTERNAL HEXAGON	3,50	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MPI	HE EXTERNAL HEXAGON	4,10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
MPI	CM INTERNAL HEXAGON	WIDE	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NEOBIOTECH	MULTI-UNIT	4,80	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
NEOBIOTECH	MULTI-UNIT	4,80	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
NEOBIOTECH	EB EXTERNAL SYSTEM	NARROW	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NEOBIOTECH	IS CONE INTERNAL	S-NARROW	41.316.094.01-2	1,60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1.593

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
NEOBIOTECH	EB EXTERNAL SYSTEM	REGULAR	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NEOBIOTECH	EB EXTERNAL SYSTEM	WIDE 3l	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NEOBIOTECH	IS CONE INTERNAL	NARROW	41.320.079.01-2	2,00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
NEOBIOTECH	IS CONE INTERNAL	REGULAR 4	41.320.079.01-2	2,00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
NEOBIOTECH	IS CONE INTERNAL	REGULAR 4,5	41.320.079.01-2	2,00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
NEOBIOTECH	IS CONE INTERNAL	WIDE 5	41.320.079.01-2	2,00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
NEOBIOTECH	EB EXTERNAL SYSTEM	WIDE B	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
NEODENT	CONO MORSE MINI-PILAR	REGULAR	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
NEODENT	II PLUS	4,30	41.316.084.01-2	1,60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1,70
NEODENT	HELIX GM/DRIVE GM/TITAMAX GM	REGULAR	41.316.084.02-2	1,60	20 Ncm	8,40	2,70	5,85	6,85	1,55	2,30	conical	30°	45° chamfer	Hex. 1,70
NEODENT	SMART HE	3,30	41.318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NEODENT	SMART HE	4,10	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
NEODENT	SMART HE	4,30	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
NEODENT	SMART HE	5,00	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
NEOSS	PROACTIVE /STRAIGHT/TAPERED	PROACTIVE	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NOBEL BIOCARE	MULTI-UNIT	REGULAR	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
NOBEL BIOCARE	STANDARD	REGULAR	41.314.045.01-2	1,40	15 Ncm	4,50	2,30	2,50	3,25	1,25	2,30	conical	31°	45° chamfer	Hex. 1,70
NOBEL BIOCARE	BRANEMARK	NARROW	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NOBEL BIOCARE	ACTIVE	NARROW	41.316.073.01-2	1,60	20 Ncm	7,30	2,20	4,71	5,56	1,74	2,30	conical	35°	45° chamfer	Hex. 1,70
NOBEL BIOCARE	MULTI-UNIT	WIDE	41.318.044.01-2	1,80	25 Ncm	4,40	2,75	3,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NOBEL BIOCARE	REPLACE	NARROW	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NOBEL BIOCARE	BRANEMARK	REGULAR	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NOBEL BIOCARE	ACTIVE	REGULAR	41.320.075.01-2	2,00	25 Ncm	7,50	2,75	5,90	6,19	1,31	2,30	conical	35°	semi-sphere	Hex. 1,70
NOBEL BIOCARE	ACTIVE	WIDE	41.320.075.01-2	2,00	25 Ncm	7,50	2,75	5,90	6,19	1,31	2,30	conical	35°	semi-sphere	Hex. 1,70
NOBEL BIOCARE	REPLACE	REGULAR	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NOBEL BIOCARE	REPLACE	WIDE	41.320.090.01-2	2,00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NOBEL BIOCARE	BRANEMARK	WIDE	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
NORIS MEDICAL	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1,40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
NORIS MEDICAL	TUFF	3,75	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
NORIS MEDICAL	TUFF TT	3,75	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
NORIS MEDICAL	ONIX	3,75	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
NORIS MEDICAL	CORTICAL	3,75	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
NORIS MEDICAL	PTERYCORE	3,75	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
NORIS MEDICAL	PTERYFIT	3,75	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
NORMON	NORMOIMPLANT HI	3,50	41.318.071.01-2	1,80	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
NORMON	NORMOIMPLANT HE	4,10	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
NOVA IMPLANTS	PSI/PCI	3,75	41.318.068.01-2	1,80	25 Ncm	6,80	4,00	5,25	5,40	1,40	2,30	conical	60°	45° chamfer	Hex. 1,630

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
NUCLEOSS	T5 SHORT IMPLANT	5,80	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NUCLEOSS	BONE LEVEL T PURE	4,20	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NUCLEOSS	BONE LEVEL T PURE	5,00	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NUCLEOSS	BONE LEVEL T PURE	6,00	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NUCLEOSS	TISSUE LEVEL T3	5,00	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
NUCLEOSS	TISSUE LEVEL T3	6,00	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSSEOLIFE	OEX	NP 3,3	41.318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OSSEOLIFE	OEX	RP 4,1	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OSSEOLIFE	OEX	WP 5	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
OSSEOLIFE	OIN	RP 4	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSSEOLIFE	OIN	WP 5	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSSTEM	TS MULTI UNIT	MINI / REGULAR	41.316.044.01-2	1,60	20 Ncm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OSSTEM	US	MINI 3,5	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OSSTEM	TS	MINI	41.316.094.01-2	1,60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
OSSTEM	TS MULTI UNIT	REGULAR	41.320.048.01-2	2,00	25 Ncm	4,80	2,70	3,30	3,40	1,40	2,30	conical	60°	45° chamfer	Hex. 1,70
OSSTEM	US	REGULAR 4,1	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OSSTEM	US	WIDE PS 5	41.320.060.01-2	2,00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OSSTEM	SS	REGULAR 4,8	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSSTEM	SS	WIDE 6	41.320.074.01-2	2,00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSSTEM	TS	REGULAR	41.320.079.01-2	2,00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
OSSTEM	US	WIDE 5,1	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
OSTEOPLUS	SHI	3,50	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
OSTEOPLUS	SHI	4,50	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
OSTEOPLUS	IRIDIUM	3,50	41.318.068.01-2	1,80	25 Ncm	6,80	4,00	5,25	5,40	1,40	2,30	conical	60°	45° chamfer	Hex. 1,70
OSTEOPLUS	ARGENTUM	4,20	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSTEOPLUS	IRIDIUM	4,20	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSTEOPLUS	PLATINIUM	4,20	41.318.075.01-2	1,80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
OSTEOPLUS	SHE	3,45	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
OSTEOPLUS	SHE	4,00	41.320.065.01-2	2,00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
OT MEDICAL	OT-F2	NARROW	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OT MEDICAL	OT-F2	REGULAR	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OT MEDICAL	OT-F3	NARROW	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OT MEDICAL	OT-F3	REGULAR	41.316.059.01-2	1,60	20 Ncm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OT MEDICAL	OT-F2	REGULAR	41.318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OT MEDICAL	OT-F2	WIDE	41.318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OT MEDICAL	OT-F3	REGULAR	41.318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
OT MEDICAL	OT-F3	WIDE	41.318.065.01-2	1,80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,630

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
OXY	EXTERNAL HEX	MINI	41.318.071.01-2	1.80	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
OXY	CONICAL CONNECTION	REGULAR	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
OXY	EXTERNAL HEX	STANDARD	41.320.070.01-2	2.00	25 N·cm	7.00	3.00	5.60	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
OXY	EXTERNAL HEX	MAXI	41.320.070.01-2	2.00	25 N·cm	7.00	3.00	5.60	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
OXY	INTERNAL CONNECTION	MINI	41.320.082.01-2	2.00	25 N·cm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
OXY	INTERNAL CONNECTION	STANDARD	41.320.082.01-2	2.00	25 N·cm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
OXY	INTERNAL CONNECTION	MAXI	41.320.082.01-2	2.00	25 N·cm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
PALTOP	ADVANCED CLASSIC	NARROW	41.316.064.01-2	1.60	20 N·cm	6.40	3.15	4.70	5.00	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	ADVANCED +	NARROW	41.316.064.01-2	1.60	20 N·cm	6.40	3.15	4.70	5.00	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	DYNAMIC	NARROW	41.316.064.01-2	1.60	20 N·cm	6.40	3.15	4.70	5.00	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	CONICAL ACTIVE	STANDARD	41.316.071.01-2	1.60	20 N·cm	7.10	2.80	5.20	5.54	1.56	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	ADVANCED CLASSIC	STANDARD	41.317.073.01-2	N1-72	25 N·cm	7.30	2.50	5.50	5.77	1.53	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	ADVANCED +	STANDARD	41.317.073.01-2	N1-72	25 N·cm	7.30	2.50	5.50	5.77	1.53	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	DYNAMIC	STANDARD	41.317.073.01-2	N1-72	25 N·cm	7.30	2.50	5.50	5.77	1.53	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	DIVA	STANDARD	41.317.073.01-2	N1-72	25 N·cm	7.30	2.50	5.50	5.77	1.53	2.30	conical	60°	45° chamfer	Hex. 1.70
PALTOP	CONICAL ACTIVE	NARROW	41.318.043.01-2	1.80	25 N·cm	4.30	2.00	2.52	2.70	1.60	2.30	conical	55°	45° chamfer	Hex. 1.70
PALTOP	UNIVERSAL	MULTI-UNIT	41.318.043.01-2	1.80	25 N·cm	4.30	2.00	2.52	2.70	1.60	2.30	conical	55°	45° chamfer	Hex. 1.70
PHIBO	TSA S3	3.70	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
PHIBO	TSH S2	3.30	41.318.065.01-2	1.80	25 N·cm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
PHIBO	TSA S4	4.70	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
PHIBO	TSA S5	6.00	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
PHIBO	TSH S5	5.00	41.320.048.01-2	2.00	25 N·cm	4.80	2.70	3.30	3.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
PHIBO	TSH S3	4.00	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
PHIBO	TSH S4	4.00	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
PROCLINIC	Aqua CM	3.50	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
PROCLINIC	CILÍNDRICO INTERNO	3.50	41.317.071.01-2	N1-72	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
PROCLINIC	CILÍNDRICO EXTERNO	3.3 MINI	41.318.065.01-2	1.80	25 N·cm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
PROCLINIC	CILÍNDRICO EXTERNO	4.1 STANDARD	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
PROCLINIC	CILÍNDRICO EXTERNO	5 MAXI	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
PROCLINIC	SP OCTA	4.80	41.320.067.01-2	2.00	25 N·cm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
PROCLINIC	SP OCTA	6.50	41.320.067.01-2	2.00	25 N·cm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
RADHEX	PHE	3.50	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
RADHEX	PHI	3.50	41.318.071.01-2	1.80	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
RADHEX	PHI	4.50	41.318.071.01-2	1.80	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
RADHEX	PHE	4.10	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
ROSTER DENT	CONICAL GOLD	3.50	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
ROSTER DENT	CONICAL GOLD P4	3.75	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.630

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
ROSTER DENT	HEX	RD	41.320.065.01-2	2.00	25 Ncm	6,50	2,70	5,00	-	1,50	2,30	straight	-	45° chamfer	Hex. 1.70
ROSTER DENT	PLUS	4.80	41.320.067.01-2	2.00	25 Ncm	6,70	2,30	3,59	5,70	1,00	2,58	conical	15°	45° chamfer	Hex. 1.70
ROSTER DENT	CONICAL GOLD	4.30	41.320.090.01-2	2.00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
ROSTER DENT	CONICAL GOLD P4	4.30	41.320.090.01-2	2.00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
SAMO	BIOSPARK / BIOROUGH TAPERED / CYLINDRICAL	4.50	41.316.094.01-2	1.60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1.70
SEWON MEDIX	IH SYSTEM	UNIVERSAL	41.314.039.01-2	1.40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1.70
SEWON MEDIX	IH2 SLA SYSTEM / IH2 RBM SYSTEM	MINI	41.316.094.01-2	1.60	20 Ncm	9,40	2,90	7,62	8,00	1,40	2,30	conical	45°	45° chamfer	Hex. 1.70
SEWON MEDIX	IH2 SLA SYSTEM / IH2 RBM SYSTEM	REGULAR	41.320.079.01-2	2.00	25 Ncm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1.70
SGS	S6/S7	3.75	41.314.039.01-2	1.40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1.70
SGS	P7N	3.00	41.314.074.01-2	1.40	15 Ncm	7,40	3,55	5,00	5,99	1,41	2,30	conical	25°	45° chamfer	Hex. 1.70
SGS	P1/P7	3.75	41.317.071.01-2	N1-72	25 Ncm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1.70
SIC INVENT	HEXAGONAL SYSTEM SICace	3.30	41.316.079.02-2	1.60	20 Ncm	7,90	3,00	6,30	-	1,60	2,30	straight	-	45° chamfer	Hex. 1.70
SIC INVENT	HEXAGONAL SYSTEM SICace	4.20	41.316.079.02-2	1.60	20 Ncm	7,90	3,00	6,30	-	1,60	2,30	straight	-	45° chamfer	Hex. 1.70
SIGNO VINCES	Mini Pilar Cónico		41.314.039.01-2	1.40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1.70
SIGNO VINCES	DUOCON	CM 3.8	41.316.076.01-2	1.60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	COMPACT	CM 3.8	41.316.076.01-2	1.60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	INFRA	CM	41.316.076.01-2	1.60	20 Ncm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	DUO	3.30	41.318.065.01-2	1.80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	INTEGRA	3.30	41.318.065.01-2	1.80	25 Ncm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	DUO	4.10	41.320.074.01-2	2.00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	INTEGRA	4.10	41.320.074.01-2	2.00	25 Ncm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	DUOCON	CM 4.6	41.320.090.01-2	2.00	25 Ncm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70
SIGNO VINCES	DUO/INTEGRA	WIDE 5	41.325.067.01-2	2,50	25 Ncm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	CONICAL ABUTMENT	Universal	41.314.039.01-2	1.40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1.70
SIN IMPLANTS	Mini Pilar	Universal	41.314.039.01-2	1.40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1.70
SIN IMPLANTS	Micro Multi Unit	Universal	41.314.039.01-2	1.40	15 Ncm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1.70
SIN IMPLANTS	SW MORSE	3.50	41.316.084.01-2	1.60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	SW MORSE	3.80	41.316.084.01-2	1.60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	SW MORSE	4.50	41.316.084.01-2	1.60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	SW MORSE	5.50	41.316.084.01-2	1.60	20 Ncm	8,40	3,50	6,80	-	1,60	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	HI SW	3.80	41.318.075.01-2	1.80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	HI SW	4.50	41.318.075.01-2	1.80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	HI SW	5.50	41.318.075.01-2	1.80	25 Ncm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	REVOLUTION COMPACT	5.00	41.320.047.01-2	2.00	25 Ncm	4,70	3,00	3,30	-	1,40	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	MULTIFUNCTIONAL ABUTMENT	Universal	41.320.047.01-2	2.00	25 Ncm	4,70	3,00	3,30	-	1,40	2,30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	TRYON	4.10	41.320.060.01-2	2.00	25 Ncm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
SIN IMPLANTS	REVOLUTION	4.10	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	STYLUS	4.10	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SIN IMPLANTS	TRYON	5.00	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
SIN IMPLANTS	REVOLUTION	3.40	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
SIN IMPLANTS	REVOLUTION	5.00	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
SIN IMPLANTS	STYLUS	3.40	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
SIN IMPLANTS	STYLUS	5.00	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	COMPACT CONICAL	4.80	41.314.039.01-2	1.40	15 N·cm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	COMPACT CONICAL	6.00	41.314.039.01-2	1.40	15 N·cm	3.90	1.80	2.10	-	1.80	2.40	straight	-	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	DEEP CONICAL	2.45	41.314.070.01-2	1.40	15 N·cm	7.00	2.30	5.39	5.65	1.35	2.30	conical	60°	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	MSC	3.00	41.316.059.01-2	1.60	20 N·cm	5.90	3.00	4.40	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SOUTHERN IMPLANT	DEEP CONICAL	2.95/3.1	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SOUTHERN IMPLANT	INTERNAL HEX	UNIVERSAL	41.317.071.01-2	N1-72	25 N·cm	7.10	2.60	5.54	5.65	1.45	2.30	conical	70°	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	TRI-NEX	3.50	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SOUTHERN IMPLANT	EXTERNAL HEXAGON	3.40	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	IT CONNECTION	4.80	41.320.067.01-2	2.00	25 N·cm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	IT CONNECTION	6.50	41.320.067.01-2	2.00	25 N·cm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	EXTERNAL HEXAGON	4.10	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	EXTERNAL HEXAGON	5.00	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	EXTERNAL HEXAGON	6.00	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SOUTHERN IMPLANT	TRI-NEX	4.30	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SOUTHERN IMPLANT	TRI-NEX	5.00	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SOUTHERN IMPLANT	TRI-NEX	6.00	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SOUTHERN IMPLANT	DEEP CONICAL	4.10	41.320.090.01-2	2.00	25 N·cm	9.00	4.00	7.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
STERI-OSS	EXTERNAL CONNECTION	3.30	41.318.065.01-2	1.80	25 N·cm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
STERN GOLD	STERN EX	4.10	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
STERN GOLD	STERN IC	4.80	41.320.067.01-2	2.00	25 N·cm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
STRAUMANN	SCREW-RETAINED NC/RC		41.314.043.01-2	1.40	15 N·cm	4.30	1.80	2.03	2.90	1.40	2.30	conical	35°	45° chamfer	Hex. 1.70
STRAUMANN	BONE LEVEL MULTI-BASE		41.314.045.01-2	1.40	15 N·cm	4.50	2.30	2.50	3.25	1.25	2.30	conical	31°	45° chamfer	Hex. 1.70
STRAUMANN	BONE LEVEL	2.90	41.314.080.01-2	1.40	15 N·cm	8.00	2.10	4.96	6.80	1.20	2.30	conical	15°	45° chamfer	Hex. 1.70
STRAUMANN	BLX	RB	41.316.066.01-2	1.60	20 N·cm	6.60	1.95	4.78	5.20	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
STRAUMANN	BLX	WB	41.316.066.01-2	1.60	20 N·cm	6.60	1.95	4.78	5.20	1.40	2.30	conical	45°	45° chamfer	Hex. 1.70
STRAUMANN	BONE LEVEL	3.30	41.316.078.01-2	1.60	20 N·cm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
STRAUMANN	BONE LEVEL	4.10	41.316.078.01-2	1.60	20 N·cm	7.84	2.20	5.51	7.04	0.80	2.30	conical	15°	45° chamfer	Hex. 1.70
STRAUMANN	STANDARD PLUS NARROW NECK 3.3	3.50	41.318.065.01-2	1.80	25 N·cm	6.50	2.83	5.00	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
STRAUMANN	EXTERNAL OCTAGON	WIDE 6.5	41.320.044.01-2	2.00	25 N·cm	4.40	2.45	2.95	3.10	1.30	2.30	conical	45°	45° chamfer	Hex. 1.70
STRAUMANN	EXTERNAL OCTAGON	REGULAR 4.8	41.320.050.01-2	2.00	25 N·cm	5.00	2.80	3.35	3.60	1.40	2.30	conical	31°	semi-sphere	Hex. 1.70
STRAUMANN	SYNOCTA	REGULAR 4.8	41.320.050.01-2	2.00	25 N·cm	5.00	2.80	3.35	3.60	1.40	2.30	conical	31°	semi-sphere	Hex. 1.70
STRAUMANN	INTERNAL OCTAGON	REGULAR 4.8	41.320.067.01-2	2.00	25 N·cm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
STRAUMANN	INTERNAL OCTAGON	WIDE 6.5	41.320.067.01-2	2.00	25 N·cm	6.70	2.30	3.59	5.70	1.00	2.58	conical	15°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	OUTLINK2	SLIM 3	41.318.051.01-2	1.80	25 N·cm	5.10	2.70	3.55	3.70	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	OUTLINK2	3.30	41.318.051.01-2	1.80	25 N·cm	5.10	2.70	3.55	3.70	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	OUTLINK2	4.1 SP	41.318.068.01-2	1.80	25 N·cm	6.80	4.00	5.25	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	PREMIUM KOHNO	3.30	41.318.068.01-2	1.80	25 N·cm	6.80	4.00	5.25	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	PREMIUM KOHNO	3.80	41.318.068.01-2	1.80	25 N·cm	6.80	4.00	5.25	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	PRAMA	3.40	41.318.068.01-2	1.80	25 N·cm	6.80	4.00	5.25	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	GLOBAL	3.80	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SWEDEN & MARTINA	GLOBAL	4.30	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SWEDEN & MARTINA	GLOBAL	4.80	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SWEDEN & MARTINA	GLOBAL	5.50	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SWEDEN & MARTINA	SHELTA	3.80	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SWEDEN & MARTINA	SHELTA	4.30	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SWEDEN & MARTINA	SHELTA	5.00	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SWEDEN & MARTINA	OUTLINK2	4.10	41.320.048.01-2	2.00	25 N·cm	4.80	2.70	3.30	3.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	OUTLINK2	5.00	41.320.048.01-2	2.00	25 N·cm	4.80	2.70	3.30	3.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
SWEDEN & MARTINA	PREMIUM KOHNO	4.25	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	PREMIUM KOHNO	5.00	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	PREMIUM KOHNO	6.00	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SWEDEN & MARTINA	STARK	REGULAR	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
SYBRON IMPLANT SOLUTIONS	ENDPORE (INNOVA)	4.10	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	ATTRACT	PLATAFORMA SWITCH	41.316.084.01-2	1.60	20 N·cm	8.40	3.50	6.80	-	1.60	2.30	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	CLASSIC-IN CR/ESTETIC-IN CR	3.60	41.316.084.01-2	1.60	20 N·cm	8.40	3.50	6.80	-	1.60	2.30	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	FIT-HEX	3.30	41.318.044.01-2	1.80	25 N·cm	4.40	2.75	3.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	CLASSIC-IN CR/ESTETIC-IN CR	4.30	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	CLASSIC-IN CR/ESTETIC-IN CR	5.00	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	CLASSIC TRADICIONAL	4.10	41.320.065.01-2	2.00	25 N·cm	6.50	2.70	5.00	-	1.50	2.30	straight	-	45° chamfer	Hex. 1.70
SYSTHEX	CLASSIC-CI/ESTETIC-CI	4.10	41.320.074.01-2	2.00	25 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	CLASSIC TRADICIONAL	5.00	41.325.054.01-2	2.50	25 N·cm	5.40	4.60	4.10	-	1.30	2.87	straight	-	semi-sphere	Hex. 1.70
SYSTHEX	CLASSIC-CI/ESTETIC-CI	5.00	41.325.067.01-2	2.50	25 N·cm	6.70	4.60	5.10	-	1.60	2.85	straight	-	semi-sphere	Hex. 1.70
TAV DENTAL	SIRIUS	NP	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
TAV DENTAL	SILHOUETTE	NP	41.316.076.01-2	1.60	20 N·cm	7.60	3.50	6.10	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
TAV DENTAL	SIRIUS	SP	41.320.082.01-2	2.00	25 N·cm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
TAV DENTAL	SILHOUETTE	SP	41.320.082.01-2	2.00	25 N·cm	8.20	4.70	6.70	-	1.50	2.40	straight	-	semi-sphere	Hex. 1.70
TBR	Z1/M	5.00	41.316.084.01-2	1.60	20 N·cm	8.40	3.50	6.80	-	1.60	2.30	straight	-	semi-sphere	Hex. 1.70
TBR	EXTERNAL HEXAGON	NARROW	41.320.060.01-2	2.00	25 N·cm	6.00	2.70	4.50	-	1.50	2.30	straight	-	semi-sphere	Hex. 1.70
TBR	CONNECT/INFINITI	3.50	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
TBR	CONNECT/INFINITI	4.00	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
TBR	CONNECT/INFINITI	5.00	41.320.068.01-2	2.00	25 N·cm	6.80	4.40	5.30	5.40	1.40	2.30	conical	60°	45° chamfer	Hex. 1.70
TBR	Z1 INTERNAL OCTAGON	3.50	41.320.071.01-2	2.00	25 N·cm	7.10	4.00	5.11	5.55	1.55	2.30	conical	20°	45° chamfer	Hex. 1.70
TBR	Z1 INTERNAL OCTAGON	4.00	41.320.071.01-2	2.00	25 N·cm	7.10	4.00	5.11	5.55	1.55	2.30	conical	20°	45° chamfer	Hex. 1.70
TBR	Z1 INTERNAL OCTAGON	5.00	41.320.071.01-2	2.00	25 N·cm	7.10	4.00	5.11	5.55	1.55	2.30	conical	20°	45° chamfer	Hex. 1.70
THOMMEN	INTERNAL HEXAGON	3.50	41.314.067.01	#N/D	#N/D	#N/D	#N/D	#N/D	#N/D	#N/D	#N/D	#N/D	#N/D	#N/D	Hex. 1.70
THOMMEN	INTERNAL HEXAGON	4.00	41.316.072.01-2	1.60	20 N·cm	7.20	3.50	5.20	5.82	1.38	2.30	conical	30°	45° chamfer	Hex. 1.70
THOMMEN	INTERNAL HEXAGON	4.50	41.316.072.01-2	1.60	20 N·cm	7.20	3.50	5.20	5.82	1.38	2.30	conical	30°	45° chamfer	Hex. 1.70
THOMMEN	INTERNAL HEXAGON	5.00	41.316.072.01-2	1.60	20 N·cm	7.20	3.50	5.20	5.82	1.38	2.30	conical	30°	45° chamfer	Hex. 1.70
THOMMEN	INTERNAL HEXAGON	6.00	41.316.072.01-2	1.60	20 N·cm	7.20	3.50	5.20	5.82	1.38	2.30	conical	30°	45° chamfer	Hex. 1.70
TITANIUM FIX	b-FIX	REGULAR	41.316.071.01-2	1.60	20 N·cm	7.10	2.80	5.20	5.54	1.56	2.30	conical	60°	45° chamfer	Hex. 1.70
TITANIUM FIX	i-FIX	REGULAR	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70
TITANIUM FIX	i-FIX	LARGA	41.318.075.01-2	1.80	25 N·cm	7.50	3.30	6.10	-	1.40	2.30	straight	-	semi-sphere	Hex. 1.70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
TITANIUM FIX	e-FIX EXT HEXAGON TORQUE INTERNO	REGULAR 4,1	41.320.068.01-2	2,00	25 N·cm	6,80	4,40	5,30	5,40	1,40	2,30	conical	60°	45° chamfer	Hex. 1,70
TITANIUM FIX	e-FIX EXTERNAL HEXAGON	REGULAR 4,1	41.320.068.01-2	2,00	25 N·cm	6,80	4,40	5,30	5,40	1,40	2,30	conical	60°	45° chamfer	Hex. 1,70
TITANIUM FIX	e-FIX SILVER	REGULAR 4,1	41.320.068.01-2	2,00	25 N·cm	6,80	4,40	5,30	5,40	1,40	2,30	conical	60°	45° chamfer	Hex. 1,70
TITANIUM FIX	c-FIX	REGULAR	41.320.074.01-2	2,00	25 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
TITANIUM FIX	b-FIX	LARGA	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
TITANIUM FIX	e-FIX SILVER	WIDE 5	41.325.054.01-2	2,50	25 N·cm	5,40	4,60	4,10	-	1,30	2,87	straight	-	semi-sphere	Hex. 1,70
TITANIUM FIX	e-FIX EXTERNAL HEXAGON	WIDE 5,1	41.325.067.01-2	2,50	25 N·cm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
TREE-OSS	MULTI UNIT	UNIVERSAL	41.314.039.01-2	1,40	15 N·cm	3,90	1,80	2,10	-	1,80	2,40	straight	-	45° chamfer	Hex. 1,70
TREE-OSS	EXTERNAL HEXAGON	3,50	41.316.059.01-2	1,60	20 N·cm	5,90	3,00	4,40	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
TREE-OSS	SIMPLE	3,75	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
TREE-OSS	TRI-CHANNEL INTERNO	3,50	41.318.075.01-2	1,80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
TREE-OSS	EXTERNAL HEXAGON	4,10	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
TREE-OSS	TRI-CHANNEL INTERNO	4,30	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
TREE-OSS	TRI-CHANNEL INTERNO	5,00	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
TREE-OSS	EXTERNAL HEXAGON	5,10	41.325.067.01-2	2,50	25 N·cm	6,70	4,60	5,10	-	1,60	2,85	straight	-	semi-sphere	Hex. 1,70
TRI IMPLANTS	TRI-VENT	3,50	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
TRI IMPLANTS	MULTI UNIT	UNIVERSAL	41.318.044.01-2	1,80	25 N·cm	4,40	2,75	3,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
TRI IMPLANTS	TRI-NARROW	3,20	41.318.080.01-2	1,80	25 N·cm	8,00	4,00	6,50	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
TRINON	Q2	4,00	41.320.060.01-2	2,00	25 N·cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
TRINON	QK	4,80	41.320.067.01-2	2,00	25 N·cm	6,70	2,30	3,59	5,70	1,00	2,58	conical	15°	45° chamfer	Hex. 1,70
UFIT	GT2	MINI	41.316.076.01-2	1,60	20 N·cm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
UFIT	GT2	MINI	41.316.076.01-2	1,60	20 N·cm	7,60	3,50	6,10	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
UFIT	GT2	REGULAR	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
UFIT	GT2	WIDE	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
UFIT	GT2	ULTRA WIDE	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
UFIT	NT2	REGULAR	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
UFIT	NT2	WIDE	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
UFIT	NT2	ULTRA WIDE	41.320.090.01-2	2,00	25 N·cm	9,00	4,00	7,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
VULKAN IMPLANTS	IN-HEX	3,75	41.317.071.01-2	N1-72	25 N·cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
WINSIX	TTi	PLATFORM SWITCHING	41.318.075.01-2	1,80	25 N·cm	7,50	3,30	6,10	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
XIVE	XIVE	3,00	41.314.076.01-2	1,40	15 N·cm	7,60	2,40	5,88	6,35	1,25	2,30	conical	45°	45° chamfer	Hex. 1,70
XIVE	MULTI UNIT	3,40	41.316.044.01-2	1,60	20 N·cm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
XIVE	MULTI UNIT	3,80	41.316.044.01-2	1,60	20 N·cm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
XIVE	MULTI UNIT	4,50	41.316.044.01-2	1,60	20 N·cm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
XIVE	MULTI UNIT	5,50	41.316.044.01-2	1,60	20 N·cm	4,40	2,50	2,90	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
XIVE	XIVE	3,40	41.316.081.01-2	1,60	20 N·cm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70

DYNAMIC SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
XIVE	XIVE	3,80	41.316.081.01-2	1,60	20 N.cm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
XIVE	XIVE	4,50	41.316.081.01-2	1,60	20 N.cm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
XIVE	XIVE	5,50	41.316.081.01-2	1,60	20 N.cm	8,10	3,00	6,35	6,73	1,37	2,30	conical	45°	45° chamfer	Hex. 1,70
YES IMPLANTS	S-SYSTEM	NARROW	41.320.079.01-2	2,00	25 N.cm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
YES IMPLANTS	S-SYSTEM	REGULAR	41.320.079.01-2	2,00	25 N.cm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
YES IMPLANTS	S-SYSTEM	WIDE	41.320.079.01-2	2,00	25 N.cm	7,90	3,30	6,32	6,50	1,40	2,30	conical	45°	45° chamfer	Hex. 1,70
ZIACOM (OSSEOLIFE)	OEX - EXTERNAL HEXAGON	NP 3,3	41.318.065.01-2	1,80	25 N.cm	6,50	2,83	5,00	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ZIACOM (OSSEOLIFE)	OEX - EXTERNAL HEXAGON	RP 4,1	41.320.060.01-2	2,00	25 N.cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ZIACOM (OSSEOLIFE)	OEX - EXTERNAL HEXAGON	WP 5	41.320.060.01-2	2,00	25 N.cm	6,00	2,70	4,50	-	1,50	2,30	straight	-	semi-sphere	Hex. 1,70
ZIACOM (OSSEOLIFE)	ION - INTERNAL HEXAGON	RP 4	41.320.074.01-2	2,00	25 N.cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ZIACOM (OSSEOLIFE)	ION - INTERNAL HEXAGON	WP 5	41.320.074.01-2	2,00	25 N.cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70
ZIMMER	EZTETIC	NP 2,9	41.316.080.01-2	1,60	20 N.cm	8,00	3,00	6,30	6,51	1,49	2,30	conical	60°	45° chamfer	Hex. 1,70
ZIMMER	TAPERED ABUTMENT MULTI-UNIT	UNIVERSAL	41.317.040.01-2	N1-72	25 N.cm	4,00	2,10	2,50	-	1,50	2,30	straight	-	45° chamfer	Hex. 1,70
ZIMMER	SCREW-VENT	3,50	41.317.071.01-2	N1-72	25 N.cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ZIMMER	SCREW-VENT	4,50	41.317.071.01-2	N1-72	25 N.cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ZIMMER	SCREW-VENT	5,70	41.317.071.01-2	N1-72	25 N.cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ZIMMER	SWISS PLUS	3,80	41.317.071.01-2	N1-72	25 N.cm	7,10	2,60	5,54	5,65	1,45	2,30	conical	70°	45° chamfer	Hex. 1,70
ZIMMER	SWISS PLUS	4,80	41.320.074.01-2	2,00	25 N.cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Hex. 1,70

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
AB	l22	STANDARD	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
AB	l5	STANDARD	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
AB	l55	STANDARD	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
AB	l10	STANDARD	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
AB	l15	STANDARD	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
AB	MULTI-UNIT D1-P14	CONICAL	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
AB	MULTI-UNIT D1-P64	CONICAL	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
AB	l2	STANDARD	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ACE	INFINITY EXTERNAL HEX	NP 3.5	40.316.008.01-2	1.60	20 Ncm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
ACE	INFINITY INTERNAL HEX	4.50	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ACE	MULTI-UNIT	CONICAL	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
ACE	INFINITY INTERNAL HEX	3.50	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	SWELL INTERNAL HEX	3.45	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	SWELL INTERNAL HEX	3.60	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	SWELL INTERNAL HEX	4.00	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	SWELL INTERNAL HEX	4.60	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	TOUAREG-S/TOUAREG-OS	3.45	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	TOUAREG-S/TOUAREG-OS	3.60	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	TOUAREG-S/TOUAREG-OS	4.00	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	TOUAREG-S/TOUAREG-OS	5.00	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ADIN	TOUAREG CLOSEFIT	UNP	40.315.008.01-2	No-80	15 Ncm	8,30	2,45	5,79	5,95	2,35	2,00	conical	65°	chamfer 45°	Unigrip
ADIN	TOUAREG CLOSEFIT	NP	40.315.008.01-2	No-80	15 Ncm	8,30	2,45	5,79	5,95	2,35	2,00	conical	65°	chamfer 45°	Unigrip
ADIN	TOUAREG CLOSEFIT	RP	40.316.008.02-2	1.60	20 Ncm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
ADIN	MULTI-UNIT	CONICAL	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
ADVAN	GFA	UNIVERSAL GTB	40.316.008.02-2	1.60	20 Ncm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
ALFAGATE	MULTI-UNIT	UNIVERSAL	40.317.005.03-2	N1-72	25 Ncm	8,50	1,60	6,58	7,00	1,50	2,20	conical	45°	chamfer 45°	Hex. 1.27
ALFAGATE	Bioactive S-line/SLA	3.75	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ALFAGATE	CONICA HEXÁGONO INTERNO	3.50	40.316.008.02-2	1.60	20 Ncm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
ALPHA BIO	DFI	UNIVERSAL	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ALPHA BIO	ATID	UNIVERSAL	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ALPHA BIO	NEO/NICE	3.50	40.316.004.03-2	1.60	20 Ncm	6,90	2,20	5,02	5,20	1,70	1,92	conical	60°	chamfer 45°	Hex. 1.25
ALPHA BIO	CONICAL STANDARD (CS)	STANDARD	40.317.004.02-2	N1-72	25 Ncm	7,52	2,20	5,12	5,77	1,75	2,10	conical	30°	chamfer 45°	Hex. 1.25
ALPHA BIO	SPI	UNIVERSAL	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ALPHA BIO	ICE	UNIVERSAL	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ANCLADEN	ANCLALOCK	3.50	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ANTHOGRYR	AXIOM REG/PX	3.40	40.316.014.01-2	1.60	20 Ncm	7,90	2,30	5,42	6,46	1,44	2,20	conical	20°	chamfer 45°	Hex. 1.19
ANTHOGRYR	AXIOM REG/PX	4.00	40.316.014.01-2	1.60	20 Ncm	7,90	2,30	5,42	6,46	1,44	2,20	conical	20°	chamfer 45°	Hex. 1.19

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
ANTHOGYR	AXIOM REG/PX	4.60	40.316.014.01-2	1.60	20 Ncm	7.90	2.30	5.42	6.46	1.44	2.20	conical	20°	chamfer 45°	Hex. 1.19
ANTHOGYR	AXIOM REG/PX	5.20	40.316.014.01-2	1.60	20 Ncm	7.90	2.30	5.42	6.46	1.44	2.20	conical	20°	chamfer 45°	Hex. 1.19
ANTHOGYR	MULTI UNIT	UNIVERSAL	40.314.014.01-2	1.40	15 Ncm	4.45	2.30	2.48	-	1.97	2.16	straight	-	chamfer 45°	Hex. 1.19
ARDS	SMART	3.75	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
ARDS	CLASSIC	3.75	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
ARDS	PREMIUM	3.75	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
ASTRA	CONE 20°	REGULAR/WIDE	40.314.003.02-2	1.40	15 Ncm	4.00	2.00	2.20	-	1.80	2.30	straight	-	chamfer 45°	Hex. 1.20
ASTRA	CONE 45°	REGULAR/WIDE	40.314.003.02-2	1.40	15 Ncm	4.00	2.00	2.20	-	1.80	2.30	straight	-	chamfer 45°	Hex. 1.20
ASTRA	YELLOW	NP	40.314.005.02-2	1.40	15 Ncm	7.50	2.40	5.50	5.70	1.80	2.10	conical	60°	chamfer 45°	Hex. 1.27
ASTRA	AQUA	RP	40.316.005.02-2	1.60	20 Ncm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
ASTRA	EVOLUTION	3.0	40.314.005.01-2	1.40	15 Ncm	7.60	3.60	5.21	6.00	1.60	2.15	conical	25°	chamfer 45°	Hex. 1.27
ASTRA	EVOLUTION	3.60	40.316.005.01-2	1.60	20 Ncm	7.50	3.60	5.44	5.90	1.60	2.13	conical	31°	chamfer 45°	Hex. 1.27
ASTRA	EVOLUTION	4.20	40.318.005.02-2	1.80	25 Ncm	7.60	3.60	5.76	6.00	1.60	2.33	conical	50°	chamfer 45°	Hex. 1.27
ASTRA	EVOLUTION EXTERNAL CONE	UNIVERSAL	40.318.005.01-2	1.80	25 Ncm	4.50	2.30	2.80	2.90	1.60	2.33	conical	70°	chamfer 45°	Hex. 1.27
AVINENT	HI	3.50	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
AVINENT	HI	4.10	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
AVINENT	TRANSEPIITAL	REGULAR	40.314.003.02-2	1.40	15 Ncm	4.00	2.00	2.20	-	1.80	2.30	straight	-	chamfer 45°	Hex. 1.20
BEGO	EXTERNAL HEXAGON	NARROW MINI (2.7/2.9/3.1)	40.316.008.01-2	1.60	20 Ncm	7.00	2.70	5.15	-	1.85	2.30	straight	-	chamfer 45°	Unigrip
BEGO	RS/RSX	3.0	40.314.004.01-2	1.40	15 Ncm	6.30	2.20	4.60	5.10	1.20	2.10	conical	25°	chamfer 30°	Hex. 1.25
BEGO	S/RI/RS/RSX	3.25/3.75	40.318.004.03-2	1.80	25 Ncm	7.65	3.30	5.17	5.75	1.90	2.40	conical	35°	chamfer 45°	Hex. 1.25
BEGO	S/RI/RS/RSX	4.10	40.318.004.03-2	1.80	25 Ncm	7.65	3.30	5.17	5.75	1.90	2.40	conical	35°	chamfer 45°	Hex. 1.25
BEGO	S/RI/RS/RSX	4.50	40.318.004.03-2	1.80	25 Ncm	7.65	3.30	5.17	5.75	1.90	2.40	conical	35°	chamfer 45°	Hex. 1.25
BEGO	S/RI/RS/RSX	5.50	40.318.004.03-2	1.80	25 Ncm	7.65	3.30	5.17	5.75	1.90	2.40	conical	35°	chamfer 45°	Hex. 1.25
BEGO	MULTI PLUS	UNIVERSAL	40.314.004.04-2	1.40	20 Ncm	4.60	2.45	2.36	2.95	1.65	2.10	conical	35°	chamfer 45°	Hex. 1.25
BIOCONCEPT	BC BONE LEVEL	NARROW	40.316.007.01-2	1.60	20 Ncm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
BIOCONCEPT	BC BONE LEVEL	REGULAR	40.316.007.01-2	1.60	20 Ncm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
BIOCONCEPT	BV TAPERED BONE LEVEL	NARROW	40.316.003.02-2	1.60	20 Ncm	10.20	2.00	7.88	8.20	2.00	2.20	conical	45°	chamfer 45°	Hex. 1.20
BIOCONCEPT	BV TAPERED BONE LEVEL	MULTI UNIT	40.314.008.01-2	1.40	15 Ncm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
BIOGENESIS	3ICON	MINI (PINK)	40.316.008.01-2	1.60	20 Ncm	7.00	2.70	5.15	-	1.85	2.30	straight	-	chamfer 45°	Unigrip
BIOGENESIS	ATICON (CONO 20°/45°)	REGULAR/WIDE	40.314.003.02-2	1.40	15 Ncm	4.00	2.00	2.20	-	1.80	2.30	straight	-	chamfer 45°	Hex. 1.20
BIOHORIZONS	EXTERNAL	3.7 (YELLOW)	40.318.003.01-2	1.80	25 Ncm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
BIOHORIZONS	INTERNAL	3 (GREY)	40.317.005.02-2	N1-72	25 Ncm	7.30	2.40	4.73	5.25	2.05	2.40	conical	45°	chamfer 45°	Hex. 1.27
BIOHORIZONS	INTERNAL	3.5 (YELLOW)	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
BIOHORIZONS	INTERNAL	4.5 (GREEN)	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
BIOHORIZONS	INTERNAL	5.7 (BLUE)	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
BIOHORIZONS	TAPERED TISSUE LEVEL	3.5 (YELLOW)	40.317.004.01-2	N1-72	25 Ncm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
BIOHORIZONS	TAPERED TISSUE LEVEL	4.5 (GREEN)	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIOHORIZONS	TAPERED TISSUE LEVEL	5.7 (BLUE)	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIOHORIZONS	MULTI-UNIT		40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
BIOMET 3i	CERTAIN	3.40	40.316.003.01-2	1.60	20 Ncm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
BIOMET 3i	CERTAIN	4.10	40.316.003.01-2	1.60	20 Ncm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
BIOMET 3i	CERTAIN	5.00	40.316.003.01-2	1.60	20 Ncm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
BIOMET 3i	LOW PROFILE	UNIVERSAL	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
BIONER	HIKELT	3.95	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIONER	HIKELT	4.90	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIONER	TRANSEPITELIAL	REGULAR	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
BIONER	TOP-DM	3.50	40.316.008.02-2	1.60	20 Ncm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
BIONER	TOP-DM	4.00	40.316.008.02-2	1.60	20 Ncm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
BIONER	TOP-DM	5.00	40.316.008.02-2	1.60	20 Ncm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
BIONER	TRANSEPITELIAL TOP DM	3.50	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
BIONER	TRANSEPITELIAL A-5M	REGULAR	40.314.008.01-2	1.40	15 Ncm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
BIONNOVATION	EXTERNAL HEXAGON	SP	40.318.003.01-2	1.80	25 Ncm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
BIONNOVATION	INTERNAL HEXAGON	SP	40.316.003.01-2	1.60	20 Ncm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
BIOTEC	SPR/CIM	3.30	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIOTEC	SPR/CIM	3.75	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIOTEC	SPR/SPTT/CIM	4.20	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIOTEC	SPR/SPTT/CIM	5.00	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIOTEC	SPR/SPTT/CIM	6.00	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BIOTECH	KONTACT	3.00	40.312.003.01-2	1.20	15 Ncm	7,85	2,70	6,20	6,55	1,30	1,90	conical	45°	chamfer 45°	Hex. 1.20
BIOTECH	KONTACT	3.6/4.2/4.8/5.4	40.314.003.03-2	1.40	15 Ncm	7,60	2,40	6,05	6,30	1,30	1,90	conical	45°	chamfer 45°	Hex. 1.20
BIOTECH	MULTI UNIT	UNIVERSAL	40.314.012.01-2	1.40	15 Ncm	4,50	1,70	2,01	2,40	2,10	2,15	conical	45°	chamfer 45°	Star 1.50
BIOTEM	AR FIXTURE R-TYPE/T-TYPE	MINI	40.316.005.02-2	1.60	20 Ncm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
BIOTEM	BR FIXTURE	MINI	40.316.008.01-2	1.60	20 Ncm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
BLUE SKY BIO	INTERNAL HEX	NARROW	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BLUE SKY BIO	INTERNAL HEX	REGULAR	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BLUE SKY BIO	INTERNAL HEX	WIDE	40.317.004.01-2	N1-72	25 Ncm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
BLUE SKY BIO	TRILOBE	3.50	40.318.008.01-2	1.80	25 Ncm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
BLUE SKY BIO	MAX	NP	40.316.005.02-2	1.60	20 Ncm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
BREDENT MEDICAL	NARROW SKY	NP 3.5	40.318.007.01-2	1.80	25 Ncm	9,10	4,25	7,32	7,45	1,65	2,18	conical	60°	chamfer 45°	Torx T6
BREDENT MEDICAL	BLUE SKY / BLUE SKY CLASSIC	4.00	40.318.007.01-2	1.80	25 Ncm	9,10	4,25	7,32	7,45	1,65	2,18	conical	60°	chamfer 45°	Torx T6
BREDENT MEDICAL	COPA SKY	4.00	40.316.007.01-2	1.60	20 Ncm	7,90	2,00	5,63	6,90	1,00	2,18	conical	15°	chamfer 45°	Torx T6

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
BTI	CONEXIÓN MULTI-IM	UNIVERSAL 4.1	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
BTI	CONEXIÓN MULTI-IM	ANCHA 5.5	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
BTI	CONEXIÓN EXTERNA	TINY 3.5	40.320.003.01-2	2.00	30 N·cm	4.70	3.00	3.33	-	1.37	2.35	straight	-	chamfer 45°	Hex. 1.20
BTI	CONEXIÓN INTERNA	ESTRECHA 3.5	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
BTI	CONEXIÓN INTERNA	UNIVERSAL 4.1	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
BTI	CONEXIÓN INTERNA	UNIVERSAL 4.1	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
BTI	CONEXIÓN INTERNA	ANCHA 5.5	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
BTI	MULTI-IM	UNIVERSAL 4.1	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
BTK	KLASSIC/KONIC EXTERNAL	3.4 EN	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
BTK	KLASSIC INTERNAL	3.5 IR	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
BTK	KLASSIC INTERNAL	4.5 IM	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
BTK	KLASSIC INTERNAL	5 IW	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
BTK	KONIC INTERNAL	3.5 IR	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
BTK	KONIC INTERNAL	4.5 IM	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
BTK	KONIC INTERNAL	5 IW	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
BTK	SAFE TL	REGULAR	40.320.007.01-2	2.00	30 N·cm	6.70	2.25	3.59	5.70	1.00	2.58	conical	15°	chamfer 45°	Torx T6
BW	EXTERNAL HEXAGON	3.30	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
BW	CONICAL INTERNAL HEXAGON	4.00	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
CAMLOG	SCREW LINE	3.30	40.316.005.04-2	1.60	20 N·cm	10.50	3.20	8.15	8.40	2.10	2.10	conical	45°	chamfer 45°	Hex. 1.27
CAMLOG	SCREW LINE	3.80	40.316.005.04-2	1.60	20 N·cm	10.50	3.20	8.15	8.40	2.10	2.10	conical	45°	chamfer 45°	Hex. 1.27
CAMLOG	SCREW LINE	4.30	40.316.005.04-2	1.60	20 N·cm	10.50	3.20	8.15	8.40	2.10	2.10	conical	45°	chamfer 45°	Hex. 1.27
CAMLOG	SCREW LINE	5.00	40.320.005.04-2	2.00	30 N·cm	10.50	3.20	8.15	8.40	2.10	2.50	conical	45°	chamfer 45°	Hex. 1.27
CAMLOG	SCREW LINE	6.00	40.320.005.04-2	2.00	30 N·cm	10.50	3.20	8.15	8.40	2.10	2.50	conical	45°	chamfer 45°	Hex. 1.27
CAMLOG	CONELOG SCREW LINE	3.80	40.316.005.07-2	1.60	20 N·cm	8.80	2.85	6.73	6.90	1.90	2.15	conical	60°	chamfer 45°	Hex. 1.27
CAMLOG	CONELOG SCREW LINE	4.30	40.316.005.07-2	1.60	20 N·cm	8.80	2.85	6.73	6.90	1.90	2.15	conical	60°	chamfer 45°	Hex. 1.27
CAMLOG	COMFOUR	5/6	40.320.003.01-2	2.00	30 N·cm	4.70	3.00	3.33	-	1.37	2.35	straight	-	chamfer 45°	Hex. 1.20
CAMLOG	ISY		40.316.005.04-2	1.60	20 N·cm	10.50	3.20	8.15	8.40	2.10	2.10	conical	45°	chamfer 45°	Hex. 1.27
CONEXAO	AR MORSE	UNIVERSAL	40.316.008.02-2	1.60	20 N·cm	7.30	2.70	5.15	5.90	1.40	2.20	conical	35°	chamfer 45°	Unigrip
CONEXAO	AR MORSE	UNIVERSAL	40.316.008.02-2	1.60	20 N·cm	7.30	2.70	5.15	5.90	1.40	2.20	conical	35°	chamfer 45°	Unigrip
CONEXAO	GRIP ACTIVES	UNIVERSAL	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
CONEXAO	GRIP ACTIVES	UNIVERSAL	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
CONEXAO	GRIP POROUS	UNIVERSAL	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
CONEXAO	GRIP POROUS	UNIVERSAL	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
CONEXAO	FLASH	UNIVERSAL	40.316.008.02-2	1.60	20 N·cm	7.30	2.70	5.15	5.90	1.40	2.20	conical	35°	chamfer 45°	Unigrip
CONEXAO	TORQ	UNIVERSAL	40.316.008.02-2	1.60	20 N·cm	7.30	2.70	5.15	5.90	1.40	2.20	conical	35°	chamfer 45°	Unigrip
CONEXAO	EXPAND	UNIVERSAL	40.316.008.02-2	1.60	20 N·cm	7.30	2.70	5.15	5.90	1.40	2.20	conical	35°	chamfer 45°	Unigrip
CORTEX	CONICAL PLATFORM	NP	40.314.005.02-2	1.40	15 N·cm	7.50	2.40	5.50	5.70	1.80	2.10	conical	60°	chamfer 45°	Hex. 1.27

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
CORTEX	CONICAL PLATFORM	RP	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
CORTEX	CONICAL PLATFORM	WP	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
CORTEX	MULTI UNIT	UNIVERSAL	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
CORTEX	SATURN	3.50	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
CORTEX	CLASSIX/DYNAMIX	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
CORTEX	MAGIX	RP	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
C-TECH	MULTI UNIT	UNIVERSAL	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
DENTAL IMPLANT KA	BONE LEVEL	NP	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
DENTAL IMPLANT KA	TISSUE LEVEL	RP	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
DENTAL IMPLANT KA	TISSUE LEVEL	WP	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
DENTAURUM	TIOLOGIC	S	40.316.005.08-2	1.60	20 N·cm	9,00	3,90	6,49	6,90	2,10	2,18	conical	45°	chamfer 45°	Hex. 1.27
DENTAURUM	TIOLOGIC	M	40.316.005.08-2	1.60	20 N·cm	9,00	3,90	6,49	6,90	2,10	2,18	conical	45°	chamfer 45°	Hex. 1.27
DENTAURUM	TIOLOGIC	L	40.316.005.08-2	1.60	20 N·cm	9,00	3,90	6,49	6,90	2,10	2,18	conical	45°	chamfer 45°	Hex. 1.27
DENTIUM	NR LINE	3.60	40.316.003.01-2	1.60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
DENTIUM	NR LINE	4.30	40.316.003.01-2	1.60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
DENTIUM	MULTI-UNIT NR LINE	4.30	40.316.003.01-2	1.60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
DENTIUM	MULTI-UNIT NR LINE	5.50	40.316.003.01-2	1.60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
DENTIUM	MULTI-UNIT SUPER LINE AND IMPLANTUM	3.6/4/4.5/5/6	40.318.005.01-2	1,80	25 N·cm	4,50	2,30	2,80	2,90	1,60	2,33	conical	70°	chamfer 45°	Hex. 1.27
DENTOFLEX	SUPREME MORSE CONE	3.30	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
DENTOFLEX	SUPREME MORSE CONE	4.50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
DIO	EXTERNAL	NARROW 3.5	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
DIO	SM SYSTEM	NARROW	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
DIO	SM SYSTEM	REGULAR/WIDE	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
DIO	UF	NARROW	40.314.003.04-2	1,40	15 N·cm	7,50	2,50	5,45	5,70	1,80	1,85	conical	45°	chamfer 45°	Hex. 1.20
DIO	UF	REGULAR	40.320.003.04-2	2,00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20
DSP BIOMEDICAL	EXTERNAL HEXAGON	3.30	40.318.003.01-2	1,80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
DSP BIOMEDICAL	EXTERNAL HEXAGON	4.10	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
DSP BIOMEDICAL	EXTERNAL HEXAGON	5.00	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
DSP BIOMEDICAL	HEXÁGONO CONE MORSE	3.30	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
DSP BIOMEDICAL	HEXÁGONO CONE MORSE	3.75	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
DSP BIOMEDICAL	HEXÁGONO CONE MORSE	4.00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
DSP BIOMEDICAL	HEXÁGONO CONE MORSE	5.00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
DSP BIOMEDICAL	CRESTAL LEVEL	3.50	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
EASY IMPLANT	MASTER C	4.50	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	MASTER C	5.00	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
EASY IMPLANT	MASTER L	3,30	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	MASTER L	3,75	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	MASTER L	4,25	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	MASTER L	4,75	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	MASTER S	3,75	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	MASTER S	4,25	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	MASTER S	4,75	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EASY IMPLANT	HEXCEL-S	3,30	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
EASY IMPLANT	HEXCEL-S	4,10	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
EASY IMPLANT	HEXCEL-S	5,00	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
EASY IMPLANT	MULTI-UNIT CONICAL ABUTMENT	UNIVERSAL	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
ECKERMANN	ALL-SPIRAL	REGULAR	40.318.004.01-2	1,80	25 N·cm	7,20	4,80	5,36	5,90	1,30	2,40	conical	30°	chamfer 45°	Hex. 1.25
ECKERMANN	TRIPLO	4,00	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
ELITE MEDICAL	EXTERNAL HEXAGON	NARROW	40.316.008.01-2	1,60	20 N·cm	7,00	2,70	5,45	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
EMFILS	NOVO COLOSSO MORSE	NARROW	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
EMFILS	NOVO COLOSSO MORSE	REGULAR	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
EMFILS	NOVO COLOSSO MORSE	WIDE	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
EMFILS	COLOSSO ECI/CRC EVOLUTION	NARROW	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
EMFILS	COLOSSO ECI/CRC EVOLUTION	REGULAR	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
EMFILS	COLOSSO ECI/CRC EVOLUTION	WIDE	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	NARROW	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	REGULAR	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
EMFILS	NOVO COLOSSO EXTERNAL HEXAGON	WIDE	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
EQUINOX	UNITI	3,30	40.314.003.03-2	1,40	15 N·cm	7,60	2,40	6,05	6,30	1,30	1,90	conical	45°	chamfer 45°	Hex. 1.20
EQUINOX	UNITI	3,70	40.316.004.02-2	1,60	20 N·cm	8,80	3,00	6,73	7,00	1,80	2,10	conical	45°	chamfer 45°	Hex. 1.25
EQUINOX	UNITI	4,30	40.316.004.02-2	1,60	20 N·cm	8,80	3,00	6,73	7,00	1,80	2,10	conical	45°	chamfer 45°	Hex. 1.25
EQUINOX	UNITI	5,60	40.316.004.02-2	1,60	20 N·cm	8,80	3,00	6,73	7,00	1,80	2,10	conical	45°	chamfer 45°	Hex. 1.25
EUROTEKNIKA	AESTHETICA	6,50	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
EUROTEKNIKA	NATURACTIS/NATURALL	3,00	40.314.005.01-2	1,40	15 N·cm	7,60	3,60	5,21	6,00	1,60	2,15	conical	25°	chamfer 45°	Hex. 1.27
EUROTEKNIKA	NATURACTIS	3,40	40.316.005.02-2	1,60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EUROTEKNIKA	NATURACTIS	3,80	40.316.005.02-2	1,60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EUROTEKNIKA	NATURACTIS	4,30	40.316.005.02-2	1,60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EUROTEKNIKA	NATURACTIS	4,80	40.316.005.02-2	1,60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EUROTEKNIKA	NATEA/NATURALL	NARROW	40.316.005.02-2	1,60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EUROTEKNIKA	NATEA/NATURALL	REGULAR	40.316.005.02-2	1,60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EUROTEKNIKA	NATEA/NATURALL	WIDE	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
EUROTEKNIKA	MULTI-UNIT TETRA	UNIVERSAL	40.320.003.04-2	2,00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
EUROTEKNIKA	MULTI-UNIT ANGULATED CONICAL	UNIVERSAL	40.320.003.04-2	2.00	30 N·cm	8.40	2.75	5.68	6.35	2.05	2.30	conical	45°	chamfer 45°	Hex. 120
EVO IMPLANTS	EVOSOFT / EVOPLUS	HEX. INT. 3.8/4.2/4.8/5.5	40.316.005.04-2	1.60	20 N·cm	10.50	3.20	8.15	8.40	2.10	2.10	conical	45°	chamfer 45°	Hex. 127
EVO IMPLANTS	EVOSHORT	HEX. INT. 4.8/5.5/6.2	40.316.005.04-2	1.60	20 N·cm	10.50	3.20	8.15	8.40	2.10	2.10	conical	45°	chamfer 45°	Hex. 127
GALIMPLANT	CONEXIÓN INTERNA	3.50	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 127
GALIMPLANT	CONEXIÓN INTERNA	4.00	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 127
GALIMPLANT	CONEXIÓN INTERNA	5.00	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 127
GALIMPLANT	PILAR MULTI-POSICION RECTO	UNIVERSAL	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
GALIMPLANT	PILAR MULTI-POSICION ANGULADO	UNIVERSAL	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
GALIMPLANT	PILAR MULTI-POSICION RECTO-ESTÉTICO	UNIVERSAL	40.320.003.01-2	2.00	30 N·cm	4.70	3.00	3.33	-	1.37	2.35	straight	-	chamfer 45°	Hex. 120
GC TECH	AADVA STANDARD / TAPERED IMPLANTS	NARROW	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 120
GC TECH	AADVA STANDARD / TAPERED IMPLANTS	REGULAR	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 120
GC TECH	AADVA STANDARD / TAPERED IMPLANTS	WIDE	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 120
GLOBAL D (TEKKA)	IN-KONE UNIVERSAL	5.00	40.318.003.03-2	1.80	25 N·cm	8.50	1.60	6.56	7.00	1.50	2.20	conical	45°	chamfer 45°	Hex. 120
GLOBAL D (TEKKA)	IN-KONE PRIMO	5.00	40.318.003.03-2	1.80	25 N·cm	8.50	1.60	6.56	7.00	1.50	2.20	conical	45°	chamfer 45°	Hex. 120
GLOBAL D (TEKKA)	IN-KONE MULTI UNIT	4.30	40.318.006.01-2	1.80	25 N·cm	6.00	3.00	3.67	4.00	2.00	2.40	conical	45°	chamfer 45°	Hex. 170
GMI	AVANTGARD	RP	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 120
GMI	FRONTIER	RP 3.3	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 120
GMI	FRONTIER	WP 4.3	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 120
GMI	PHOENIX	NARROW 3.3	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 120
GMI	UNIVERSAL	PS-RP 4.8	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
GT MEDICAL	BEST-FIT EXTERNAL HEXAGON	NARROW	40.316.008.01-2	1.60	20 N·cm	7.00	2.70	5.15	-	1.85	2.30	straight	-	chamfer 45°	Unigrip
GT MEDICAL	BEST-FIT INTERNAL HEXAGON	WIDE	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 127
HAHN IMPLANT	TAPERED IMPLANT	3 (Grey)	40.314.005.01-2	1.40	15 N·cm	7.60	3.60	5.21	6.00	1.60	2.15	conical	25°	chamfer 45°	Hex. 127
HAHN IMPLANT	TAPERED IMPLANT	3.5 (Purple)	40.316.008.02-2	1.60	20 N·cm	7.30	2.70	5.15	5.90	1.40	2.20	conical	35°	chamfer 45°	Unigrip
HAHN IMPLANT	MULTI UNIT ABUTMENT SYSTEM	UNIVERSAL	40.314.005.01-2	1.40	15 N·cm	7.60	3.60	5.21	6.00	1.60	2.15	conical	25°	chamfer 45°	Hex. 127
HIOSSSEN	ETII SA / ETIII SA	MINI	40.316.003.02-2	1.60	20 N·cm	10.20	2.00	7.88	8.20	2.00	2.20	conical	45°	chamfer 45°	Hex. 120
HIOSSSEN	ETII SA / ETIII SA	REGULAR	40.320.003.04-2	2.00	30 N·cm	8.40	2.75	5.68	6.35	2.05	2.30	conical	45°	chamfer 45°	Hex. 120
HIOSSSEN	ETIII BA	MINI	40.316.003.02-2	1.60	20 N·cm	10.20	2.00	7.88	8.20	2.00	2.20	conical	45°	chamfer 45°	Hex. 120
HIOSSSEN	ETIII BA	REGULAR	40.320.003.04-2	2.00	30 N·cm	8.40	2.75	5.68	6.35	2.05	2.30	conical	45°	chamfer 45°	Hex. 120
HI-TEC	TAPERED SELF THREAD	3.50	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 127
HI-TEC	TAPERED SELF THREAD	4.50	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 127
HI-TEC	LOGIC PLUS	3.70	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 127
HI-TEC	LOGIC PLUS	3.90	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 127
IBS	MAGICORE	4/4.5	40.316.002.01-2	1.60	20 N·cm	7.00	2.79	4.86	5.44	1.56	2.30	conical	45°	chamfer 45°	Cuad. 130
IBS	MAGICORE	5.50	40.316.002.01-2	1.60	20 N·cm	7.00	2.79	4.86	5.44	1.56	2.30	conical	45°	chamfer 45°	Cuad. 130
IBS	MAGIC FC	3.80	40.320.003.04-2	2.00	30 N·cm	8.40	2.75	5.68	6.35	2.05	2.30	conical	45°	chamfer 45°	Hex. 120

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
IBS	N.R. FIX	3,80	40.320.003.04-2	2,00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20
IDO IMPLANT	I DO IMPLANT	UNIVERSAL	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
IHDE DENTAL (IMBIODENT)	BONE LEVEL PLUS	3,30	40.316.007.01-2	1,60	20 N·cm	7,90	2,00	5,63	6,90	1,00	2,18	conical	15°	chamfer 45°	Torx T6
IHDE DENTAL (IMBIODENT)	BONE LEVEL PLUS	4,10	40.316.007.01-2	1,60	20 N·cm	7,90	2,00	5,63	6,90	1,00	2,18	conical	15°	chamfer 45°	Torx T6
IHDE DENTAL (IMBIODENT)	BONE LEVEL PLUS	4,80	40.316.007.01-2	1,60	20 N·cm	7,90	2,00	5,63	6,90	1,00	2,18	conical	15°	chamfer 45°	Torx T6
IMPLACIL DE BORTOLI	EXTERNAL HEXAGON	3,50	40.316.008.01-2	1,60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
IMPLACIL DE BORTOLI	INTERNAL HEXAGON	3,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLACIL DE BORTOLI	INTERNAL HEXAGON	4,00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
IMPLACIL DE BORTOLI	INTERNAL HEXAGON	5,00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
IMPLACIL DE BORTOLI	CONE MORSE	3,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLACIL DE BORTOLI	CONE MORSE	4,00	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLACIL DE BORTOLI	CONE MORSE	4,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLACIL DE BORTOLI	CONE MORSE	5,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLANT DIRECT	REACTIVE	3,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLANT DIRECT	REPLANT	3,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLANT DIRECT	REPLUS	3,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
IMPLANT DIRECT	LEGACY	3,50	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
IMPLANT DIRECT	LEGACY	4,50	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
IMPLANT DIRECT	LEGACY	5,70	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
IMPLANT DIRECT	SWISHPLANT	4,80	40.320.007.01-2	2,00	30 N·cm	6,70	2,25	3,59	5,70	1,00	2,58	conical	15°	chamfer 45°	Torx T6
IMPLANT DIRECT	SWISHPLANT	6,50	40.320.007.01-2	2,00	30 N·cm	6,70	2,25	3,59	5,70	1,00	2,58	conical	15°	chamfer 45°	Torx T6
IMPLANT DIRECT	SWISH ACTIVE	3,00	40.316.008.02-2	1,60	20 N·cm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
IMPLANT DIRECT	INTERACTIVE	3,00	40.316.008.02-2	1,60	20 N·cm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
IMPLANT DIRECT	PPA	5,00	40.318.006.01-2	1,80	25 N·cm	6,00	3,00	3,67	4,00	2,00	2,40	conical	45°	chamfer 45°	Hex. 1.70
IMPLANT GENESIS	AKTIVE SYSTEM	STANDARD	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
IMPLASA HÖSCHT	DOUBLEGOTIC/GOTIC/NORDIC	3,50	40.318.003.01-2	1,80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
INTRA LOCK	INTRAHEX	3,75/4	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
INTRA LOCK	INTRAHEX	4,75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
INTRA LOCK	INTRAHEX MULTI UNIT	UNIVERSAL	40.316.008.01-2	1,60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
INTRA LOCK	UNIHEX	WIDE	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
JDENTALCARE	JDEVOLUTION/JDEVOLUTION PLUS	3,70	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
JDENTALCARE	JDEVOLUTION/JDEVOLUTION PLUS	4.00	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
JDENTALCARE	JDEVOLUTION/JDEVOLUTION PLUS	5.00	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
JDENTALCARE	JDEVOLUTION S	3.20	40.316.008.02-2	1.60	20 N·cm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
JDENTALCARE	CONICAL ABUTMENT	UNIVERSAL	40.314.008.01-2	1.40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
JDENTALCARE	JD ICON PLUS	3.70	40.318.003.02-2	1.80	25 N·cm	8,00	3,60	6,00	-	2,00	2,10	straight	-	chamfer 45°	Hex. 1.20
JDENTALCARE	JD ICON PLUS	4.00	40.318.003.02-2	1.80	25 N·cm	8,00	3,60	6,00	-	2,00	2,10	straight	-	chamfer 45°	Hex. 1.20
JDENTALCARE	JD ICON PLUS	4.80	40.318.003.02-2	1.80	25 N·cm	8,00	3,60	6,00	-	2,00	2,10	straight	-	chamfer 45°	Hex. 1.20
JDENTALCARE	JD ICON PLUS T	3.50	40.318.003.02-2	1.80	25 N·cm	8,00	3,60	6,00	-	2,00	2,10	straight	-	chamfer 45°	Hex. 1.20
JDENTALCARE	JD ICON PLUS T	3.50	40.318.003.02-2	1.80	25 N·cm	8,00	3,60	6,00	-	2,00	2,10	straight	-	chamfer 45°	Hex. 1.20
JDENTALCARE	JD ICON PLUS T	3.50	40.318.003.02-2	1.80	25 N·cm	8,00	3,60	6,00	-	2,00	2,10	straight	-	chamfer 45°	Hex. 1.20
KENTRON	INTERNAL HEXAGON	NP 3.8	40.316.003.01-2	1.60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
KENTRON	INTERNAL HEXAGON	RP 4.5	40.316.003.01-2	1.60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
KEYSTONE	INTERNAL TILOBE PRIMACONNEX	3.50	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
KEYSTONE	INTERNAL TILOBE PRIMACONNEX	4.10	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
KEYSTONE	INTERNAL TILOBE PRIMACONNEX	5.00	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
KEYSTONE	RESTORE	RD 4.1	40.320.003.03-2	2.00	30 N·cm	7,90	3,70	5,55	6,05	1,85	2,40	conical	45°	chamfer 45°	Hex. 1.20
KEYSTONE	RESTORE	SD 3.3	40.318.006.01-2	1.80	25 N·cm	6,00	3,00	3,67	4,00	2,00	2,40	conical	45°	chamfer 45°	Hex. 1.70
KLOCKNER	ESSENTIAL CONE	4.50	40.318.012.01-2	1.80	25 N·cm	7,25	2,40	4,93	5,25	2,00	2,15	conical	45°	chamfer 45°	Star 1.50
KLOCKNER	ESSENTIAL CONE	6.00	40.318.012.01-2	1.80	25 N·cm	7,25	2,40	4,93	5,25	2,00	2,15	conical	45°	chamfer 45°	Star 1.50
KLOCKNER	VEGA	NV	40.316.012.01-2	1.60	20 N·cm	8,00	2,65	5,53	6,00	2,00	2,15	conical	45°	chamfer 45°	Star 1.50
KLOCKNER	NK2-SK2	4.20	40.318.006.01-2	1.80	25 N·cm	6,00	3,00	3,67	4,00	2,00	2,40	conical	45°	chamfer 45°	Hex. 1.70
KLOCKNER	KL	NARROW	40.316.008.02-2	1.60	20 N·cm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
KLOCKNER	KL	REGULAR	40.320.003.02-2	2.00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
KLOCKNER	VEGA	RV	40.318.012.02-2	1.80	25 N·cm	8,00	2,60	5,68	6,00	2,00	2,15	conical	45°	chamfer 45°	Star 1.50
KLOCKNER	MULTI UNIT PERMANENT	UNIVERSAL	40.314.012.01-2	1.40	15 N·cm	4,50	1,70	2,01	2,40	2,10	2,15	conical	45°	chamfer 45°	Star 1.50
KLOCKNER	ESSENTIAL CONE OCTACONE 12°	4.50	40.318.006.01-2	1.80	25 N·cm	6,00	3,00	3,67	4,00	2,00	2,40	conical	45°	chamfer 45°	Hex. 1.70
KLOCKNER	ESSENTIAL CONE OCTACONE 25°	4.50	40.318.006.01-2	1.80	25 N·cm	6,00	3,00	3,67	4,00	2,00	2,40	conical	45°	chamfer 45°	Hex. 1.70
LASAK	BIONIQ	2.90	40.314.004.02-2	1.40	15 N·cm	8,40	2,50	5,99	6,70	1,70	2,00	conical	35°	chamfer 45°	Hex. 1.25
LASAK	BIONIQ	3.5/4/5	40.316.004.01-2	1.60	20 N·cm	8,60	2,70	6,16	6,90	1,70	2,00	conical	30°	chamfer 45°	Hex. 1.25
LASAK	MULTI UNIT	UNIVERSAL	40.314.004.03-2	1.40	15 N·cm	4,30	1,80	2,30	-	2,00	2,00	straight	-	chamfer 45°	Hex. 1.25
LEADER	TIXOS INTERNAL HEXAGON	3.50	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
LEADER	TIXOS INTERNAL HEXAGON	4.00	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
LEADER	TIXOS INTERNAL HEXAGON	5.00	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
LEADER	TIXOS INTERNAL HEXAGON	6.00	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
LEADER	TIXOS EXTERNAL HEXAGON	4.10	40.320.003.02-2	2.00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
LEADER	TIXOS EXTERNAL HEXAGON	5.00	40.320.003.02-2	2.00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
LEADER	TIXOS MORSE CONEXION	3.60	40.316.003.01-2	1.60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
LIFECORE	EXTERNAL HEXAGON	NARROW	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
LIFECORE	EXTERNAL HEXAGON	REGULAR	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
MEDENTIKA	MULTI-UNIT	UNIVERSAL	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
MEDENTIS	ICX TEMPLANT	3.45	40.316.007.01-2	1.60	20 N·cm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
MEDENTIS	ICX TEMPLANT	3.75	40.316.007.01-2	1.60	20 N·cm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
MEDENTIS	ICX TEMPLANT	4.10	40.316.007.01-2	1.60	20 N·cm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
MEDENTIS	ICX TEMPLANT	4.80	40.316.007.01-2	1.60	20 N·cm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
MEGAGEN	ANYONE EXTERNAL	SMALL 3.5	40.316.008.01-2	1.60	20 N·cm	7.00	2.70	5.15	-	1.85	2.30	straight	-	chamfer 45°	Unigrip
MEGAGEN	ANYONE EXTERNAL	WIDE 5	40.320.003.03-2	2.00	30 N·cm	7.90	3.70	5.55	6.05	1.85	2.40	conical	45°	chamfer 45°	Hex. 1.20
MEGAGEN	ANYONE EXTERNAL	SUPERWIDE 5.5	40.320.003.03-2	2.00	30 N·cm	7.90	3.70	5.55	6.05	1.85	2.40	conical	45°	chamfer 45°	Hex. 1.20
MEGAGEN	ANYRIDGE	SMALL	40.318.003.02-2	1.80	25 N·cm	8.00	3.60	6.00	-	2.00	2.10	straight	-	chamfer 45°	Hex. 1.20
MEGAGEN	ANYRIDGE	REGULAR	40.318.003.02-2	1.80	25 N·cm	8.00	3.60	6.00	-	2.00	2.10	straight	-	chamfer 45°	Hex. 1.20
MEGAGEN	ANYRIDGE	WIDE	40.318.003.02-2	1.80	25 N·cm	8.00	3.60	6.00	-	2.00	2.10	straight	-	chamfer 45°	Hex. 1.20
MEGAGEN	MINI	3/3.25	40.314.003.04-2	1.40	15 N·cm	7.50	2.50	5.45	5.70	1.80	1.85	conical	45°	chamfer 45°	Hex. 1.20
MEGAGEN	EZ PLUS	SMALL	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 1.20
MEGAGEN	MULTI UNIT N TYPE	UNIVERSAL	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
MEGAGEN	MULTI UNIT S TYPE	UNIVERSAL	40.314.003.01-2	1.40	15 N·cm	3.90	1.91	2.10	-	1.80	2.40	straight	-	chamfer 45°	Hex. 1.20
MEGAGEN	CONE ABUTMENT	3.80	40.320.003.05-2	2.00	30 N·cm	4.80	3.30	3.65	3.90	0.90	2.45	conical	45°	chamfer 45°	Hex. 1.20
MICRODENT	GENIUS	3.50	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
MICRODENT	GENIUS	4.00	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
MICRODENT	GENIUS	4.50	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
MICRODENT	TRYLOGIC	3.80	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
MICRODENT	TRYLOGIC	4.30	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
MICRODENT	TRYLOGIC	5.10	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
MICRODENT	EKTOS	3.50	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
MIS	MULTI-UNIT	UNIVERSAL	40.316.005.06-2	1.60	20 N·cm	3.80	1.80	2.20	-	1.60	2.10	straight	-	chamfer 45°	Hex. 1.27
MIS	LANCE	WIDE	40.320.003.01-2	2.00	30 N·cm	4.70	3.00	3.33	-	1.37	2.35	straight	-	chamfer 45°	Hex. 1.20
MIS	SEVEN/M4	NARROW	40.316.005.05-2	1.60	20 N·cm	7.60	2.70	5.21	5.50	2.10	2.10	conical	60°	chamfer 45°	Hex. 1.27
MIS	SEVEN/M4	STANDARD	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
MIS	SEVEN/M4	WIDE	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
MIS	C1	NARROW	40.316.005.05-2	1.60	20 N·cm	7.60	2.70	5.21	5.50	2.10	2.10	conical	60°	chamfer 45°	Hex. 1.27
MIS	C1	STANDARD	40.317.005.01-2	N1-72	25 N·cm	7.60	2.15	5.17	5.40	2.20	2.20	conical	60°	chamfer 45°	Hex. 1.27
MIS	C1	WIDE	40.317.005.01-2	N1-72	25 N·cm	7.60	2.15	5.17	5.40	2.20	2.20	conical	60°	chamfer 45°	Hex. 1.27
MIS	V3	STANDARD	40.317.005.01-2	N1-72	25 N·cm	7.60	2.15	5.17	5.40	2.20	2.20	conical	60°	chamfer 45°	Hex. 1.27
MOZO GRAU	MG INHEX	2.3 MINI	40.314.005.02-2	1.40	15 N·cm	7.50	2.40	5.50	5.70	1.80	2.10	conical	60°	chamfer 45°	Hex. 1.27
MOZO GRAU	MG INHEX	2.8 STANDARD	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
MOZO GRAU	MG INHEX	3.8 MAXI	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
MPI	CM INTERNAL HEXAGON	REGULAR	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
MPI	CM INTERNAL HEXAGON	WIDE	40.320.005.03-2	2,00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
NEOBIOTECH	IS CONE INTERNAL	S-NARROW	40.316.003.02-2	1.60	20 N·cm	10,20	2,00	7,88	8,20	2,00	2,20	conical	45°	chamfer 45°	Hex. 1.20
NEOBIOTECH	IS CONE INTERNAL	NARROW	40.320.003.04-2	2,00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20
NEOBIOTECH	IS CONE INTERNAL	REGULAR 4	40.320.003.04-2	2,00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20
NEOBIOTECH	EB EXTERNAL SYSTEM	NARROW	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
NEOBIOTECH	MULTI-UNIT	4.80	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
NEODENT	SMART HE	3.30	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
NEODENT	SMART HE	4.10	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
NEODENT	SMART HE	4.30	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
NEODENT	CONO MORSE MINI-PILAR	REGULAR	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
NEOSS	PROACTIVE /STRAIGHT/TAPERED	PROACTIVE	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
NOBEL BIO CARE	ACTIVE	3.00	40.314.008.02-2	1,40	15 N·cm	6,70	2,50	4,87	5,30	1,40	1,80	conical	35°	chamfer 45°	Unigrip
NOBEL BIO CARE	ACTIVE	NARROW	40.316.008.02-2	1.60	20 N·cm	7,30	2,70	5,15	5,90	1,40	2,20	conical	35°	chamfer 45°	Unigrip
NOBEL BIO CARE	BRANEMARK	NARROW	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
NOBEL BIO CARE	MULTI-UNIT	REGULAR	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
NOBEL BIO CARE	MULTI-UNIT	WIDE	40.318.006.01-2	1.80	25 N·cm	6,00	3,00	3,67	4,00	2,00	2,40	conical	45°	chamfer 45°	Hex. 1.70
NOBEL BIO CARE	REPLACE	NARROW	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
NOBEL BIO CARE	STANDARD	REGULAR	40.314.003.01-2	1,40	15 N·cm	3,90	1,91	2,10	-	1,80	2,40	straight	-	chamfer 45°	Hex. 1.20
NORIS MEDICAL	TUFF	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
NORIS MEDICAL	TUFF TT	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
NORIS MEDICAL	ONIX	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
NORIS MEDICAL	CORTICAL	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
NORIS MEDICAL	PTERYCORE	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
NORIS MEDICAL	PTERYFIT	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
NORIS MEDICAL	MULTI UNIT	UNIVERSAL	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
NORMON	NORMOIMPLANT HI	3.50	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
NOVA IMPLANTS	PSI/PCI	3.75	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
NUCLEOSS	T5 SHORT IMPLANT	5.80	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
NUCLEOSS	BONE LEVEL T PURE	4.20	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
NUCLEOSS	BONE LEVEL T PURE	5.00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
NUCLEOSS	BONE LEVEL T PURE	6.00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
NUCLEOSS	TISSUE LEVEL T3	5.00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
NUCLEOSS	TISSUE LEVEL T3	6.00	40.320.007.02-2	2,00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
OSSEOLIFE	OEX	NP 3.3	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
OSSEOLIFE	OEX	RP 4.1	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
OSSEOLIFE	OEX	WP 5	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
OSSEOLIFE	OIN	RP 4	40.320.007.02-2	2.00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
OSSEOLIFE	OIN	WP 5	40.320.007.02-2	2.00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
OSSTEM	TS	MINI	40.316.003.02-2	1.60	20 N·cm	10,20	2,00	7,88	8,20	2,00	2,20	conical	45°	chamfer 45°	Hex. 1,20
OSSTEM	TS	REGULAR	40.320.003.04-2	2.00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1,20
OSSTEM	TS MULTI UNIT	REGULAR	40.320.003.01-2	2.00	30 N·cm	4,70	3,00	3,33	-	1,37	2,35	straight	-	chamfer 45°	Hex. 1,20
OSSTEM	SS	REGULAR 4,8	40.320.007.02-2	2.00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
OSSTEM	SS	WIDE 6	40.320.007.02-2	2.00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
OSSTEM	US	MINI 3,5	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
OSSTEM	US	WIDE PS 5	40.320.003.03-2	2.00	30 N·cm	7,90	3,70	5,55	6,05	1,85	2,40	conical	45°	chamfer 45°	Hex. 1,20
OSTEOPLUS	ARGENTUM	4,20	40.318.006.01-2	1.80	25 N·cm	6,00	3,00	3,67	4,00	2,00	2,40	conical	45°	chamfer 45°	Hex. 1,70
OSTEOPLUS	IRIDIUM	3,50	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1,20
OSTEOPLUS	IRIDIUM	4,20	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
OSTEOPLUS	SHI	3,50	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1,27
OSTEOPLUS	SHI	4,50	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1,27
OSTEOPLUS	PLATINIUM	4,20	40.318.002.02-2	1.80	25 N·cm	8,30	2,60	6,60	-	1,70	2,45	straight	-	chamfer 45°	Cuad. 1,30
OT MEDICAL	OT-F2	NARROW	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
OT MEDICAL	OT-F2	REGULAR	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
OT MEDICAL	OT-F2	REGULAR	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1,20
OT MEDICAL	OT-F2	WIDE	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1,20
OT MEDICAL	OT-F3	NARROW	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
OT MEDICAL	OT-F3	REGULAR	40.316.008.01-2	1.60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
OT MEDICAL	OT-F3	REGULAR	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1,20
OT MEDICAL	OT-F3	WIDE	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1,20
OXY	EXTERNAL HEXAGON	MINI	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1,20
OXY	CONICAL CONNECTION	REGULAR	40.318.003.02-2	1.80	25 N·cm	8,00	3,60	6,00	-	2,00	2,10	straight	-	chamfer 45°	Hex. 1,20
PALTOP	ADVANCED CLASSIC	NARROW	40.316.004.04-2	1.60	20 N·cm	8,70	2,15	6,46	7,00	1,70	2,20	conical	45°	chamfer 45°	Hex. 1,25
PALTOP	ADVANCED CLASSIC	STANDARD	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1,27
PALTOP	ADVANCED +	NARROW	40.316.004.04-2	1.60	20 N·cm	8,70	2,15	6,46	7,00	1,70	2,20	conical	45°	chamfer 45°	Hex. 1,25
PALTOP	ADVANCED +	STANDARD	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1,27
PALTOP	DYNAMIC	NARROW	40.316.004.04-2	1.60	20 N·cm	8,70	2,15	6,46	7,00	1,70	2,20	conical	45°	chamfer 45°	Hex. 1,25
PALTOP	DYNAMIC	STANDARD	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1,27
PALTOP	DIVA	STANDARD	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1,27
PALTOP	CONICAL ACTIVE	NARROW	40.316.005.01-2	1.60	20 N·cm	7,50	3,60	5,44	5,90	1,60	2,13	conical	31°	chamfer 45°	Hex. 1,27
PALTOP	CONICAL ACTIVE	STANDARD	40.316.005.01-2	1.60	20 N·cm	7,50	3,60	5,44	5,90	1,60	2,13	conical	31°	chamfer 45°	Hex. 1,27
PHIBO	TSA S3	3,70	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1,27
PHIBO	TSA S4	4,70	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
PHIBO	TSA S5	6,00	40.318.002.01-2	1.80	25 N·cm	7,00	3,20	5,20	-	1,80	2,50	straight	-	chamfer 45°	Cuad. 1,30

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
PHIBO	TSH S2	3,30	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
PHIBO	TSH S3	4.00	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
PHIBO	TSH S4	4.00	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
PHIBO	TSH S5	5.00	40.320.003.01-2	2.00	30 N·cm	4.70	3.00	3.33	-	1.37	2.35	straight	-	chamfer 45°	Hex. 1.20
PROCLINIC	CILÍNDRICO EXTERNO	3,3 MINI	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
PROCLINIC	CILÍNDRICO EXTERNO	4,1 STANDARD	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
PROCLINIC	CILÍNDRICO EXTERNO	5 MAXI	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
PROCLINIC	CILÍNDRICO INTERNO	3,50	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
PROCLINIC	Aqua CM	3,50	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
RADHEX	PHE	3,50	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
RADHEX	PHI	3,50	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
RADHEX	PHI	4,50	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
ROSTER DENT	HEX	RD	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
ROSTER DENT	PLUS	4,80	40.320.007.01-2	2.00	30 N·cm	6.70	2.25	3.59	5.70	1.00	2.58	conical	15°	chamfer 45°	Torx T6
ROSTER DENT	CONICAL GOLD	3,50	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
ROSTER DENT	CONICAL GOLD	4,30	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
ROSTER DENT	CONICAL GOLD P4	3,75	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
ROSTER DENT	CONICAL GOLD P4	4,30	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
SAMO	BIOSPARK / BIOROUGH TAPERED / CYLINDRICAL	4,50	40.316.005.04-2	1.60	20 N·cm	10.50	3.20	8.15	8.40	2.10	2.10	conical	45°	chamfer 45°	Hex. 1.27
SEWON MEDIX	IH2 SLA SYSTEM / IH2 RBM SYSTEM	MINI	40.316.003.02-2	1.60	20 N·cm	10.20	2.00	7.88	8.20	2.00	2.20	conical	45°	chamfer 45°	Hex. 1.20
SEWON MEDIX	IH2 SLA SYSTEM / IH2 RBM SYSTEM	REGULAR	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
SEWON MEDIX	IH SYSTEM	UNIVERSAL	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
SGS	P1/P7	3,75	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
SGS	P7N	3.00	40.314.005.01-2	1.40	15 N·cm	7.60	3.60	5.21	6.00	1.60	2.15	conical	25°	chamfer 45°	Hex. 1.27
SGS	S6/S7	3,75	40.314.003.02-2	1.40	15 N·cm	4.00	2.00	2.20	-	1.80	2.30	straight	-	chamfer 45°	Hex. 1.20
SIGNO VINCES	DUOCON	CM 3,8	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
SIGNO VINCES	DUOCON	CM 4,6	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
SIGNO VINCES	COMPACT	CM 3,8	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
SIGNO VINCES	DUO	3,30	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
SIGNO VINCES	DUO	4,10	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
SIGNO VINCES	INTTEGRA	3,30	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
SIGNO VINCES	INTTEGRA	4,10	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
SIGNO VINCES	Mini Pilar Cónico		40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
SIGNO VINCES	INFRA	CM	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
SIN IMPLANTS	SW MORSE	3,50	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	SW MORSE	3,80	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 1.20

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
SIN IMPLANTS	SW MORSE	4,50	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	SW MORSE	5,50	40.316.003.01-2	1,60	20 N·cm	8,40	2,50	6,60	-	1,80	2,00	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	HI SW	3,80	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
SIN IMPLANTS	HI SW	4,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
SIN IMPLANTS	HI SW	5,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
SIN IMPLANTS	TRYON	4,10	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	TRYON	5,00	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	REVOLUTION	3,40	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	REVOLUTION	4,10	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	REVOLUTION	5,00	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	STYLUS	3,40	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	STYLUS	4,10	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	STYLUS	5,00	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	REVOLUTION COMPACT	5,00	40.320.003.01-2	2,00	30 N·cm	4,70	3,00	3,33	-	1,37	2,35	straight	-	chamfer 45°	Hex. 1.20
SIN IMPLANTS	Pilar Cónico	Universal	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
SIN IMPLANTS	Mini Pilar	Universal	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
SIN IMPLANTS	Micro Multi Unit	Universal	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
SIN IMPLANTS	Pilar Multifuncional	Universal	40.320.003.01-2	2,00	30 N·cm	4,70	3,00	3,33	-	1,37	2,35	straight	-	chamfer 45°	Hex. 1.20
SOUTHERN IMPLANT	EXTERNAL HEXAGON	3,40	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SOUTHERN IMPLANT	EXTERNAL HEXAGON	5,00	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SOUTHERN IMPLANT	EXTERNAL HEXAGON	6,00	40.320.003.02-2	2,00	30 N·cm	7,00	3,25	5,00	-	2,00	2,40	straight	-	chamfer 45°	Hex. 1.20
SOUTHERN IMPLANT	MSC	3,00	40.316.008.01-2	1,60	20 N·cm	7,00	2,70	5,15	-	1,85	2,30	straight	-	chamfer 45°	Unigrip
SOUTHERN IMPLANT	IT CONNECTION	4,80	40.320.007.01-2	2,00	30 N·cm	6,70	2,25	3,59	5,70	1,00	2,58	conical	15°	chamfer 45°	Torx T6
SOUTHERN IMPLANT	IT CONNECTION	6,50	40.320.007.01-2	2,00	30 N·cm	6,70	2,25	3,59	5,70	1,00	2,58	conical	15°	chamfer 45°	Torx T6
SOUTHERN IMPLANT	TRI-NEX	3,50	40.318.008.01-2	1,80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
SOUTHERN IMPLANT	INTERNAL HEX	UNIVERSAL	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
SOUTHERN IMPLANT	COMPACT CONICAL	4,80	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
SOUTHERN IMPLANT	COMPACT CONICAL	6,00	40.314.008.01-2	1,40	15 N·cm	3,60	1,80	2,10	-	1,50	2,10	straight	-	chamfer 45°	Unigrip
SOUTHERN IMPLANT	DEEP CONICAL	2,45	40.314.005.02-2	1,40	15 N·cm	7,50	2,40	5,50	5,70	1,80	2,10	conical	60°	chamfer 45°	Hex. 1.27
SOUTHERN IMPLANT	DEEP CONICAL	2,95/3,1	40.316.005.02-2	1,60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
SOUTHERN IMPLANT	DEEP CONICAL	4.10	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
STERI-OSS	CONEXIÓN EXTERNA	3.30	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
STERNGOLD	STERN EX	4.10	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
STERNGOLD	STERN IC	4.80	40.320.007.01-2	2.00	30 N·cm	6.70	2.25	3.59	5.70	1.00	2.58	conical	15°	chamfer 45°	Torx T6
STRAUMANN	STANDARD PLUS NARROW NECK 3.3	3.50	40.318.003.01-2	1.80	25 N·cm	6.80	3.30	5.20	-	1.60	2.30	straight	-	chamfer 45°	Hex. 1.20
STRAUMANN	BONE LEVEL	2.90	40.314.007.02-2	1.40	15 N·cm	7.00	2.10	6.20	2.25	0.80	2.10	conical	15°	chamfer 45°	Torx T6
STRAUMANN	BONE LEVEL	3.30	40.316.007.01-2	1.60	20 N·cm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
STRAUMANN	BONE LEVEL	4.10	40.316.007.01-2	1.60	20 N·cm	7.90	2.00	5.63	6.90	1.00	2.18	conical	15°	chamfer 45°	Torx T6
STRAUMANN	BONE LEVEL MULTI-BASE		40.314.003.02-2	1.40	15 N·cm	4.00	2.00	2.20	-	1.80	2.30	straight	-	chamfer 45°	Hex. 1.20
STRAUMANN	INTERNAL OCTAGON	REGULAR 4.8	40.320.007.01-2	2.00	30 N·cm	6.70	2.25	3.59	5.70	1.00	2.58	conical	15°	chamfer 45°	Torx T6
STRAUMANN	INTERNAL OCTAGON	WIDE 6.5	40.320.007.01-2	2.00	30 N·cm	6.70	2.25	3.59	5.70	1.00	2.58	conical	15°	chamfer 45°	Torx T6
STRAUMANN	EXTERNAL OCTAGON	REGULAR 4.8	40.320.002.01-2	2.00	30 N·cm	4.90	3.26	3.26	3.50	1.40	2.49	conical	45°	chamfer 45°	Cuad. 1.30
STRAUMANN	EXTERNAL OCTAGON	WIDE 6.5	40.320.002.01-2	2.00	30 N·cm	4.90	3.26	3.26	3.50	1.40	2.49	conical	45°	chamfer 45°	Cuad. 1.30
STRAUMANN	SCREW-RETAINED NC/RC		40.314.007.01-2	1.40	15 N·cm	4.00	1.80	2.01	2.80	1.20	2.20	conical	35°	chamfer 45°	Torx T6
SWEDEN & MARTINA	OUTLINK2	SLIM 3	40.318.006.01-2	1.80	25 N·cm	6.00	3.00	3.67	4.00	2.00	2.40	conical	45°	chamfer 45°	Hex. 1.70
SWEDEN & MARTINA	OUTLINK2	3.30	40.318.006.01-2	1.80	25 N·cm	6.00	3.00	3.67	4.00	2.00	2.40	conical	45°	chamfer 45°	Hex. 1.70
SWEDEN & MARTINA	OUTLINK2	4.1 SP	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	PREMIUM KOHNO	3.30	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	PREMIUM KOHNO	3.80	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	PREMIUM KOHNO	4.25	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
SWEDEN & MARTINA	PREMIUM KOHNO	5.00	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
SWEDEN & MARTINA	PREMIUM KOHNO	6.00	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
SWEDEN & MARTINA	STARK	REGULAR	40.320.007.01-2	2.00	30 N·cm	6.70	2.25	3.59	5.70	1.00	2.58	conical	15°	chamfer 45°	Torx T6
SWEDEN & MARTINA	GLOBAL	3.80	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	GLOBAL	4.30	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	GLOBAL	4.80	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	GLOBAL	5.50	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	SHELTA	3.80	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
SWEDEN & MARTINA	SHELTA	4.30	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	SHELTA	5.00	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SWEDEN & MARTINA	PRAMA	3.40	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SYSTHEX	ATTRACT	PLATAFORMA SWITCH	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 1.20
SYSTHEX	CLASSIC-IN CR/ESTETIC-IN CR	3.60	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 1.20
SYSTHEX	CLASSIC-IN CR/ESTETIC-IN CR	4.30	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SYSTHEX	CLASSIC-IN CR/ESTETIC-IN CR	5.00	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
SYSTHEX	CLASSIC-CI/ESTETIC-CI	4.10	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
SYSTHEX	FIT-HEX	3.30	40.318.006.01-2	1.80	25 N·cm	6.00	3.00	3.67	4.00	2.00	2.40	conical	45°	chamfer 45°	Hex. 1.70
SYSTHEX	CLASSIC TRADICIONAL	4.10	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
TAV DENTAL	SIRIUS	NP	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
TAV DENTAL	SIRIUS	SP	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
TAV DENTAL	SILHOUETTE	NP	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
TAV DENTAL	SILHOUETTE	SP	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
TBR	Z1/M	5.00	40.316.003.01-2	1.60	20 N·cm	8.40	2.50	6.60	-	1.80	2.00	straight	-	chamfer 45°	Hex. 1.20
TBR	CONNECT/INFINITI	4.40	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
TBR	CONNECT/INFINITI	5.00	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
THOMMEN	HEXAGONO INTERNO	3.50	40.314.003.02-2	1.40	15 N·cm	4.00	2.00	2.20	-	1.80	2.30	straight	-	chamfer 45°	Hex. 1.20
THOMMEN	HEXAGONO INTERNO	4.00	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
THOMMEN	HEXAGONO INTERNO	4.50	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
THOMMEN	HEXAGONO INTERNO	5.00	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
THOMMEN	HEXAGONO INTERNO	6.00	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
TITANIUM FIX	b-FIX	REGULAR	40.316.005.02-2	1.60	20 N·cm	8.25	3.00	6.25	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
TITANIUM FIX	b-FIX	LARGA	40.320.005.03-2	2.00	30 N·cm	10.30	3.50	8.30	-	2.00	2.33	straight	-	chamfer 45°	Hex. 1.27
TITANIUM FIX	i-FIX	REGULAR	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
TITANIUM FIX	i-FIX	LARGA	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
TITANIUM FIX	e-FIX EXT HEXAGON TORQUE INTERNO	REGULAR 4.1	40.320.003.02-2	2.00	30 N·cm	7.00	3.25	5.00	-	2.00	2.40	straight	-	chamfer 45°	Hex. 1.20
TITANIUM FIX	e-FIX SILVER	REGULAR 4.1	40.320.003.01-2	2.00	30 N·cm	4.70	3.00	3.33	-	1.37	2.35	straight	-	chamfer 45°	Hex. 1.20
TITANIUM FIX	c-FIX	REGULAR	40.320.007.02-2	2.00	30 N·cm	7.40	3.30	6.00	-	1.40	2.30	straight	-	semi-sphere	Torx T6
TREE-OSS	SIMPLE	3.75	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
TREE-OSS	HEXAGONO EXTERNO	3.50	40.316.008.01-2	1.60	20 N·cm	7.00	2.70	5.15	-	1.85	2.30	straight	-	chamfer 45°	Unigrip
TREE-OSS	TRI-CHANNEL INTERNO	3.50	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip
TREE-OSS	MULTI UNIT	UNIVERSAL	40.314.008.01-2	1.40	15 N·cm	3.60	1.80	2.10	-	1.50	2.10	straight	-	chamfer 45°	Unigrip
TRI IMPLANTS	TRI-VENT	3.50	40.317.004.01-2	N1-72	25 N·cm	7.60	2.80	5.60	5.77	1.83	2.30	conical	70°	chamfer 45°	Hex. 1.27
TRI IMPLANTS	TRI-NARROW	3.20	40.318.008.01-2	1.80	25 N·cm	8.30	2.50	6.50	-	1.80	2.45	straight	-	chamfer 45°	Unigrip

STRAIGHT SCREWS TECHNICAL SPECIFICATIONS

BRAND	MODEL	PLATFORM	REFERENCE	METRIC	TORQUE	TOTAL LENGTH	THREAD LENGTH	A LENGTH	B LENGTH	C LENGTH	HEAD DIAMETER	SEAT	ANGLE	THREAD ENTRY	CONNECTION
TRI IMPLANTS	MULTI UNIT	UNIVERSAL	40.318.005.01-2	1.80	25 N·cm	4,50	2,30	2,80	2,90	1,60	2,33	conical	70°	chamfer 45°	Hex. 1.27
UFIT	GT2	Mini	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
UFIT	GT2	REGULAR	40.320.005.03-2	2.00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
UFIT	GT2	WIDE	40.320.005.03-2	2.00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
UFIT	GT2	ULTRA WIDE	40.320.005.03-2	2.00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
UFIT	GT2	Mini	40.316.005.02-2	1.60	20 N·cm	8,25	3,00	6,25	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
UFIT	NT2	REGULAR	40.320.005.03-2	2.00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
UFIT	NT2	WIDE	40.320.005.03-2	2.00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
UFIT	NT2	ULTRA WIDE	40.320.005.03-2	2.00	30 N·cm	10,30	3,50	8,30	-	2,00	2,33	straight	-	chamfer 45°	Hex. 1.27
VULKAN IMPLANTS	IN-HEX	3.75	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
WINSIX	TTi	PLATFORM SWITCHING	40.318.008.01-2	1.80	25 N·cm	8,30	2,50	6,50	-	1,80	2,45	straight	-	chamfer 45°	Unigrip
XIVE	XIVE	3.00	40.314.003.03-2	1.40	15 N·cm	7,60	2,40	6,05	6,30	1,30	1,90	conical	45°	chamfer 45°	Hex. 1.20
XIVE	XIVE	3.40	40.316.004.02-2	1.60	20 N·cm	8,80	3,00	6,73	7,00	1,80	2,10	conical	45°	chamfer 45°	Hex. 1.25
XIVE	XIVE	3.80	40.316.004.02-2	1.60	20 N·cm	8,80	3,00	6,73	7,00	1,80	2,10	conical	45°	chamfer 45°	Hex. 1.25
XIVE	XIVE	4.50	40.316.004.02-2	1.60	20 N·cm	8,80	3,00	6,73	7,00	1,80	2,10	conical	45°	chamfer 45°	Hex. 1.25
XIVE	XIVE	5.50	40.316.004.02-2	1.60	20 N·cm	8,80	3,00	6,73	7,00	1,80	2,10	conical	45°	chamfer 45°	Hex. 1.25
YES IMPLANTS	S-SYSTEM	NARROW	40.320.003.04-2	2.00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20
YES IMPLANTS	S-SYSTEM	REGULAR	40.320.003.04-2	2.00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20
YES IMPLANTS	S-SYSTEM	WIDE	40.320.003.04-2	2.00	30 N·cm	8,40	2,75	5,68	6,35	2,05	2,30	conical	45°	chamfer 45°	Hex. 1.20
ZIACOM (OSSEOLIFE)	OEX - EXTERNAL HEXAGON	NP 3.3	40.318.003.01-2	1.80	25 N·cm	6,80	3,30	5,20	-	1,60	2,30	straight	-	chamfer 45°	Hex. 1.20
ZIACOM (OSSEOLIFE)	ION - INTERNAL HEXAGON	RP 4	40.320.007.02-2	2.00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
ZIACOM (OSSEOLIFE)	ION - INTERNAL HEXAGON	WP 5	40.320.007.02-2	2.00	30 N·cm	7,40	3,30	6,00	-	1,40	2,30	straight	-	semi-sphere	Torx T6
ZIMMER	EZTETIC	NP 2.9	40.316.004.04-2	1.60	20 N·cm	8,70	2,15	6,46	7,00	1,70	2,20	conical	45°	chamfer 45°	Hex. 1.25
ZIMMER	SCREW-VENT	3.50	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ZIMMER	SCREW-VENT	4.50	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ZIMMER	SCREW-VENT	5.70	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ZIMMER	SWISS PLUS	3.80	40.317.004.01-2	N1-72	25 N·cm	7,60	2,80	5,60	5,77	1,83	2,30	conical	70°	chamfer 45°	Hex. 1.27
ZIMMER	TAPERED ABUTMENT MULTI-UNIT	UNIVERSAL	40.317.005.03-2	N1-72	25 N·cm	8,50	1,60	6,58	7,00	1,50	2,20	conical	45°	chamfer 45°	Hex. 1.27

SCREWDRIVERS




STRAIGHT SCREWS

Straight screws cover all the thread metrics available on the market. We have several lengths for each metric to make the adaptation to the milled structures easier.



DYNAMIC MILLING TOOL SPECIFICATIONS

MAIN COMPATIBILITY	REFERENCE	CUTTING DIAMETER	SEAT	CUTTING LENGTH	USEFUL LENGTH (max. drilling depth)	STEM CUTTING DIAMETER	SUPPORT DIAMETER (SHANK)	TOTAL LENGTH
		ϕ_c		L_c	L_u	ϕ_v	ϕ_s	L_t
BEGO RS/RSX 3° ASTRA EVOLUTION 3.0° *Only for titanium and soft materials	33.325.472.01-2	1.4	25	0.4	4.7	0.5	3	50
	33.425.472.01-2	1.4	25	0.4	4.7	0.5	4	50
	33.625.472.01-2	1.4	25	0.4	4.7	0.5	6	50
STRAUMANN BONE LEVEL NP STRAUMANN BONE LEVEL RP MEDENTIS ICX TEMPLANT 4.1 STRAUMANN SYNOCTA 3.5 MEDENTIS ICX NARROW	33.315.804.01-2	1.6	15	0.7	8	0.65	3	50
	33.415.804.01-2	1.6	15	0.7	8	0.65	4	50
	33.615.804.01-2	1.6	15	0.7	8	0.65	6	50
ANTHOGYR AXIOM RG/PX XNP ANTHOGYR AXIOM RG/PX RP ANTHOGYR AXIOM RG/PX WP	33.320.704.01-2	1.6	20	0.7	7	0.8	3	50
	33.420.704.01-2	1.6	20	0.7	7	0.8	4	50
	33.620.704.01-2	1.6	20	0.7	7	0.8	6	50
ASTRA EVOLUTION 3.6 ANKYLOS ALPHABIO CONICAL STANDARD CONNECTION LASAK BIONIQ OR NEODENT GM ANKYLOS BALANCE BASE	33.330.734.01-2	1.6	30	0.7	7.3	0.8	3	50
	33.430.734.01-2	1.6	30	0.7	7.3	0.8	4	50
	33.630.734.01-2	1.6	30	0.7	7.3	0.8	6	50
NOBEL BIO CARE ACTIVE NP NOBEL BIO CARE ACTIVE 3.0 LASAK BIONIQ QN	33.335.754.01-2	1.6	35	0.7	7.5	0.65	0.5	50
	33.435.754.01-2	1.6	35	0.7	7.5	0.65	4	50
	33.635.754.01-2	1.6	35	0.7	7.5	0.65	6	50
OSSTEM TS NP CAMLOG SCREW LINE 3.8 NP CAMLOG SCREW LINE 4.3 RP KLOCKNER VEGA NV XIVE S 3.4 BIOTECH DENTAL KONTACT XNP BIOTECH DENTAL KONTACT RP DIO UF NP CAMLOG SCREW-LINE 3.3	33.345.804.01-2	1.6	45	0.7	8	0.65	3	50
	33.445.804.01-2	1.6	45	0.7	8	0.65	4	50
	33.645.804.01-2	1.6	45	0.7	8	0.65	6	50
MIS C1 NP MIS M4 NP CONOLOG 3.8 CONOLOG 4.3 ASTRA YELLOW ALPHABIO CONICAL HEX CONNECTION	33.360.754.01-2	1.6	60	0.7	7.5	0.65	3	50
	33.460.754.01-2	1.6	60	0.7	7.5	0.65	4	50
	33.660.754.01-2	1.6	60	0.7	7.5	0.65	6	50
BIOMET 3i CERTAIN NP ASTRA AQUA	33.390.754.01-2	1.6	90	0.7	7.5	0.65	3	50
	33.490.754.01-2	1.6	90	0.7	7.5	0.65	4	50
	33.690.754.01-2	1.6	90	0.7	7.5	0.65	6	50
ASTRA EVOLUTION 4.2	33.350.775.01-2	1.7	50	0.7	7.7	0.8	3	50
	33.450.775.01-2	1.7	50	0.7	7.7	0.8	4	50
	33.650.775.01-2	1.7	50	0.7	7.7	0.8	6	50



DYNAMIC MILLING TOOL SPECIFICATIONS

MAIN COMPATIBILITY	REFERENCE	CUTTING DIAMETER	SEAT	CUTTING LENGTH	USEFUL LENGTH (max. drilling depth)	STEM CUTTING DIAMETER	SUPPORT DIAMETER (SHANK)	TOTAL LENGTH
		Øc		Lc	Lu	Øv	Øs	Lt
BIOMET 3i CERTAIN RP NOBEL BIOACARE BRANEMARK NP NOBEL BIOACARE REPLACE NP MEGAGEN ANYRIDGE RP BIOMET 3i CERTAIN WP	33.390.805.01-2	1,7	90	0,7	8	0,65	3	50
	33.490.805.01-2	1,7	90	0,7	8	0,65	4	50
	33.690.805.01-2	1,7	90	0,7	8	0,65	6	50
"BEGO S/RI 3.25-3.75 BEGO S/RI 4.1 BEGO S/RI 4.5 BEGO S/RI 5.50 STRAUMANN SCREW-RETAINED NC/RC BEGO MULTI-PLUS"	33.335.676.01-2	1,8	35	1	6,7	0,9	3	50
	33.435.676.01-2	1,8	35	1	6,7	0,9	4	50
	33.635.676.01-2	1,8	35	1	6,7	0,9	6	50
"KLOCKNER ESSENTIAL CONE 4.5 DIRECTO IMPLANTE KLOCKNER ESSENTIAL CONE 4.5 OCTACONE 12° KLOCKNER ESSENTIAL CONE 4.5 OCTACONE 25° KLOCKNER VEGA RV XIVE S 3.8 XIVE S 4.5 BIOHORIZONS 3.0 STRAUMANN SYNOCCTA 6.5 STRAUMANN BLX RB STRAUMANN BLX WB STRAUMANN TLX NT STRAUMANN TLX RT STRAUMANN TLX WT DENTIMU MU SUPERLINE"	33.345.856.01-2	1,8	45	1	8,5	0,9	3	50
	33.445.856.01-2	1,8	45	1	8,5	0,9	4	50
	33.645.856.01-2	1,8	45	1	8,5	0,9	6	50
"MIS C1 RP PALTOP UNIVERSAL MULTI UNIT MIS C1 WP S&M PREMIUM KHONO 3.3 S&M PREMIUM KHONO 3.8 S&M OUTLINK 3.3 S&M OUTLINK 4.1 S&M PREMIUM KHONO 4.25 BREDENT SKY NP BREDENT SKY RP ADIN TOUAREG/CLOSEFIT NP ADIN TOUAREG/CLOSEFIT UNP CAMLOG CONELOG 3.3 GLOBAL D (TEKKA) EASY IMPLANT MINI ROOT R"	33.360.756.01-2	1,8	60	1	7,5	0,9	3	50
	33.460.756.01-2	1,8	60	1	7,5	0,9	4	50
	33.660.756.01-2	1,8	60	1	7,5	0,9	6	50
"ZIMMER SCREW-VENT 3.5 ZIMMER SCREW-VENT 4.5 ASTRA EVOLUTION UNIT ABUTMENT ZIMMER TYPE 5.7 OXY FIXO	33.370.716.01-2	1,8	70	1	7,1	0,9	3	50
	33.470.716.01-2	1,8	70	1	7,1	0,9	4	50
	33.670.716.01-2	1,8	70	1	7,1	0,9	6	50

DYNAMIC MILLING TOOL SPECIFICATIONS

MAIN COMPATIBILITY	REFERENCE	CUTTING DIAMETER	SEAT	CUTTING LENGTH	USEFUL LENGTH (max. drilling depth)	STEM CUTTING DIAMETER	SUPPORT DIAMETER (SHANK)	TOTAL LENGTH
		Øc		Lc	Lu	Øv	Øs	Lt
NOBEL BIO CARE BRANEMARK RP NOBEL BIO CARE MULTI-UNIT RP BIOMET 3i OSSEOTITE NP BTI EXTERNAL CONNECTION NP BTI INTERNAL CONNECTION RP MIS MULTI-UNIT ST KEYSTONE PRIMA NP KEYSTONE PRIMA RP KEYSTONE PRIMA WP NEOSS PROACTIVE 3.4 NEOSS PROACTIVE 4.1 BIOMET 3i OSSEOTITE WP BTI EXTERNAL CONNECTION WP BTI MULTI-IM UNIVERSAL RP ANTHOGYRD MULTI-UNIT 4.8 BEGO MINI BTI INTERNAL WP LASAK MULTI-UNIT QN/OR SIC SICACE 3.3 SIC SICACE 4.2 IMPLANT DIRECT	33.390.716.01-2	1,8	90	1	7,1	0,9	3	50
	33.490.716.01-2	1,8	90	1	7,1	0,9	4	50
	33.690.716.01-2	1,8	90	1	7,1	0,9	6	50
*STRAUMANN INTERNAL OCTAGON RP STRAUMANN INTERNAL OCTAGON 6.5"	33.315.708.01-2	1,8	15	1	7	1	3	50
	33.415.708.01-2	2	15	1	7	1	4	50
	33.615.708.01-2	2	15	1	7	1	6	50
STRAUMANN SYNOCTA RP	33.330.708.01-2	2	30	1	7	1	3	50
	33.430.708.01-2	2	30	1	7	1	4	50
	33.630.708.01-2	2	30	1	7	1	6	50
NOBEL BIO CARE ACTIVE RP NOBEL BIO CARE ACTIVE WP	33.335.758.01-2	2	35	1	7,5	1	3	50
	33.435.758.01-2	2	35	1	7,5	1	4	50
	33.635.758.01-2	2	35	1	7,5	1	6	50
OSSTEM TS RP CAMLOG SCREW-LINE 5.0 CAMLOG SCREW-LINE 6.0	33.345.808.01-2	2	45	1	8	1	3	50
	33.445.808.01-2	2	45	1	8	1	4	50
	33.645.808.01-2	2	45	1	8	1	6	50
NOBEL BIO CARE REPLACE RP ASTRA LILAC NOBEL BIO CARE REPLACE WP ASTRA EVOLUTION 4.8 NOBEL BIO CARE BRANEMARK WP ASTRA EVOLUTION 5.4 NOBEL BIO CARE REPLACE 6.0	33.390.958.01-2	2	90	1	9,5	1	3	50
	33.490.958.01-2	2	90	1	9,5	1	4	50
	33.690.958.01-2	2	90	1	9,5	1	6	50

SCREWDRIVER ADAPTOR

Screwdriver for the Dynamic Scanbody System

Ref. 43.621.410.01-2

Screwdriver with manual handle
Standard length: 21 mm



Ref. 43.624.410.01-2

Contra-angle
Length: 24 mm



Ref. 43.621.415.01-2

Tiny Screwdriver with manual handle
Length: 21 mm



Ref. 43.620.411.01-2

Multi-Unit
Contra-angle
Length: 20 mm



COMPLEMENTS

Manual handle

Made of stainless steel.
They are used to connect screwdriver bits with the contra-angle connection



Large manual handle for laboratory

Ref. 49.601.000.03-2

Ideal to manipulate models in the laboratory
Length: 55.65 mm



Manual handle for clinic

Ref. 49.601.000.01-2

Clinic handle: used to position the prosthesis in the mouth prior to torque control in the clinic.
Length: 15.65 mm



Universal manual torque wrench prosthetic

Ref. 11.990.990.07-2

Torque wrench
4 mm square connection.
Torque 10-35N.c



Dynamic Screw Transfer

Ref. 49.413.000.01-2

Manual torque wrench adapter prosthetic

Piece to connect the screwdriver with contra-angle connection to the torque wrench.



Universal Manual torque wrench adapter
Ref. 49.604.000.05-2
4 mm Square connection



Straumann Manual torque wrench adapter
Ref. 49.604.000.07-2
Straumann connection



Nobel Biocare Manual torque wrench adapter
Ref. 49.604.000.08-2



MIS Manual torque wrench adapter
Ref. 49.604.000.09-2

DYNAMIC SCREWDRIVERS

Screwdriver with hexalobular head, exclusively to the 3.0 Dynamic Abutment System.
Lengths: 18, 24, 32 mm

Hexalobular 1,70 mm. Length: 18 mm

Ref. 43.618.201.01-2



Hexalobular 1,70 mm. Length: 24 mm

Ref. 43.624.201.01-2



Hexalobular 1,70 mm Length: 32 mm

Ref. 43.632.201.01-2



DAS MU SYSTEM COMPONENTS



Ratchet
49.409.000.01-2



Analog
22.612.209.01-2



MU Scanbody 8 mm
53.422.209.02-2 (Non-engaging*)



Screwdriver
43.321.316.01-2



Titanium Abutment
35.312.209.21-2 (Non-engaging*)
35.322.209.21-2 (Engaging*)



MU Dynamic Scanbody
52.408.137.01-2



43.322.316.01-2



Mu ANG Insertion tool
49.422.000.01-2



Digital Analog
34.312.209.01-2



Dynamic Scanbody Adaptor
50.312.209.01-2



Healing Cap Regular
40.320.003.88-2



Digital Analog Positioner
49.309.000.01-2



Screwdriver Adaptor
43.621.410.01-2
43.624.410.01-2



Healing Cap Wide
40.320.003.89-2



MU ScAnalog
23.412.209.01-2



Impression coping
29.301.000.11-2 (Non-engaging*)
29.301.000.10-2 (Engaging*)



MU Scanbody 4.5 mm
53.412.209.01-2



Reference Scanbody

54.322.209.31-2 (Non-engaging*)
54.312.209.31-2 (Engaging*)

Peek Pins



TYPE A
49.414.000.01-2 (6 mm)
49.415.000.01-2 (9 mm)
49.416.000.01-2 (13 mm)



TYPE B
49.414.000.02-2 (6 mm)
49.415.000.02-2 (9 mm)
49.416.000.02-2 (13 mm)



TYPE C
49.414.000.03-2 (6 mm)
49.415.000.03-2 (9 mm)
49.416.000.03-2 (13 mm)



TYPE D
49.416.000.04-2 (13 mm)



TYPE E
49.416.000.05-2 (13 mm)

CAPS

Regular



49.418.000.01-2 (3,8 mm)
49.419.000.01-2 (6 mm)
49.420.000.01-2 (8 mm)

Wide



49.418.000.02-2 (3,8 mm)
49.419.000.02-2 (6 mm)
49.420.000.02-2 (8 mm)



MU Dynamic TiBase

31.312.209.01-2 (Non-engaging*)
31.312.209.01-2 (Engaging*)



MU Dynamic 3TiBase

31.322.209.21-2 (Non-engaging*)



Dynamic Screw

41.320.040.01-2



Provisional (temporary) Dynamic Screw

41.320.050.02-2



Straight Screw

40.320.003.06-2



Dynamic Screwdriver

43.618.201.01-2 (18 mm)
43.624.201.01-2 (24 mm)
43.632.201.01-2 (32 mm)



Screwdriver Hex.1,2

43.601.103.02-2



MU DMTone

33.391.716.01-2 Shank Ø3
33.491.716.01-2 Shank Ø4
33.691.716.01-2 Shank Ø6

TALLADIUM GUARANTEE

TERMS AND CONDITIONS

These guarantee terms and conditions ("T&C") cover the entire range of Talladium products ("Products"), manufactured by TALLADIUM ESPAÑA S.L. and distributed by Geoda Medical S.L. or official dealers. The guarantee described in these T&C is exclusively in benefit of the clinician ("Clinician") and of the dental technician ("Technician") and not for the benefit of third parties or institutions, including patients.

GUARANTEE PERIOD

TALLADIUM ESPAÑA S.L. offers a lifelong guarantee for its entire range of products starting from the date of issue of the invoice.

GUARANTEE SCOPE

Subject to the limitations and exceptions described in these T&C, TALLADIUM ESPAÑA S.L. will offer the following benefits:

QUALITY: If there are defects in the materials or in the manufacturing of the Product, TALLADIUM ESPAÑA S.L. will replace the Product with no additional cost.

SAFETY: If, having complied with all the product indications, the prosthesis should have to be made again, due to a fault in the Dynamic Abutment or Dynamic Titanium Base system, TALLADIUM ESPAÑA S.L. will replace the abutments and screws necessary to remake the prosthesis, as well as the costs derived from its manufacturing.

In case of having used our products and having complied with all the product indications, the implants suffer any damage, TALLADIUM ESPAÑA S.L. will pay the cost of the implants. This coverage will only be valid during the first 6 months after the collocation of the prosthesis which includes our products.

CLAIM REQUIREMENTS AND PROCEDURE

To receive the benefits indicated in these T&C, the treating Clinician must satisfy the following requirements:

- a) The claim must be notified to TALLADIUM ESPAÑA S.L. within (30) days since the date the claimed defect was detected.
- b) This requires that the Clinician or Technician must contact the customer service department by telephone or by e-mail to make the claim.
- c) A claim form will be completed, which, together with a document or report which justifies the faulty Product and the faulty Product itself, will be sent by the customer to TALLADIUM ESPAÑA S.L. offices, within the previously indicated period.
- d) Clinicians or Technicians presenting a claim in agreement with these T&C must be up to date in any payments owing to TALLADIUM ESPAÑA S.L. or to any of its subsidiaries, at the time when the claim form is presented.
- e) All the use procedures of our Products must be carried out in agreement with the instructions of TALLADIUM ESPAÑA S.L. as well as in accordance with commonly accepted dentistry practices.
- f) The expenses derived from this procedure will be assumed by the customer. The return shipping costs will be assumed by TALLADIUM ESPAÑA S.L. in all those cases covered by these T&C. Regardless of the guarantee rights, claims should be notified as soon as possible in order to comply with regulatory requirements.

GENERAL LIMITATIONS OF THIS GUARANTEE

With the exception of the guarantee described in these T&C, neither TALLADIUM ESPAÑA S.L. nor its representatives, nor third parties manufacturing or distributing the Products, represent or offer a guarantee, agreement or any other express or implicit, oral or written, commitment, with respect to the Products (without limitation), including guarantees involved in the marketing, durability or suitability for individual uses or purposes. In addition and within the maximum extent permitted by the relative law, TALLADIUM ESPAÑA S.L. rejects (on its own behalf, and on behalf of its representatives and third parties that manufacture or distribute Products) any responsibility with respect to any direct or indirect damage caused, which may result from or be a consequence of the design, composition of the dental prosthesis into which the Products are integrated.


GUARANTEE EXCLUSIONS

TALLADIUM ESPAÑA S.L. limits this guarantee to:

- Transformed abutments that form part of the dental prosthesis. But not the screws used to anchor them.
- Those products that are not used with the accessories and parts marketed by Talladium España

AMENDMENT OR SUSPENSION OF THE GUARANTEE

TALLADIUM ESPAÑA S.L. reserves the right to amend or withdraw these T&C at any time and without prior notification. Any modification or suspension shall not affect products already placed in patients.



This printed version is probably
outdated and may not include all
products or compatibilities.

**Scan and download the
latest version of this catalogue**

www.dynamicabutment.com



**DYNAMIC
ABUTMENT
SOLUTIONS**

Virginia Woolf, 17
25005 · Lleida (Spain)

INTERNATIONAL +34 873 450 709
das@dynamicabutment.com

SPAIN +34 973 289 580
spain@dynamicabutment.com



www.dynamicabutment.com

